

# Beer's- Lambert Law and Standard Curves

# Spectrophotometer:

- Spectrophotometer can be used to measure the amount of **light absorbed or transmitted** by a solution. **How?**
  - It consist of two parts: spectrometer and photometer.
    - a. **Spectrometer** is used for producing light of any selected color (wavelength)
    - b. **Photometer** is used for measuring the intensity of light.
  - By using the spectrophotometer, we can **quantitatively measure absorbance**, and this information can be used to determine the concentration of the absorbing molecule.
  - **More concentrated** solution will absorb **more light** and transmits less [why?]
- So:
- The more concentrated solution → **high** absorbance value.
  - And less concentrated solution → **less** absorbance value.


# Spectrophotometer cont.:

□ Wavelength in this instrument divided into:

1. Invisible range-ultraviolet- (from 100 to 360 nm)
2. Visible range (above 360 nm -700 nm)

💡 What are the type of cuvettes used in each range ? Why ?

□ **Blank:** contain everything except the compound to be measured



# Beer-Lambert law:

## □ Principle:

The absorption of light by a solution is described by the Beer-Lambert law as:

- There is **linear** relationship between **absorbance** and **concentration** of an absorbing species.

$$A = a_m \times c \times l$$

## □ Where:

- **A** = is the absorbance of the solution [Ab].
- **a<sub>m</sub>** = the molar extinction (absorption) coefficient.
- **l** = length of the light path through the solution.
- **c** = concentration of the absorbing substance.

# Beer-Lambert law cont.:

□ From the law we observe :

1. **Direct** relationship between **C** and **A**.



$$A = a_m \times c \times l$$

2. **Direct** proportional between **I** and **A**.

➤ If we plot absorbance against concentration, we get a straight line.

➤ The **linear relationship** between **concentration** and **absorbance** is both simple and straightforward.

➔ Which is why we prefer to express the **Beer-Lambert law** using **absorbance** as a measure of the absorption rather than %T ( % transmittance).




# Standard curve for concentrations:

- It is a **graph** that shows the relationship between different known concentrations of a substance and the absorbance at a specific wavelength.
- **Standard curves** are most commonly used to determine the concentration of a substance, using serial dilution of solutions (standard solutions) of known concentrations.

## ➔ What is a standard solution?

- Is a solution containing a **precisely known concentration** of an element or a substance.
- A series of known standard solutions can be prepared by diluting the stock known solution.

(We should calculate the concentration of the diluted solution by the formula:  $C_1 \times V_1 = C_2 \times V_2$ )



# Some points to consider:

1. The absorbance of a solution with “unknown concentration” preferred to be **lower** than the highest Ab value in the standard curve.  
  
→ [ the absorbance of solution with “unknown concentration”, is within the range of absorbance values of solution with “known concentration solutions”].
2. If your unknown sample had an absorbance higher than the highest absorbance recorded by the standard [out of the range], **how would you determine its concentration correctly?**

## Two choices:

- I. **Increase** the concentration of standard solutions.
  - II. **Dilute** the solution with “unknown concentration” → measure the Ab after dilution → determine its concentration from the curve → then multiply the value by the **dilution factor**.
3. If there is insufficient volume the light may pass over the solution instead of going through it.
  4. Absorbance has **no units**.

# Determining the unknown concentration of a solution by known absorbance value



## From standard curve

- The standard is constructed by plotting the **absorbance** values vs **concentration** and solving for the best straight line which is given in terms of  $y = mx + b$ 
  - where; y is absorbance, x concentration, b is the y-intercept and m is the slope.
- Measure the absorbance of the “solution with unknown concentration” to determine the concentration.

## Beer-Lambert law

- Using available information of any standard solution to determine the “ $\epsilon$ ”
- Then using these information to get the unknown concentration using:  $A = \epsilon lc$
- Note: “ $\epsilon$ ” will changed when the wavelength changed.



# Example:

Determine the concentration of a solution with an absorbance value of **0.041** from the standard curve.

□  $y = mx + b$  (Where:  $y$  is absorbance,  $x$  concentration,  $b$  is the y-intercept and  $m$  is the slope)

Based on the curve, the slope ( $m$ ) = 0.8 and ( $b$ ) = 0

Then:

$$0.041 = 0.8(x) + 0$$

$$X = 0.041/0.8 = 0.051M$$

Or use the TREND formula




Figure 1. Standard curve of concentration of 0.1 M of copper sulphate

# Practical Part

# Objectives:

---

- To understand the concept of Beer-Lambert law and its application.
- Getting familiar with the standard curve.
- Determination of an unknown concentration for two solutions.

# Method:

How to determine the **concentrations of the standard solutions** ?

You are provided by:

1. Standard solution (Stock solution) of Copper Sulfate with 0.1 M [known concentration].
2. Solution with Unknown concentration "A".
3. Solution with Unknown concentration "B".

Set up 8 test tubes, as following table:

Tube	0.1M Copper Sulfate Standard Solution (ml)	Distal water (ml)	Solutions with unknown concentration (ml)
Blank	-	5 ml	-
A	1ml	4 ml	-
B	2 ml	3 ml	-
C	3 ml	2 ml	-
D	4 ml	1 ml	-
E	5 ml		-
Solution "A"	-	-	2 ml
Solution "B"	-	-	2 ml

Mix the contents using the vortex.

Measure the absorbance of each tube at 600 nm against the blank.

Wear chemical splash goggles and chemical resistant clothing such as gloves and aprons when handling copper sulphate

# Results:

- Calculate the concentrations of standard solutions (Tube A-E).
- Plot the standard curve (Absorbance vs. Concentration).
- Determine the concentration of Solution "A" and "B" from the standard curve.
- Calculate the concentration of Solution "A" and "B", from the law using ( $0.942 \text{ M}^{-1}\text{cm}^{-1}$ ).
- Calculate the extinction coefficient of your Copper Sulfate solution (using one of the known tubes).

<b>Tube</b>	<b>Absorbance at 600nm</b>	<b>Concentration M</b>
<b>A</b>		
<b>B</b>		
<b>C</b>		
<b>D</b>		
<b>E</b>		
<b>Solution "A"</b>		From the curve=
<b>Solution "B"</b>		From the curve=