

**King Saud University
Faculty of Science
Chemistry Department**

**General Practical Chemistry
EXPERIMENTS REPORTS
101 Chem & 104 Chem**

**Text Book:
Practical General Chemistry
By
Dr. Ahmad Al-Owais & Dr. Abdulaziz Al-Wassil**



LIST OF EXPERIMENTS

Week	Date	Experiment	Book Page
1		Determination of a Liquid Density	131 - 133
2		Preparation of a Standard Solution of Sodium Carbonate	63 - 86
3		Determination of Organic Indicators for Acid Base Titrations	95 - 98
4		Determination of Sodium Hydroxide Concentration By Titrations With A Standard Solution of Hydrochloric Acid +Quiz 1	99 - 109
5		Determination of Acetic Acid Concentration By Titrations With A Standard Solution of Sodium Hydroxide	99 - 109
6		Determination of Hydrochloric Acid Concentration By Titrations With A Standard Solution of Sodium Carbonate	87 - 93
7		First Exam	
8		Measurement of Gas Diffusion (Graham's Law of Diffusion)	137 - 140
9		Determination of Critical Solution Temperature+ Quiz 2	177 - 180
10		Hess's Law	197 - 203
11		Effect of Concentration on Reaction Rate	183 - 190
12		Determination of The Molar Mass of An Organic Compound by The Depression of Its Freezing Point	205 - 216
13		Determination of The Molar Mass of An Organic Compound by The Steam Distillation+ Quiz 3	191 - 196
14		Second Exam	



EXPERIMENT (1)

Determination of a Liquid Density

DATE:

STUDENT'S NAME:

STUDENT'S NUMBER:

Symbols


Mass of the empty beaker in g = m_1

Mass of the beaker and the liquid in g = m_2

Mass of the liquid ($m_2 - m_1$) in g = m

Volume of the liquid in $\text{cm}^3 = V$

Density of liquid ($\frac{m}{V}$) in $\text{g} / \text{cm}^3 = d$



Notes:

The liquid used in this experiment is:

Results and Calculations:

- Calculate the liquid's density in all cases and put in the following table:

V (cm³)								
m₂ (g)								
m₁ (g)								
m (g)								
d (g / cm³)								

- Plot the relationship between the mass of the liquid (m) on the Y-axis versus its volume (V) on the X-axis, and find the liquid's density from the slope.

$$\text{slope} = d = \frac{\Delta m}{\Delta V} \frac{-}{-} = \frac{-}{-} = \underline{\hspace{2cm}} =$$

Experiment (1)
Graphical relation between mass of liquid (m) and its volume (V)



EXPERIMENT (2)

Preparation of a Standard Solution of Sodium Carbonate

DATE:

STUDENT'S NAME:

STUDENT'S NUMBER:

Molar masses (g mol^{-1}): C = 12 , O = 16 , Na = 23



Results & calculation:

Molarity of sodium carbonate standard solution = $C_{\text{Na}_2\text{CO}_3} =$ mol L^{-1}

Molar mass of $\text{Na}_2\text{CO}_3 = M_{\text{Na}_2\text{CO}_3} =$ g mol^{-1}

Volume of solution in L = $V =$ L

Number of mole of $\text{Na}_2\text{CO}_3 = n_{\text{Na}_2\text{CO}_3} = C_{\text{Na}_2\text{CO}_3} \times V = \times =$ mol

Mass of $\text{Na}_2\text{CO}_3 = m_{\text{Na}_2\text{CO}_3} = n_{\text{Na}_2\text{CO}_3} \times M_{\text{Na}_2\text{CO}_3} = \times =$ g

EXPERIMENT (3)

Determination of Organic Indicators for Acid Base Titrations

DATE:

STUDENT'S NAME:

STUDENT'S NUMBER:

A) Titration of a strong acid (HCl) with a strong base (NaOH) using the pH meter

Molar masses (g mol^{-1}): H = 1 , O = 16 , Na = 23

Results & calculation:

Volume of base added (V_{base})	0	5	10	15	20	22.5	24	24.5	24.8	25	26	28	30
Calculated pH	1.2	1.3	1.4	1.6	1.9	2.3	2.8	3.3	3.6	9.7	11	11.4	11.6
Measured pH													



HCl molarity = $M =$ molar

HCl volume = $V =$ mL

NaOH molarity = $M' =$ molar

NaOH volume (from diagram) = $V' =$ mL

- Knowing that the pH range for methyl orange indicator (M.O.) is from 3.1 to 4.4, and for phenol phthaline indicator (Ph.Ph.) is from 8 to 10, plot pH (on the Y-axis) against V_{base} (on the X-axis).



From the graph:

- The pH range at the equivalent point is from () to ().
- The suitable indicator for this titration is ().

2. Calculation of the base molarity :

3. Calculation of the base concentration in g L^{-1} :

4. Calculation of pH, pOH, $[\text{H}^+]$ and $[\text{OH}^-]$:

	HCl solution	NaOH solution
pH		
pOH		
$[\text{H}^+]$		
$[\text{OH}^-]$		

Experiment (3)

**A) Graphical relation between the (pH) and the volume of base added (V)
(Titration of a strong acid with a strong base)**



EXPERIMENT (3)

Determination of Organic Indicators for Acid Base Titrations

DATE:

STUDENT'S NAME:

STUDENT'S NUMBER:

B) Titration of a weak acid (CH₃COOH) with a strong base (NaOH) using the pH meter

Molar masses (g mol⁻¹): H = 1 , C = 12 O = 16
--

Results & calculation:

Volume of base added (V _{base})	0	2	4	6	8	10	12	14	16	17	17.5	18	18.5
Calculated pH	1.2	1.3	1.4	1.6	1.9	2.3	2.8	3.3	3.6	9.7	11	11.4	11.6
Measured pH													

CH₃COOH molarity = M = molar



CH₃COOH volume = V = mL

NaOH molarity = M' = molar

NaOH volume (from diagram) = V' = mL

Requirements:

- Knowing that the pH range for methyl orange indicator (M.O.) is from 3.1 to 4.4, and for phenol phethaline indicator (Ph.Ph.) is from 8 to 10, plot pH (on the Y-axis) against V_{base} (on the X-axis).



Created with

From the graph:

- The pH range at the equivalent point is from () to ().
- The suitable indicator for this titration is ().

5. Calculation of the base molarity :

6. Calculation of the base concentration in g L^{-1} :

7. Calculation of pH, pOH, $[\text{H}^+]$ and $[\text{OH}^-]$:

	HCl solution	NaOH solution
pH		
pOH		
$[\text{H}^+]$		
$[\text{OH}^-]$		

Experiment (3)

B) Graphical relation between the (pH) and the volume of base added (V)
(Titration of a weak acid with a strong base)



EXPERIMENT (4)

Determination of Sodium Hydroxide Concentration By Titrations With A Standard Solution of Hydrochloric Acid

DATE:

STUDENT'S NAME:



STUDENT'S NUMBER:

Molar masses (g mol^{-1}): H = 1 , O = 16 Na = 23

Results:

FIRST: Volume of NaOH using Ph.Ph. as indicator:

Exp.	Initial reading	Final reading	Volume (V) mL	Average
1				
2				
3				



Calculations:

1. Volume of NaOH = $V =$ mL

2. Volume of HCl = $V' =$ mL

3. Molarity of HCl = $M' =$ mol L⁻¹

4. The reaction equation is:


5. Calculation of the base molarity:

- Ph.Ph. indicator used is ().
- pH range of indicator is from () to ().
- At the end point the color of indicator changed from () to ().
- From the reaction equation using Ph .Ph. as indicator:
n =
n' =

6 Calculation of the base concentration in g L⁻¹:

SECOND: Volume of NaOH using M.O. as indicator:

Exp.	Initial reading	Final reading	Volume (V) mL	Average
1				
2				
3				



Calculations:

1. Volume of NaOH = $V =$ mL

2. Volume of HCl = $V' =$ mL

3. Molarity of HCl = $M' =$ mol L⁻¹

4. The reaction equation is:

5. Calculation of the base molarity:

- Indicator used is ().
- pH range of indicator is from () to ().
- At the end point the color of indicator changed from () to ().
- From the reaction equation using M.O. as indicator:
 $n =$
 $n' =$

6. Calculation of the base concentration in g L⁻¹:

EXPERIMENT (5)

Determination of Acetic Acid Concentration By Titrations With A Standard Solution of Sodium Hydroxide

DATE:

STUDENT'S NAME:


STUDENT'S NUMBER:

Molar masses (g mol^{-1}): H = 1 , C = 12 O = 16
--

Results:

FIRST: Volume of NaOH using M.O. as indicator:

Exp.	Initial reading	Final reading	Volume (V) mL	Average
1				
2				
3				




Calculations:

1. Volume of NaOH = $V =$ mL
2. Volume of $\text{CH}_3\text{COOH} = V' =$ mL
3. Molarity of NaOH = $M' =$ mol L^{-1}
4. The reaction equation is:
5. Calculation of the acid molarity:
 - Indicator used is ().
 - pH range of indicator is from () to ().
 - At the end point the color of indicator changed from () to ().
 - From the reaction equation using M.O. as indicator:
 $n =$
 $n' =$
6. Calculation of the acid concentration in g L^{-1} :

SECOND: Volume of NaOH using Ph.Ph. as indicator:

Exp.	Initial reading	Final reading	Volume (V) mL	Average
1				
2				
3				



Calculations:

1. Volume of NaOH = $V =$ mL
2. Volume of CH₃COOH = $V' =$ mL
3. Molarity of NaOH = $M' =$ mol L⁻¹
4. The reaction equation is:
5. Calculation of the acid molarity:
 - Indicator used is ().
 - pH range of indicator is from () to ().

- At the end point the color of indicator changed from () to ().
- From the reaction equation using M.O. as indicator:
 $n =$
 $n' =$

6. Calculation of the acid concentration in g L^{-1} :

EXPERIMENT (6)

Determination of Hydrochloric Acid Concentration By Titrations With A Standard Solution of Sodium Carbonate

DATE:

STUDENT'S NAME:


STUDENT'S NUMBER:

Molar masses (g mol^{-1}): H = 1 , Cl = 35.45
--

Results:

Volume of HCl using M.O. as indicator:

Exp.	Initial reading	Final reading	Volume (V) mL	Average
1				
2				
3				



Calculations:

1. Volume of HCl = $V =$ mL
2. Volume of $\text{Na}_2\text{CO}_3 = V' =$ mL
3. Molarity of $\text{Na}_2\text{CO}_3 = M' =$ mol L^{-1}
4. The reaction equation is:
5. Calculation of the acid molarity:
 - Indicator used is ().
 - pH range of indicator is from () to ().
 - At the end point the color of indicator changed from () to ().
 - From the reaction equation using M.O. as indicator:
 $n =$
 $n' =$
6. Calculation of the acid concentration in g L^{-1} :

EXPERIMENT (7)


Measurement of Gas Diffusion (Graham's Law of Diffusion)

STUDENT'S NAME:

STUDENT'S NUMBER:

Molar masses (g mol^{-1}): H = 1 , N = 14 Cl = 35.45

$$\frac{r_{\text{NH}_3}}{r_{\text{HCl}}} = \frac{\sqrt{d_{\text{HCl}}}}{\sqrt{d_{\text{NH}_3}}} = \frac{\sqrt{M_{\text{HCl}}}}{\sqrt{M_{\text{NH}_3}}} \quad (\text{Graham's law})$$



$$\frac{r_{\text{NH}_3}}{r_{\text{HCl}}} = \frac{\frac{L_{\text{NH}_3}}{\text{time}}}{\frac{L_{\text{HCl}}}{\text{time}}}$$

$$\frac{L_{\text{NH}_3}}{L_{\text{HCl}}} = \frac{\sqrt{M_{\text{HCl}}}}{\sqrt{M_{\text{NH}_3}}}$$

Results:

1. Distance moved by HCl gas (L_{HCl}) = _____ cm
2. Distance moved by NH_3 gas (L_{NH_3}) = _____ cm
3. Reaction equation:

Calculations:

1. The theoretical ratio between the molar masses of the two gases (Y):
2. The measured ratio between the molar masses for the two gases (X):
3. The practical molar mass of one of the two gases (M_Y) knowing the theoretical molar mass of the other gas and the values of L_{HCl} and L_{NH_3} using Graham's law:

$$\frac{L_{NH_3}}{L_{HCl}} = \frac{\sqrt{M_{HCl}}}{\sqrt{M_{NH_3}}}$$

- 4 Calculation of the theoretical molar mass of the same gas using the molar masses of its atoms (M_X):
- 5 Error percentage:
 - First method:

$$\text{Error percentage} = \pm \frac{\text{difference between theoretical and practical ratios}}{\text{theoretical ratio}} \times 100$$

$$\text{Error percentage} = \pm \frac{Y - X}{Y} \times 100$$

- Second method

$$\text{Error percentage} = \pm \frac{\text{difference between theoretical and practical molar masses}}{\text{theoretical molar mass}} \times 100$$

$$\text{Error percentage} = \pm \frac{M_Y - M_X}{M_X} \times 100$$

EXPERIMENT (8)

Determination of Critical Solution Temperature

DATE:

STUDENT'S NAME:

STUDENT'S NUMBER:

Results:

Experimental results and calculations:

EXP. No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14
Mass of phenol (g)	4	4	4	4	4	4	4	4	4	4	4	4	4	4
Mass of water (g)	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Mass of solution (g)	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Mass % of water														
Mass % of phenol														
Miscibility temperature (°C)														

Calculations:

1. A SAMPLE CALCULATION OF A MASS PERCENTAGE

Mass % of water = $\frac{\text{mass of water}}{\text{mass of mixture}} \times 100$	Mass % of phenol = $\frac{\text{mass of phenol}}{\text{mass of mixture}} \times 100$

2. From the graphical relation between miscibility temperature (Y-axis) and the mass percentage of phenol (X-axis):

• The critical solution temperature (C.S.T) = °C

• Mass % of phenol =

• Mass % of water =

Experiment (8)

Graphical relation between miscibility temperature and the mass percentage of phenol



EXPERIMENT (9)


Hess's Law

DATE:

STUDENT'S NAME:

STUDENT'S NUMBER:

Diagrammatic illustration of Hess's law:



General Notes:

1. Density of NaOH solution = $d = 1 \text{ g / cm}^3$
2. Specific heat of NaOH solution = $\rho_{\text{solution}} = 4.18 \text{ J / g } ^\circ\text{C}$
3. Specific heat for calorimeter (glass) = $\rho_{\text{calorimeter}} = 0.836 \text{ J / g } ^\circ\text{C}$
4. Volume of NaOH solutions used in all experiment = $V = 50 \text{ mL}$
5. Mass of NaOH solutions in g = $m_{\text{solution}} = V \times d = 50 \times 1 = 50 \text{ g}$
6. Mass of calorimeter (glass tube) in g = $m_{\text{calorimeter}}$
7. Initial temperature in $^\circ\text{C} = t_1$
8. Final temperature in $^\circ\text{C} = t_2$
9. Temperature change in $^\circ\text{C} = \Delta t$
10. Heat gained by solution in J = $q_1 = \rho_{\text{solution}} \times m_{\text{solution}} \times \Delta t$
11. Heat gained by calorimeter in J = $q_2 = \rho_{\text{calorimeter}} \times m_{\text{calorimeter}} \times \Delta t$
12. Total heat gained in J = $Q = q_1 + q_2$
13. Number of moles of NaOH used in:
 - experiment 1 = $(n_{\text{NaOH}})_1 = \frac{(m_{\text{NaOH}})_1}{M_{\text{NaOH}}}$
 - experiment 2 = $(n_{\text{NaOH}})_2 = \frac{(m_{\text{NaOH}})_2}{M_{\text{NaOH}}}$
 - experiment 3 = $(n_{\text{NaOH}})_3 = (\text{molarity})_{\text{NaOH}} \times V_{\text{NaOH}}$

14. $\Delta H = \frac{Q}{n_{\text{NaOH}}}$

Calculations and results:

	Experiment 1	Experiment 2	Experiment 3
t_1 (°C)			
t_2 (°C)			
Δt (°C)			
q_1 (J)			
q_2 (J)			
Q (J)			
n_{NaOH} (mol)			
ΔH (kJ mol ⁻¹)			

Verification of Hess's law using thermochemical equations:

- 1)
- 2)
- 3)

EXPERIMENT (10)

Effect of Concentration on Reaction Rate

(Determination of the order of the sodium thiosulphate and hydrochloric acid reaction)

DATE:

STUDENT'S NAME:


STUDENT'S NUMBER:

Reaction equation:

Rate law:

Arrhenius equation is:

The graphical plot:



Determination of the reaction order with respect to sodium thiosulphate:

- Symbols

1. Volume of $\text{Na}_2\text{S}_2\text{O}_3$ before dilution = $V_{\text{Na}_2\text{S}_2\text{O}_3}$
2. Volume of $\text{Na}_2\text{S}_2\text{O}_3$ After dilution = $V'_{\text{Na}_2\text{S}_2\text{O}_3} = 29 \text{ mL}$
3. Molarity of $\text{Na}_2\text{S}_2\text{O}_3$ before dilution = $M_{\text{Na}_2\text{S}_2\text{O}_3} = 0.15 \text{ mol L}^{-1}$
4. Molarity of $\text{Na}_2\text{S}_2\text{O}_3$ after dilution = $M'_{\text{Na}_2\text{S}_2\text{O}_3} = \frac{M_{\text{Na}_2\text{S}_2\text{O}_3} \times V_{\text{Na}_2\text{S}_2\text{O}_3}}{V'_{\text{Na}_2\text{S}_2\text{O}_3}}$
5. Reaction time in seconds = t
6. Reaction rate in seconds⁻¹ = $\frac{1}{t}$

- **Calculatoinis**

1. Calculation of $M'_{Na_2S_2O_3}$ in the reactions of $Na_2S_2O_3$ with HCl in 29 mL solution:

- $(M'_{Na_2S_2O_3})_1 =$

- $(M'_{Na_2S_2O_3})_2 =$

- $(M'_{Na_2S_2O_3})_3 =$

- $(M'_{Na_2S_2O_3})_4 =$

- $(M'_{Na_2S_2O_3})_5 =$

Exp.	$V_{Na_2S_2O_3}$	V_{H_2O}	V_{HCl}	$M'_{Na_2S_2O_3}$	t	$\frac{1}{t}$	$-\log M'_{Na_2S_2O_3}$	$-\log \frac{1}{t}$
1	25	0	4					
2	20	5	4					
3	15	10	4					
4	10	15	4					
5	5	20	4					

2. Obtaining the order, n, from the plot of $\log \frac{1}{t}$ versus $\log M'_{Na_2S_2O_3}$ according to:

$$\log \frac{1}{t} = \log k + n \log M'_{Na_2S_2O_3}$$

slope =

n =

Determination of the reaction order with respect to hydrogen chloride:

- **Symbols**

1. Volume of HCl before dilution = V_{HCl}
2. Volume of HCl after dilution = $V'_{HCl} = 15 \text{ mL}$
3. Molarity of HCl before dilution = $M_{HCl} = 1 \text{ mol L}^{-1}$
4. Molarity of HCl after dilution = $M'_{HCl} = \frac{M_{HCl} \times V_{HCl}}{V'_{HCl}}$
5. Reaction time in seconds = t
6. Reaction rate in seconds⁻¹ = $\frac{1}{t}$

- **calculatins**

1. Calculation of M'_{HCl} in the reactions of HCl with $Na_2S_2O_3$ in 15 mL solution:

- $(M'_{\text{HCl}})_1 =$
- $(M'_{\text{HCl}})_2 =$
- $(M'_{\text{HCl}})_3 =$
- $(M'_{\text{HCl}})_4 =$
- $(M'_{\text{HCl}})_5 =$

Exp.	V_{HCl}	$V_{\text{H}_2\text{O}}$	$V_{\text{Na}_2\text{S}_2\text{O}_3}$	M'_{HCl}	t	$\frac{1}{t}$	$-\log M'_{\text{HCl}}$	$-\log \frac{1}{t}$
1	5	0	10					
2	4	1	10					
3	3	2	10					
4	2	3	10					
5	1	4	10					

2. Obtaining the order, n, from the plot of $\log \frac{1}{t}$ versus $\log M'_{\text{HCl}}$ according to the equation:

$$\log \frac{1}{t} = \log k + m \log \log M'_{\text{HCl}}$$

slope =

n =

Rate law:

Rate constant:

k =

Experiment (10)

Graphical relation between $\log \frac{1}{t}$ and $\log M_{\text{Na}_2\text{S}_2\text{O}_3}$



Experiment (10)

Graphical relation between $\log \frac{1}{t}$ and $\log M'_{HCl}$



EXPERIMENT (11)

Determination of the Molar Mass of An Organic Compound By The Depression of Its Freezing Point

DATE:

STUDENT'S NAME:

STUDENT'S NUMBER:

Results:

	Unknown A	Unknown B
Mass of solvent m_1 (g)		
Mass of solute m_2 (g)		
t_{solvent} ($^{\circ}\text{C}$)		
t_{solution} ($^{\circ}\text{C}$)		
$\Delta t_f = t_{\text{solvent}} - t_{\text{solution}}$ ($^{\circ}\text{C}$)		

Calculations:

- Molal freezing point depression constant of solvent = $K_f = 1.86^{\circ}\text{C molal}^{-1}$
- Molar mass of solute in $\text{g mol}^{-1} = M_2 = K_f \frac{m_2 \times 1000}{\Delta t_f \times m_1}$

$$(M_2)_A =$$

$$(M_2)_B =$$

EXPERIMENT (12)


Determination of the Molar Mass of An Organic Compound By The Steam Distillation

DATE:

STUDENT'S NAME:

STUDENT'S NUMBER:

Results:



1. Volume of water after distillation = $V_{H_2O} =$ cm^3
2. Density of water = $d_{H_2O} = 1 \text{ g cm}^{-3}$
3. Volume of unknown liquid after distillation $V_B =$ cm^3
4. Density of unknown liquid = $d_B = 1.106 \text{ g cm}^{-3}$
5. Atmosphere pressure in Riyadh $P^\circ = 720 \text{ mmHg}$
- 6.

P°_{water} (mmHg)	489.8	504.7	526.0
$T_{b, \text{water}}$ ($^\circ\text{C}$)	88	89	90

7. $P_{\text{total}} = P^\circ = 720 \text{ mmHg}$

Calculations:

1. Calculation of the unknown vapor pressure:

2. Calculation of the unknown molar mass: