



BCH 312

Biochemical Calculations [Practical]

Nora Aljebrin


Office: Building 5, 3rd floor, 5T267

Website: www.fac.ksu.edu.sa/naljebrin

E.mail: naljebrin@ksu.edu.sa




Mark Distribution



- Report [5 Marks]
- Quiz [7 Marks]
- Final [13 Marks, 8 for practical and 5 for theoretical]


- Final exam date:
 - Monday 30/2/1435H - 22/12/2014



+ Identification of Glassware, pH meter, and Spectrophotometer

+ Objectives

- Identify different glassware and the accuracy of them
- To be familiar with the use of pipetting technique
- To learn how to handle the pH meter and to measure pH values
- To learn how to handle the spectrophotometer




+ Glassware



■ Volumetric flask

- It is used for preparing large amounts of standard solutions and reagents that require highly accurate concentrations.
- The most common volumes are 250, 500, 1000, and 2000 ml. However smaller volumes are available such as 100, 50, 25, and even 1, 2, and 5 ml.





Ehrlenmeyer flasks


- This type of flask is a conical container with multiple volume markings that serve as an estimate of volume and is available in a variety of sizes, such as 250, 500 and 1000 ml.
- It is used in solution preparation but are less accurate than volumetric flasks.
- The purpose of an Ehrlenmeyer flask is to help dissolve a solid solute into solution before transfer to a volumetric flask for final volume adjustment





■ Graduated cylinder

- It is a tall flask with multiple volume gradation.
- The purpose of this flask is to aliquot volumes of a fluid in making up a reagent dilution or to determine the volume of an unknown fluid.




■ Beakers

- They are used to hold stock solutions or diluents for short-term storage or during laboratory procedures.
- They are available in a variety of total volume sizes, but have the least accuracy of the three types of flasks.



+



Burette and its stand



Glass rod



Funnel



FOLD FILTER
PAPER IN HALF



FOLD INTO QUARTERS WITH TOP
SECTION SMALLER THAN BOTTOM



TEAR OFF CORNER OF
SMALLER SECTION



OPEN CONE

Filter paper



Pipettes

- There are two main type of pipettes are used in biochemical laboratory:

a) Volumetric or transfer pipettes

- They are designed to deliver a fixed volume of liquid
- It is consists of a cyclindrical bulb joined at both ends to narrowed glass tubing.
- Non-blown out
- More accurate than measuring pipettes


b) Graduated or measuring pipettes

- These consists of a plain narrow tube drawn out to a tip and graduated uniformly along its length.
- Some are blown out
- There are two types of graduated pipettes are available:
 - 1-Mohr, graduated between two marks
 - 2-Serological pipettes with graduation marks down to the tip.

+ Pipette Types




Volumetric
pipette



Graduated
pipettes


+ How to use pipettes




- The pipette is first washed with water ,then rinsed several times with a little of the solution to be used and finally filled to just above the mark , the liquid is allowed to fall to the mark .
- The solution is allowed to drain into the appropriate vessel with the jet of the pipette touching the wall of the vessel .
- After the flow of the liquid has stopped, the jet is held against the wall for some times and then removed .





1. Use mechanical suction



2. Wipe off outside of pipette with gauze




3. Adjust the meniscus; read meniscus at the bottom of the curved liquid






Smallest division of graduated pipette



- How to know the smallest division of a pipette:



Smallest division:

0.1 ml

Smallest division:

0.2 ml

Smallest division:

0.02 ml

Smallest division:

0.02 ml



Method



Part 1

- Examine the 3 pipette placed on your laboratory bench .
- Record their types and the volume of their smallest division .

Pipette	Type	Smallest division
A		
B		
C		



Part 2

- By using distilled water , pipette into weighted beaker
 - with 5ml graduate pipette(Mohr) =5ml water
 - with 5ml measuring cylinders =5ml water

	Weight of the beaker	Weight of beaker +water	Weight of water
graduate pipette			
measuring cylinder			

+ Laboratory Equipment


■ pH Meter:



- Hydrogen ion concentration of many solution is low and difficult to measure accurately.
- So, the term pH introduced as a way of expressing hydrogen ion concentration
- pH define as the negative logarithm of the hydrogen ion concentration .

→ **pH = - log₁₀ [H⁺]**


- pH range value (0 – 14), the higher pH number, the lower the hydrogen ion concentration and vice versa





Glass electrode

- The glass electrode consists of a very thin bulb, blown onto a hard glass tube.
- The bulb is made of high conductivity glass which is sensitive to pH.
- The bulb contains a solution of hydrochloric acid (0.1N) and is connected to a platinum lead via silver - silver chloride electrode which is reversible with respect to hydrogen ions.
- The glass electrode is very sensitive and readily responds to changes in hydrogen ion concentration.






Method


- Standardize the PH meter by placing the electrode in a solution of known PH(PH 4 , 7 , 9) .
- Wash the electrode with distilled water and dry with a tissue then put it into sample solution A & B , read PH .
- **Note:** After use the electrode you should storage it in distilled water and never be allowed to dry out . If the electrode get dry it will required reactivation.

Solution	pH Value
Standard 4	
Standard 7	
Standard 9	
A	
B	

+ Spectrophotometer

- It is an instrument used to measure the intensity of light at a given wavelength that is transmitted or absorbed by a sample.
- It consists of two parts:
 - **Spectrometer** is designed to emit the light at different wavelengths
 - **Photometer** contains a photoelectric cell and the potentials are recorded on a scale which read out as absorbance or transmittance.





Wavelength in this instrument is divided into:

- Invisible range (ultraviolet) from 100 to 360 nm → [Quartz cuvette are used]
- Visible range (400 -700 nm) → [Glass or plastic cuvette are used]
- Blank: It contains everything except the compound to be measured.

+ Method

- Adjust the spectrophotometer to zero using blank solution in the cuvette and read the absorbance of standard solution and the solution of unknown concentration at 280 nm.
- Read your result in the table below:

NO.	Solution	Absorbance
1	Standard solution (0.5 gm/100 ml of BSA)	
2	Solution of Unknown concentration	



- Calculate the concentration of unknown solution from the following formula:

$$C_u = (A_u / A_s) \times C_s$$

Where:

A_u = Absorbance of the solution of unknown concentration

A_s = Absorbance of standard solution

C_s = Concentration of standard solution

C_u = Concentration of the solution of unknown concentration