

Preparation of Biological Solutions and Serial Dilutions



- Objective:

1- To learn how to prepare solutions.

2-To get familiar with solution dilutions.



- Introduction:


- It is very important to understand how to prepare solutions and make dilutions and it is an essential skill for biochemists which is necessary knowledge needed for doing any experiment.
- A simple solution is basically two substances that are mixed together.
- One of them is called the solute and the other is the solvent. A **solute** is the substance to be dissolved (sugar). A **solvent** is the one doing the dissolving (water) forming a homogenous mixture.



DISSOLVE THE SOLUTE
IN THE SOLVENT



SOLUTE



SOLVENT



A. Preparation of biological solutions:

- There are many units for concentration:

1. Molarity.

2. W/V %.

3. W/W %.



1.Molarity:

- It is the number of moles of solute dissolved in one liter of solution.

Molar = no. of mole/vol. in L.

- *Remember that:*

- No of mole = weight (g) / molecular weight

- **Example: 0.5 Molar (M) solution:**

- That mean there are 0.5 mole dissolved in 1000ml (1L).

Note: Number of mole (No of mole)



Example:

- **How many gram to Prepare 2 M of NaCl in 100 ?**
- 2 mole of NaCl present in 1000 ml [or 1Liter] of solvent (dis.H₂O)
- And we know that

No of mole = weight (g) / molecular weight.

So, $[2 \text{ mole} = \text{weight (g)} / 58.5] \longrightarrow \text{weight (g)} = 2 \times 58.5 = 117 \text{ g.}$

- This weight needed if 1000 ml is required to be prepared. Since we need to prepare only 100 ml.

117 g ~~====>~~ 1000 ml.
? g ~~====>~~ 100 ml.

$$[(100 \times 117)/1000] = 11.7 \text{ g}$$

- **11.7 g of NaCl** dissolved in small volume of dis.H₂O, then complete the volume up to 100 ml.

Note: The MW of NaCl is 58.44 =(35.5+23)



- Practically how to prepare 2M NaCl:

1. Place a beaker in a balance and zero the balance.
2. Weight 11.7 grams of NaCl , in the beaker and dissolve in little water (less than 100 ml).
3. Once the solid is dissolved the volume is transferred to 100 ml volumetric flask.
4. Brought up to a final volume 100 ml.



2. W/V %:Weight/Volume Percentage Concentration:

- The number of grams of solute dissolved in 100 mL of solution.

For example:

- Prepare 3% of NaOH:

It means 3 grams of NaOH is dissolved in 100 ml of the solution.

- Example:

Prepare 50 ml of 4% NaOH

4g-----> 100 ml

?----->50 ml

- The weight in grams of NaOH needed to prepare 4% NaOH is
 $= (4 \times 50)/100 = 2 \text{ g.}$

So,

- 2 grams of NaOH is dissolved in little water and the volume made up to 50 ml



3. W/W %:

- The number of grams of solute dissolved in 100 gram of solution.
- The concentrations of many commercial acids are giving in terms of w/w%.
- In order to calculate the volume of the stock solution required for a given preparation the density (**specific gravity**) of stock solution should be provided.
- **Weight (wt) = volume (ml) x SG x w/w% (as decimal).**
- To calculate $w/w\%$ as decimal = $(w/w)/100$.
- For example: $w/w\% = 20\% \longrightarrow 20 / 100 = 0.20$



- Prepare 100 ml with 0.4 M HCl solutions starting with the concentrated HCl solution you are provided with. (w/w% = 36% , S.Gr = 1.15).

- **Weight= volume(ml) x SGr x w/w% (as decimal)**

• **Important Note!!!:** the volume in this formula is not the required volume in the question, it is the volume of the concentrated HCl that you must add.

- First we must calculate the weight by the following:

Mole= Molarity x volume in liter = $0.4 \times 0.1 = 0.04$ mole

Weight= mole x MWt $\longrightarrow 0.04 \times 36.5 = 1.46$ g


$1.46 = \text{volume} \times 1.15 \times 0.36 \longrightarrow \text{Volume} = \mathbf{3.53\text{ ml}}$

- **So, 3.53 ml of stock** (i.e. concentrated HCl) solution is needed and the volume made up to **100 ml by the addition of water.**



B .Dilution of Solution :

- It means to add more solvent without the addition of more solute →
To make it less concentrated.



1. Volume to volume dilutions (ratio).
2. Preparing dilutions by using the $V_1 \times C_1 = V_2 \times C_2$ formula.
3. Serial Dilutions

1:4 dilution ?

vol. Of solute

Total vol.

1 ml from solute + 3 ml from **solute** = total volume 4



1) Volume to volume dilutions (ratio):

- This type of dilutions describes the ratio of the solute to the final volume of the dilute solution.
- **For example**, to make 1:10 dilution of a 1.0 M NaCl solution:
one part of the 1.0 M NaCl solution, should be mixed with nine parts of water, for a total of ten parts,
- Therefore 1 :10 dilution means 1 part + 9 parts of water .
- **Thus if 10 ml of the 1:10 dilution was needed**, then 1ml of 1.0 M NaCl should be mixed with 9 ml of water.



- If 100 ml of 1:10 dilution was needed, then 10 ml of the 1.0 M NaCl should be mixed with 90 ml of water. The final concentration of NaCl in both cases will be 0.1 M.

Example:

- Prepare 2:10 dilution of solution (A) with 7M , but the total volume is 20ml not 10 ml?

- Dilution factor (D.F) = $10/2 = \underline{5}$


2 ml \rightarrow 10 ml

?? \rightarrow 20 ml

$$= (2 \times 20) / 10 = 4 \text{ ml}$$

So,

- 4 ml from solution(A) of 7M is needed and complete volume up to 20 ml (adding 16 ml water).

Note: [16 ml water = 4 ml - 20 ml]. 

(2) Preparing dilutions by using the $V_1C_1=V_2C_2$ formula:

- Sometimes it is necessary to use one solution to make a specific amount of a more dilute solution .
- To do this the following formula can be used: $V_1 \times C_1 = V_2 \times C_2$.

Where:

V_1 = Volume of starting solution needed to make the new solution.

C_1 = Concentration of starting solution.

V_2 = Final volume of new solution.

C_2 = Final concentration of new solution.



- **For example:**

- **Make 5ml of 0.25M solution from a 1.0M solution:**

- **Since:** $V_1 \times C_1 = V_2 \times C_2$.

Where: $V_1 = ?$, $C_1 = 1M$, $V_2 = 5ml$, $C_2 = 0.25M$

$$(V_1) \times (1M) = (5ml) \times (0.25M).$$

$$V_1 = [(5 \times 0.25)/1] = 1.25 \text{ ml}$$

- **So,** 1.25ml of the 1M solution is needed (starting solution) then complete the volume up to 5 ml.

- **Note:**

(since the diluted solution should have a final volume of 5ml)

- $(V_1 - V_2) = 5ml - 1.25ml = 3.75ml$).

- 3.75ml of diluent (generally water) should be added to the 1.25ml of starting solution.



3) Serial Dilutions:

- It is a stepwise dilution of a solution, where the **dilution factor is constant at each step.**

- The source of dilution material for each step comes from the diluted material of the previous step.

- **Example:**

- Starting with a 2.0 M stock solution of hydrochloric acid, prepare four standard solutions by serial dilution of the following Molarity respectively 1 M, 0.5 M, 0.25 M, 0.125 M. [with 1:2 dilution]

- **Dilution factor (D.F)** = final volume / aliquot volume = $2/1 = 2$

- **To prepare standard solution 1,**

- 1 ml of the stock 2.0M solution is needed and volume made up to 2 ml with distilled water (never forget to mix properly).



-To prepare standard solutions 2-4, 1 ml of the previously diluted solution is taken and volume is made up to a final volume of 2 ml by the addition of distilled water.

- How to calculate the concentration of the diluted solutions if they unknown ?

- **First:** find the D.F:

- **Dilution factor (D.F)** = final volume / aliquot volume = $2/1 = 2$

- **Second:** divide the previous solution concentration by the D.F:

- Concentration of solution 1 = $2.0 \text{ M stock solution} / 2 = 1 \text{ M}$

- Concentration of solution 2 = $1\text{M}/2 = 0.5 \text{ M}$

- Concentration of solution 3 = $0.5\text{M}/2 = 0.25 \text{ M}$

- Concentration of solution 4 = $0.25/2 = 0.125 \text{ M}$



Thank you

