# Quantitative amino acids estimation by Ninhydrin method

BCH303 [Practical]

#### **Amino acid quantitation:**

- Free amino acids, as well as amino acids released from macromolecules such as peptides, proteins or glycoproteins.
- Rapid and accurate.
- Importance?

#### Amino acid analysis:

• Refers to the methodology used to determine the amino acid composition or content of proteins, peptides, and other pharmaceutical preparations.

• Methods?

#### Ninhydrin:

- Detecting amino acids.
- Microgram amounts.
- Non-selectivity?
- Other convenient reagents are available which can react with the alpha amino group to form colored or fluorescent derivatives. These include fluorescamine, dansyl chloride, dabsyl chloride, etc., <u>used in the detection of trace amounts of amino acids at the nanogram level.</u>

## Quantitative estimation of amino acid using Ninhydrin reagent:

- Ruhemann's purple (RP) was discovered by Siegfried Ruhemann in 1910.
- Ruhemann's purple formed by the reaction at 570 nm is measured, whereas for **imino** acids, the absorption happens at 440 nm.
- The color intensity produced is proportional to the amino acid concentration.
- So ?

#### **Standard curve:**

- Direct relationship between color and concentration → direct relationship between concentration and absorbance.
- The standard curve (also called calibration curve): is a type of graph used as a quantitative research technique that shows the relationship between different known concentrations of a substance and the absorbance at a specific wave length.
- Is most commonly used to determine the concentration of a substance (unknown), using serial dilution of solutions (Standard solutions) of known concentrations.



#### **Standard curve:**




Figure 1. A standard curve showing the relation between the absorbance of different concentrations of a substance.

### **Practical part**

#### Aim:

- Determination of amino acids quantity using ninhydrin reaction.
- Understanding and contructing a standard curve.

#### **Principle:**

- At nuetral pH, ninhydrin destroys each primary α-amino acid to form Ruhemann's purple.
- Maximum absorption at about 570 nm.
- Imino acids yields a yellow- orange product at neutral pH.
- The intensity of the color resulted is <u>linearly proportional</u> to the concentration of the amino acids present in the solution

#### **Results:**

Table 1. Concentration od standard amino acid solution and their absorbance at 570 nm.

absorbance at 570 mm.		
Test tube	Amino acid concentration	Absorbance at 570 nm
	[µg/ml]	
Blank		
A		
В		
С		
D		
E		
Unknown		
sample		




Figure 1. standard curve of amino acid using ninhydrin method.

#### **Home Work:**

1. Search for a method used in amino acid analysis aimed to identify the identity or composition of individual amino acids.