BCH312 [Practical]

#### **Preparation of Different Buffer Solutions**



- All biochemical reactions occur under <u>strict conditions</u> of the concentration of hydrogen ion.
- Biological life cannot withstand large changes in hydrogen ion concentrations which we measure as the pH.
- Those solutions that have the <u>ability to resist changes</u> in pH upon the addition of **limited amounts** of acid or base are called **BURFERS**.

# **Two types of Buffers**



## **Mechanism of Action:**

**How buffers can resist the change in pH?** 

-Example using [HA/A<sup>-</sup>] buffer:
→ Where: HA is Weak acid and A- is conjugated base [its salt].

```
HA \rightleftharpoons H^+ + A^-
```

If  $\mathbf{H}^+$  (acid) is added to this buffer system  $\rightarrow \mathbf{H}^+$  will react with <u>conjugated base</u>  $\rightarrow$  to give conjugate acid.



If **OH**<sup>-</sup> (base) is added to this buffer system  $\rightarrow$  **OH**<sup>-</sup> will react with conjugated acid  $\rightarrow$  to give <u>conjugate base</u> and H<sub>2</sub>O.



#### **Mechanism of Action**



# Mechanism of Action cont':

#### **Example:**

- **Buffer system:**  $CH_3COOH / CH_3COO^-$
- □ When acid [H<sup>+</sup>] added:

 $CH_3COO^- + H^+ \longrightarrow CH_3COOH$ 

■ When base [OH<sup>-</sup>] added:

 $CH_3COOH + OH^- \longrightarrow CH_3COO^- + H_2O$ 

• NOTE: It resists pH changes when it's two components are present in specific proportions.

→ Thus the buffer is effective as long <u>as it does not run out</u> of one of its components. (There are enough conjugated base and conjugated acid to absorb the  $H^+$  ions or  $OH^-$  ions added to the system respectively).

conjugated baseconjugated acidCH3COO-CH3COOH

### Henderson-Hasselbalch equation:

#### It is often used to perform:

- 1. To calculate the pH of the Buffer.
- 2. To preparation of Buffer.

$$pH = pK_a + \log \frac{\left[A^{-}\right]}{\left[HA\right]}$$

- It relates the Ka [dissociation constant] of a weak acid, [HA] concentration of weak acid component, [A-] concentration of conjugate base [salt of the weak acid] component and the pH of the buffer.
- □ The equation is derived from the acid dissociation constant.

### Henderson-Hasselbalch equation cont':

■ A buffer is **best used close to its pKa** [to act as a good buffer the pH of the solution must be within one pH unit of the pKa].

→ The buffer capacity is optimal when the ratio of the weak acid to its salt is 1:1; that is, when pH = pKa

# **Calculating the pH:**



[Pkw : number of dissociation constant of H<sub>2</sub>O].

#### Preistieel Pert



#### □ To understand the behaviour and nature of buffers solutions.

□ To learn how to prepare buffers.

## A) Nature of buffers:

#### **Method:**

- 1. You are provided with: 0.2M solution of  $CH_3COOH$ , 0.2M solution of  $CH_3COONa$ .
- → Determine which is the weak acid and which is the conjugated base [or its salt].

2. Calculate the volume that you must take from  $CH_3COOH$  and  $CH_3COONa$  to prepare the following mixtures with final **volume of the solution =20 ml** :

- 1. 100% [HA]
- 2. 75% [HA], 25% [A<sup>-</sup>]
- 3. 50% [HA], 50% [A<sup>-</sup>]
- 4. 25% [HA], 75% [A<sup>-</sup>]

3. Calculate the pH for each solution [pKa of  $CH_3COOH = 4.76$ ].

4. Follow the table.

### **Calculations:**

1. To Calculate the volume that you must take from CH3COOH and CH3COONa to prepare the previous mixtures with final volume of the solution =20 ml:

#### (A) <u>100% [HA]:</u>

The final volume is 20ml, So: 20 x 100% = (20 x 100)/100 = 20 mlTake 20ml of HA and measure the pH.

#### (**B**) <u>75% [HA], 25% [A<sup>-</sup>]:</u>

From HA=  $20 \times 75\% = (75 \times 20) / 100 = 15$  ml

From  $A^{-} = 20 \ge 25\% = (25 \ge 20)/100 = 5 = 5$  ml

→ Mix 15ml HA and 5 ml A<sup>-</sup> and measure the pH (measured PH) note that the total volume is 20 ml [15ml + 5ml = 20ml]

□ The same way for other mixtures ...

Note: HA : as  $CH_3COOH$ . A<sup>-</sup> : as  $CH_3COONa$ .

#### **Calculations cont':**

2. To Calculate the pH for the previous mixtures with pKa of CH<sub>3</sub>COOH = 4.76 :

(A) <u>100% [HA]:</u>

$$pH = (pKa + p[HA]) \longrightarrow p[HA] = -log \ 0.2 = 0.69 \implies pH = (4.76 + 0.69) = 2.72$$

(**B**) <u>75% [HA]</u>, <u>25% [A<sup>-</sup>]</u>:

pH = pka + log [A<sup>-</sup>] / [HA] → pH= 4.76 + log [A<sup>-</sup>]/[HA] → [HA] =  $C_1 \times V_1 = C_2 \times V_2$ =0.2 × 15 =  $C_2 \times 20 = C2 = 0.15M$ 

→ 
$$[A^-] = C_1 X V_1 = C_2 X V_2$$
  
= 0.2 X 5 = C<sub>2</sub> X 20 = C2 = 0.05 M

So, pH =  $4.76 + \log 0.05 / 0.15 \rightarrow \text{pH} = \frac{4.282}{2}$ 

#### **Calculations cont':**

#### (C) <u>50%[HA], 50%[A]:</u>

pH = Pka + log [A<sup>-</sup>] / [HA] → pH= 4.76 + log [A<sup>-</sup>]/[HA] → [HA] =  $C_1 \times V_1 = C_2 \times V_2$ =0.2 × 10 = C2 × 20 = C2 = 0.1M

→ [A-] = 
$$C_1 \times V_1 = C_2 \times V_2$$
  
= 0.2 × 10 =  $C_2 \times 20 = C2 = 0.1$  M

So, pH =  $4.76 + \log 0.1/0.1 \rightarrow pH = 4.76 + 0 = 4.76$  [pH=pka]

#### (D) <u>25% [HA], 75% [A-] :</u>

pH = pka + log [A-] / [HA] → pH= 4.76 + log [A-]/[HA] → [HA] =  $C_1 X V_1 = C_2 X V_2$ =0.2 X 5 =  $C_2 X 20 = C2 = 0.05M$ 

→ 
$$[A^-] = C_1 X V_1 = C_2 X V_2$$
  
= 0.2 X 15 = C<sub>2</sub> X 20 = C2 = 0.15 M

So, pH =  $4.76 + \log 0.15 / 0.05 \rightarrow pH = 5.24$ 



Solutions	HA	$\mathbf{A}^{\text{-}}$	Final volume	Calculated	Measured	2M HCl	Measured	The difference
	CH <sub>3</sub> COOH (ml)	CH <sub>3</sub> COONa (ml)	( <b>ml</b> )	рН	рН	( <b>ml</b> )	рН	
100%[HA]	20 ml	0	20 ml	2.729		0.1		
75%[HA],25%[A <sup>-</sup> ]	15 ml	5 ml	20 ml	4.28		0.1		
50%[HA],50%[A <sup>-</sup> ]	10 ml	10 ml	20 ml	4.76		0.1		
25%[HA],75%[A <sup>-</sup> ]	5 ml	15 ml	20 ml	5.24		0.1		

### **B)** Preparation of buffer:

You are provided with 0.2M acetic acid and solid sodium acetate.
 Prepare 50ml of a 0.19M acetate buffer pH =4.86 if you know that (pKa=4.7).

#### **Calculations:**



# **Calculations cont' (first method):**

**To prepare a buffer Henderson-Hasselbalch equation is used:** 

**pH = pka +log [A-] / [HA]** 

First calculate the concentration of the weak acid and its conjugated base that make up the buffer with 0.19M: Assume  $[A^-] = y$  and [HA] = 0.19 - y➔ So:  $4.86 = 4.76 + \log \frac{y}{0.19 - y}$  $0.1 = \log \frac{y}{0.19-y}$   $\rightarrow$  by taking the "Anti log for both sides"  $\rightarrow$   $1.26 = \frac{y}{0.19-y}$ y = 1.26 x (0.19-y)  $\rightarrow$  y = 0.24 - 1.26 y  $\rightarrow$  y + 1.26 y = 0.24  $\rightarrow$  2.26 y = 0.24y=0.11 M [which is the concentration of [A<sup>-</sup>] in the buffer ] So, 0.11 + 0.08 = 0.19M[HA] = 0.19 - 0.11= 0.08 M [which is the concentration of [HA] in the buffer ]

## **Calculations cont' (second method):**

**To prepare a buffer Henderson-Hasselbalch equation is used:** 

**pH = pka +log [A-] / [HA]** 

First calculate the concentration of the weak acid and its conjugated base that make up the buffer with 0.19M:

First: 
$$4.86 = 4.76 + \log \frac{[A-]}{[HA]}$$
  
 $0.1 = \log \frac{[A-]}{[HA]} \Rightarrow by taking the "Anti log for both sides" \Rightarrow \frac{[A-]}{[HA]} = 1.26 = \frac{1.26}{1}$   
SO:  $\frac{1.26}{2.26}$  of total = [A-] and  $\frac{1}{2.26}$  of total = [HA]  
[A-] =  $\frac{1.26}{2.26} \ge 0.11$  M [which is the concentration of [A-] in the buffer ]  
[HA-] =  $\frac{1}{2.26} \ge 0.19 = 0.08$  M [which is the concentration of [HA] in the buffer ]  
 $0.11 + 0.08 = 0.19$ M

### **Calculations cont':**

• To calculate the volume needed from [HA] to prepare the buffer, No. of mole of [HA] should be calculated first:

No. of mole = Molarity x Volume of solution in L =  $0.08 \times 0.05 = 0.004$ mole So, M of stock = no. of mole / Volume in Liter 0.2 = 0.004 / V

→ V = 0.02 L = 20 ml

**To calculate the weight needed from [A-] to prepare the buffer, No. of mole of [A-] should be calculated first:** 

**No. of mole = Molarity x volume of solution in L** 

= 0.11 X 0.05 = 0.0055 mole

weight in (g) of [A-] = No. of moles x MW

→ =0.0055 x 82 = 0.451 g

Problem 1-27, p39 Problem 1-28, p40



■ Now take 20 ml from 0.2M acetic acid and 0.451 g from solid sodium acetate and then complete the volume up to 50 ml by addition of water.

# **C) Testing for buffering behaviour:**

- □ In one beaker add 10ml of 0.19M acetate buffer that you have prepared, and in another beaker add 10ml of 0.2M KCl.
- □ Measure the pH.
- □ Add 0.1ml from 2M HCl to for both solutions.
- □ Measure the pH after the addition.

Solution	Measured pH	Add 2M HCl	Measured pH
0.19M acetate buffer		0.1 ml	
o.2M KCl		0.1 ml	