



BCH 312
Experiment (3)

Titration of a
weak acid with
strong base



+ Objectives



- 1) To study titration curves.
- 2) To determine the pK_a value of a weak acid.
- 3) To reinforce the understanding of buffers.

+ Weak Acid

- Weak acids or bases do not dissociate completely, therefore an equilibrium expression with K_a must be used.
- The K_a is a **quantitative measure of the strength of an acid in solution.**
- K_a is usually expressed as pK_a .
 - **$pK_a = -\log K_a$**
- As an acid gets weaker, its K_a gets smaller and pK_a gets larger and vice versa, ***for example:***
 - HCl is a strong acid , it has a large value of K_a (low pK_a)
 - CH_3COOH is a weak acid , it has low value of K_a (high pK_a)
- pK_a values of weak acids can be determined mathematically or practically by the use of titration curves.



Titration curves

- There are many uses of titration, one of them is to indicate the pK_a value of the weak acid by using the titration curve.
- Titration Curves are produced by monitoring the pH of a given volume of a sample solution after successive addition of acid or alkali.
- The curves are usually **plots of pH against the volume of titrant added.**

+ How to calculate the pH at each point of the titration curve

■ [1] Before any addition of strong base the (starting point):


- The weak acid is in the full protonation form $[\text{CH}_3\text{COOH}]$ (electron donor).
- In this point pH of weak acid $< \text{pK}_a$
- We can calculate the pH from:

$$\text{pH} = (\text{pK}_a + P[\text{HA}]) / 2$$

■ [2] When certain amount of strong base added (any point before the middle of titration):

- The weak acid is starting to dissociate, $[\text{CH}_3\text{COOH}] > [\text{CH}_3\text{COO}^-]$ (Donor > Acceptor).
- In this point the pH of weak acid $< \text{pK}_a$
- We can calculate the pH from:

$$\text{pH} = \text{pK}_a + \log[\text{A}^-] / [\text{HA}]$$



+Continue

■ [3] At middle of titration:

- $[\text{CH}_3\text{COOH}] = [\text{CH}_3\text{COO}^-]$ (Donor=Acceptor) ,
- $\text{pH} = \text{pK}_a$.
- The component of weak acid work as a **Buffer** (A solution that can resistant the change of PH) .
- Buffer capacity= $\text{pK}_a \pm 1$


Note: pK_a is defined as the pH value at middle of titration at which they will be $[\text{donor}] = [\text{acceptor}]$.

■ [4] any point after mid of titration:

- $[\text{CH}_3\text{COOH}] < [\text{CH}_3\text{COO}^-]$, (Donor < Acceptor) .
- In this point the $\text{pH} > \text{pK}_a$
- We can calculate the pH from:
 - $\text{pH} = \text{pK}_a + \log \frac{[\text{A}^-]}{[\text{HA}]}$

■ [5] At the end point

- The weak acid is fully dissociated $[\text{CH}_3\text{COO}^-]$ (electron acceptor) .
- pH of weak acid $> \text{pK}_a$
- Approximatly, all the solution contains CH_3COONa so we first must calculate pOH , then the pH:
 - $\text{pOH} = (\text{pK}_b + \text{p}[\text{A}^-]) / 2$
 - $\text{pH} = \text{pK}_w - \text{pOH}$



+ Notes:

The pH calculated by different way :

- [1] At starting point $\text{pH} = (\text{pK}_a + \text{p}[\text{HA}]) / 2$
- [2] At any point within the curve (after , in or after middle titration)
 $\text{pH} = \text{pK}_a + \log[\text{A}^-]/[\text{HA}]$ (Henderson-Hasselbalch equation)
- [3] At end point $\text{pOH} = (\text{pK}_b + \text{p}[\text{A}^-]) / 2$
 $\text{pH} = \text{pK}_w - \text{pOH}$

Henderson-Hasselbalch equation is an equation that is often used to :

- To calculate the pH of the buffer
- Used in buffer preparation.
- To calculate the pH in any point within the **titration curve** (Except starting and ending point)

+ Example

- Determine the pH value of 500 ml weak acid (0.1M) , titrated with 0.1M KOH (Pka=5) ?

After addition 100 ml , 250 ml , 375 and 500 ml of KOH??

[1] pH after the addition of 100 ml of KOH?

- $\text{pH} = \text{pK}_a + \log \frac{A^-}{HA}$



- No. of moles of HA remaining = No. of moles of HA originally – No. of moles of KOH
- No. of moles of KOH = $M \times V \text{ (in L)} = 0.1 \times 0.1 = \underline{0.01 \text{ moles}} = \text{No. of moles of } A^-$
- No. of moles of HA originally = $M \times V \text{ (in L)} = 0.1 \times 0.5 = 0.05 \text{ moles}$
- No. of moles of HA remaining = $0.05 - 0.01 = \underline{0.04 \text{ moles}}$
- $\text{pH} = 5 + \log 0.01 / 0.04$
- $\text{pH} = 4.4$



[2] pH after the addition of 250 ml of KOH??

- No. of moles of HA remaining = No. of moles of HA originally – No. of moles of KOH
- No. of moles of KOH = $0.1 \times 0.25 = \underline{0.025 \text{ moles}}$ = **No. of moles of A⁻**
- No. of moles of HA originally = $0.1 \times 0.5 = 0.05 \text{ moles}$
- No. of moles of HA remaining = $0.05 - 0.025 = \underline{0.025 \text{ mole}}$
- $\text{pH} = 5 + \log 0.025/0.025$
- **pH = 5 = Pka** (at mid point , The component of weak acid work as a Buffer , has a buffering capacity 5 ± 1)



[3] pH after the addition of 375 ml of KOH??

- No. of moles of HA remaining = No. of moles of HA originally – No. of moles of KOH
- No. of moles of KOH = $0.1 \times 0.375 = \underline{0.0375 \text{ moles}}$ = **No. of moles of A⁻**
- No. of moles of HA originally = $0.1 \times 0.5 = 0.05 \text{ moles}$
- No. of moles of HA remaining = $0.05 - 0.0375 = \underline{0.0125 \text{ mole}}$
- $\text{pH} = 5 + \log 0.0375/0.0125$
- $\text{pH} = 5.48$



✚ [4] pH after the addition of 500 ml of KOH??

Note: 500 ml the same volume of weak acid that mean the all weak acid as $[CH_3COO^-]$

- $pOH = (pK_b + p[A^-]) / 2$ $pK_b = pK_w - pK_a \rightarrow pK_b = 14 - 5 = 9$
- $p[A^-] = -\log [A^-] \rightarrow [A^-] = ??$
- No of moles of KOH = $0.1 \times 0.5 = 0.05$ moles = **No. of moles of A^-**
- Total volume = $500 + 500 = 1000 = 1L$
- $[A^-] = 0.05/1 = 0.05 M$
- $p[A^-] = -\log 0.05 = 1.3$
- $pOH = (9 + 1.3) / 2 = 5.15$
- $pH = pK_w - pOH$
- $pH = 14 - 5.15 = 8.85$ (at end point)



Method

- You are provided with 10 ml of a 0.1M CH_3COOH weak acid solution, titrate it with 0.1m NaOH adding the base drop wise mixing, and recording the pH after each 0.5 ml NaOH added until you reach a pH=10.

Volume of 0.1 NaOH (ml)	pH
0	
0.5	
1	
1.5	



Discussion:

- 1- Plot a Curve of pH versus ml of NaOH added, determine the pK_a value from the curve.
- 2- Calculate the pH of the weak acid HA solution after the addition of 3ml, 5ml, and 10ml of NaOH.
- 3- Compare your calculated pH values with those obtained from Curve.
- 4- At what pH-range did the acid show buffering behavior? What are the chemical species at that region, what are their proportions?