

**THE EFFECT OF  
INCUBATION TIME ON  
THE RATE OF AN ENZYME  
CATALYZED REACTION**


The progress of an enzyme catalyzed reaction may be followed by:

Measuring either the quantity of substrate used up

OR

the quantity of product formed and plotting against time.


Typically, a curve of the following type is obtained



**The rate of the reaction is highest at time zero and decreases with increasing time, eventually falling to zero itself, where the above curve reaches a plateau.**

This usually occurs either **when all the substrate is used up or when equilibrium is reached.**

- The initial rate of reaction  $v_i$ , measured as the **tangent to the above curve at the origin (time= 0)**, is used in the study of enzyme kinetics and is affected by factors which include enzyme concentration, substrate concentration, temperature and pH.



# Principle

**Under acid conditions**, the enzyme catalyzes the hydrolysis of p-nitrophenyl phosphate (pNPP) to inorganic phosphate and p-nitrophenol.

**If base is added** to the mixture after the completion of the reaction, the p-nitrophenol is converted to a colored form which absorbs lights at 405 nm.

- Assuming a path length of 1.0 cm and an extinction coefficient of  $18.8 \times 10^3$  liter mol<sup>-1</sup> cm<sup>-1</sup>, you can use the absorbance at 405 nm to calculate the number of micromoles of p-nitrophenol released. Since this is a fixed-time assay that is stopped after 5 minutes, the velocity of the reaction ( $\mu\text{moles of p-nitrophenol/minute}$ ) can be easily computed.

**Set up a series of identical enzyme reaction tubes each of which is allowed to incubate for a different period of time (zero through 30 minutes).**

**The results should indicate how long the reaction is linear under the given conditions of substrate and enzyme concentration.**

- 1-To demonstrate the validity of the assumption that the relationship between product yield and time has been linear throughout.
- 2- To Calculate the initial velocity of reaction,  $v_i$ .

## **Objectives**



# Materials

## Chemicals:

- Sodium acetate buffer
- Magnesium chloride
- p-nitrophenyl phosphate
- Potassium hydroxide
- Stock solution of crude wheat germ Acid Phosphatase

## Equipments:

- Water bath
- Spectrophotometer
- Test tubes Pipettes

## Method

- Prepare a series of seven reaction tubes labeled 0 through 30 minutes at 5-minute intervals (0, 5, 10 ... minutes).
- To each of these tubes add:  
0.5 ml of sodium acetate buffer (pH 5.7),  
0.5 ml of  $\text{MgCl}_2$ ,  
0.5 ml of p-nitrophenyl phosphate, (**Except blank**)  
5 ml of distilled water.
- Place all the tubes in a test rack situated in a water bath at 37 °C and let the temperature equilibrate for 5 minutes.

Add 0.5ml of the enzyme to each tube, mix, start the stopwatch, and let the reaction proceed for 5 minutes before adding the KOH to terminate the reaction.

Total incubation time (min)	Clock time (min)	
	Start reaction (add enzyme)	Stop reaction (add KOH)
0	0	0
5	0	5
10	2	12
15	4	19
20	6	26
25	8	33
30	10	40

## Results

Incubation time(min)	Absorbance at 405nm
0 (blank)	
5	
10	
15	
20	
25	
30	

- $\text{Velocity} = A \times 10^6 / \epsilon \times \text{Time}$

Use the extinction coefficient ( $18.8 \times 10^3 \text{ liter mol}^{-1} \text{ cm}^{-1}$ ) for p-nitrophenol to calculate the micromoles of product released at each time point. Prepare a graph, plotting  $\mu\text{moles}$  of p-nitrophenol released against time.

THANK YOU