

Estimation of Serum Creatinine, Urine Creatinine and Creatinine Clearance

-Kidney functions:

- The kidneys serve three essential functions:

1. They function as filters, removing metabolic products and toxins from the blood and excreting them through the urine.
2. They regulate the body's fluid status, electrolyte balance, and acid-base balance.
3. The kidneys produce or activate hormones that are involved in erythropoiesis, Ca^{2+} metabolism, and the regulation of blood pressure and blood flow.



-Renal function tests (RFT):

- Are used to detect the presence of **renal diseases** and **assess their progress**.
- The most widely used test is to measure the glomerular filtration rate (GFR), that is, the rate of filtrate formation by the kidneys.

- Glomerular Filtration Rate:

- Under normal conditions, approximately 625 mL of plasma flow through the kidneys **each minute** and the volume of plasma filtered is **125 ml/ min** which is called the glomerular filtration rate.
- **Glomerular filtration rate (GFR):** is the volume of plasma filtered by the kidneys per unit of time.
- **GFR** is an important and the best overall measurement in the evaluation of kidney function.

-Measuring the GFR:


- Accurate measurement of the GFR by clearance tests requires determination of the concentration, in **plasma and urine**, of a substance is known to be completely filtered from the plasma at the glomerulus.
- This substance **must not be reabsorbed nor secreted by renal tubules, broken down, or accumulated by the tubules** and must remain at a **constant concentration** in the plasma throughout the period of urine collection.
- **It's clearance is given by: $\text{clearance} = \frac{U \cdot V}{P}$, where:**
 - ➔ U= concentration of any substance in urine.
 - P= concentration of the same substance in plasma .
 - V= volume of urine(ml/min).

Substances clearance used for Measuring GFR

	Inulin Clearance	Creatinine Clearance	Urea Clearance
Source	Non-toxic fructose polymer.	End-product of skeletal muscle creatine metabolism.	End product of protein metabolism.
Advantages	Not reabsorbed or secreted.	An <u>endogenous</u> product of muscle metabolism; near constant production.	An <u>endogenous</u> product of protein.
Disadvantages	Not made by body; must be injected (exogenous).	Small amount is secreted .	-Partially reabsorbed. -synthesis varies with diet.

- Creatinine:

- Creatinine is derived from “**creatine**” which is synthesized in **the liver, kidney and pancreas** it moves through the circulation and is taken up entirely by muscles.
- Creatinine is a substance that, in health, is easily excreted by the kidney.
- It is the byproduct of muscle energy metabolism and is produced at a constant rate according to the muscle mass of the individual.
- Endogenous creatinine production is **constant** as long as the **muscle mass remains constant**.



- Serum Creatinine:

- High plasma creatinine:

- Plasma creatinine tends to be higher in subjects with a **large muscle mass**.

- Other **non-renal** causes of increased plasma creatinine include the following :

1. A high meat intake can cause a **temporary increase**.
2. Transient, small increases may occur **after vigorous exercise** .

- Urine Creatinine:

- **Decreased urine creatinine is found in:**

- Advanced renal disease.
- Renal stenosis.
- Hyperthyroidism.

Increased urine creatinine is found in:

- Diabetes mellitus.
- Hypothyroidism.

- Creatinine clearance:

- A measure of the amount of creatinine eliminated (filtered) from the blood by the kidneys.
- Creatinine is cleared from the body fluids **almost entirely** by glomerular filtration (small amount is secreted by kidney tubules).
- Therefore, the clearance of creatinine can be used to **assess GFR.**
- Because measurement of creatinine clearance **does not** require intravenous infusion into the patient, **this method is much more widely used than inulin clearance** for estimating GFR clinically.

- Clinical Implications:

1. Decreased creatinine clearance is found in any condition that decreases renal blood flow:

- a. Impaired kidney function.
- b. Shock, dehydration.
- c. Hemorrhage.

2. Increased creatinine clearance is found in:

- a. Pregnancy.

-Reference Values:

- Urine creatinine: 1- 2 g/ 24h
- (Serum) creatinine: 0.6–1.2 mg/dL
- Normal creatinine clearance= 100-130 ml/min/1.73m²

- Note:

- What 1.73 m² means?

- Kidney function is proportional to kidney size, which is proportional to body surface area.
A of 1.73 m² is the normal mean value for young adults.
- Adjustment for body surface area is necessary when comparing a patient's estimated **GFR** to normal values or to the levels defining the stages of **Chronic kidney disease (CKD)**.

Chart 2 - Chronic kidney disease staging

Stage	Description	GF (ml/min/1.73m ²)
I	Kidney lesion with normal or increased GF	≥ 90
II	Kidney lesion with mild GF decrease	60-89
III	Kidney lesion with moderate GF decrease	30-59
IV	Kidney lesion with marked GF decrease	15-29
V	Functional kidney failure or undergoing SRT	< 15

SRT- substitutive renal therapy. Source: National Kidney Foundation, 2002.

Practical Part

- Objective:

- 1- To estimate creatinine in serum and urine.
- 2- To calculate creatinine clearance value.

- Principle:

(Jaffe's method):

- Colorimetric estimation of creatinine using the alkaline picrate method:

Creatinine + picric acid \longrightarrow **Creatinine picrate (orange)**

- Absorbance at **520 nm**.

-Method:

1- Set up a series of 8 test tube as follows:

Chemical	Standard (3mg/dl) (serum)		Test (serum) (C)		Standard (0.75mg/dl) (Urine)		Test (urine) (G)		Blank
	(A)	(B)	(C)	(D)	(E)	(F)	(G)	(H)	
Water	1.5 ml	1.5 ml	1.5 ml	1.5 ml	1.5 ml	1.5 ml	1.5 ml	1.5 ml	2 ml
Standard (serum)	0.5 ml	0.5 ml	-	-	-	-	-	-	-
Serum Sample	-	-	0.5 ml	0.5 ml	-	-	-	-	-
Standard (Urine)	-	-	-	-	0.5 ml	0.5 ml	-	-	-
Urine Sample	-	-	-	-	-	-	0.5 ml	0.5 ml	-
Picric acid	6 ml	6 ml	6 ml	6 ml	6 ml	6 ml	6 ml	6 ml	6 ml

2-Immerse the Tubes carefully in the boiling water bath **for 40 seconds.**

4- Pipette 0.6 ml of NaOH to all tube.

5- Let the tubes stand **for 20 min.**

6- Read the absorbance at **520 nm.**

-Results:

Tube	Standard (serum)		Test (serum)		Standard (urine)		Test(Urine)	
	(A)	(B)	(C)	(D)	(A)	(B)	(C)	(D)
Absorbance at 520 nm								
Average(Mean of Absorbance)								

- Calculation:

Patient information: 24h urine volume = 100ml, gender: women, body surface: 1.6m², DF=100.

1-Serum creatinine =

(Mean Absorbance of serum test ÷ Mean Absorbance of Standard) X concentration of Serum standard = mg / dl

2-Urine creatinine =

(Mean Absorbance of Urine test ÷ Mean Absorbance of Standard) X concentration of Urine standard X DF= mg / dl

(To compare with normal range, convert from mg/dl to g/24 h)

3- Creatinine Clearance :

$$=U.V/ P$$

$$= [(Urinary creatinine (mg/dl)) / (plasma creatinine (mg/dl))] \times Urine volume(ml/min) = B$$

B -----> 1.6 m² (person surface area)

? -----> 1.73 m²

-Corrected for surface area= ml/min/1.73 m²

-Example:

- Find the Creatinine Clearance if you know that the Urine creatinine $U = 488 \text{ mg/dl}$, Serum creatinine $P = 2.32 \text{ mg/dl}$, Volume of urine in 24 h $= 100 \text{ ml}$ and A "surface area" $= 1.6 \text{ m}^2$?

→ Creatinine Clearance: $= U.V / P$

$$= (488 \text{ mg/dl} \div 2.32 \text{ mg/dl}) \times (100 \div 1440^*) = 14.6 \text{ ml/min}$$

14.6 ml/ min in 1.6 m^2 , find the creatinine clearance for 1.73 m^2 surface area :

$$= (14.6 \times 1.73) \div 1.6 = 15.8 \text{ ml /min/}1.73\text{m}^2$$

OR

→ Creatinine Clearance: $= (U \times V \times 1.73) / (P \times 1440 \times A)$

$$= (488 \text{ mg/dl} \times 100 \times 1.73) / (2.32 \times 1440 \times 1.6)$$

$$= 15.8 \text{ ml/min /}1.73\text{m}^2$$

* To convert 24 hour to min ($24 \times 60 = 1440$)

-Discussion:

- Comment on the concentration of creatinine in serum.
- Comment on the concentration of creatinine in urine.
- Comment on the value of **Creatinine Clearance** .