



# Preparation of Buffer Solutions by Different Laboratory Ways

# Dissociation of triprotic acid:


- **Triprotic acid** is acid that contain three hydrogens ions.
- It dissociates in solution in three steps, with three  $K_a$  values.
- **phosphoric acid** is an example of triprotic acid .
- It dissociates in solution as following:



# Preparation of buffer by several ways:

- For example if you was asked to prepare sodium phosphate buffer [  $\text{NaH}_2\text{PO}_4$  /  $\text{Na}_2\text{HPO}_4$  ]: you can prepare it by.....
  1. By mixing  $\text{NaH}_2\text{PO}_4$  (conjugate acid ) and  $\text{Na}_2\text{HPO}_4$  (conjugate base) in the proper proportions.
  2. By starting with  $\text{H}_3\text{PO}_4$  and converting it to  $\text{NaH}_2\text{PO}_4$  plus  $\text{Na}_2\text{HPO}_4$  by adding the proper amount of **NaOH**.
  3. By starting with  $\text{NaH}_2\text{PO}_4$  and converting a portion of it to  $\text{Na}_2\text{HPO}_4$  by adding **NaOH**.
  4. By starting with  $\text{Na}_2\text{HPO}_4$  and converting a portion of it to  $\text{NaH}_2\text{PO}_4$  by adding a strong acid such as **HCL**.
  5. By starting with  $\text{Na}_3\text{PO}_4$  and converting it to  $\text{Na}_2\text{HPO}_4$  plus  $\text{NaH}_2\text{PO}_4$  by adding **HCL**.
  6. By mixing  $\text{Na}_3\text{PO}_4$  and  $\text{NaH}_2\text{PO}_4$  in the proper proportions.

**HCl**  
'donate H<sup>+</sup>'



**NaOH**  
'accept H<sup>+</sup>'

**Example:** Prepare 0.1 liters of 0.045 M sodium phosphate buffer, pH=7.5, [pKa1= 2.12, pKa2 = 7.21 and pKa3 = 12.30]:

a) From concentrated (15M) H<sub>3</sub>PO<sub>4</sub> and solution of 1.5 M NaOH .


b) From solid NaH<sub>2</sub>PO<sub>4</sub> and solid NaOH.

Regardless of which method is used , the first step involves determine the buffer ionic species, calculating number of moles and amounts of the two ionic species in the buffer.


**Calculations:**

**1<sup>st</sup> →** Write the equations of phosphoric acid dissociation and the pKa of corresponding ones:

Because phosphoric acid [H<sub>3</sub>PO<sub>4</sub>] is **triprotic acid** it has 3 dissociation phases so:



**2<sup>nd</sup> →** Choose the pKa value which is near the pH value of the required buffer, to be able to know the ionic species involved in your buffer:



→ The pH of the required buffer [pH=7.5] is near the value of pKa<sub>2</sub> , consequently , the two major ionic species present are H<sub>2</sub>PO<sub>4</sub><sup>-</sup> ( conjugate acid ) and HPO<sub>4</sub><sup>2-</sup> ( conjugate base), with the HPO<sub>4</sub><sup>2-</sup> **predominating** {since the pH of the buffer is slightly basic}.

## Calculations cont':

**3<sup>rd</sup>** → calculate No. of moles for the two ionic species in the buffer:

$$\text{pH} = \text{pKa}_2 + \log \left[ \frac{\text{HPO}_4^{2-}}{\text{H}_2\text{PO}_4^-} \right] \quad \rightarrow \text{Note that: } [\text{A}^-] = \text{HPO}_4^{2-}, [\text{HA}] = \text{H}_2\text{PO}_4^-$$

- Since the buffer concentration is **0.045M**, so assume  $[\text{A}^-] = y$ ,  $[\text{HA}] = 0.045 - y$ :

$$7.5 = 7.2 + \log ( y / 0.045 - y )$$

$$7.5 - 7.2 = \log ( y / 0.045 - y )$$

$$0.3 = \log( y / 0.045 - y ) \rightarrow \text{antilog for both sides}$$

$$\rightarrow 2 = ( y / 0.045 - y ) \rightarrow y = 0.09 - 2y \rightarrow 3y = 0.09 \rightarrow y = 0.09/3 = \underline{\underline{0.03\text{M}}} \rightarrow \text{conc. of } [\text{HPO}_4^{2-}] = [\text{A}^-] = y$$

$$\text{So, conc. of } [\text{H}_2\text{PO}_4^-] = [\text{HA}] = 0.045 - y = 0.045 - 0.03 = \underline{\underline{0.015\text{M}}}$$

- **Now find the number of mole for the two ionic species in the buffer:**

$$\text{- No. of moles of } \text{HPO}_4^{2-} (\text{A}^-) = M \times V = 0.03 \times 0.1 = \boxed{0.003 \text{ moles.}}$$

$$\text{- No. of moles of } \text{H}_2\text{PO}_4^- (\text{HA}) = M \times V = 0.015 \times 0.1 = \boxed{0.0015 \text{ moles.}}$$

→ Note that Total no. of moles of phosphate buffer  
=  $M \times V = 0.045 \times 0.1 = 0.0045$  moles.

Now, to prepare the required buffer:

a) From concentrated (15M)  $\text{H}_3\text{PO}_4$  and solution of 1.5 M NaOH .

Calculations:

Start with **0.0045 mole** of  $\text{H}_3\text{PO}_4$  add **0.0045 moles** of NaOH to convert  $\text{H}_3\text{PO}_4$  completely to  $\text{H}_2\text{PO}_4^-$  (HA) , then add **0.003 moles** of NaOH to convert  $\text{H}_2\text{PO}_4^-$  to give  $\text{HPO}_4^{2-}$  (A<sup>-</sup>):

No. of moles needed of NaOH=  $0.0045+0.003=$  0.0075 moles



→ Volume of NaOH needed= no.of moles / M =  $0.0075/ 1.5 = 0.005 \text{ L} =$  5 ml

→ Volume of  $\text{H}_3\text{PO}_4$  needed =no.of moles / M =  $0.0045/ 15 =0.0003 \text{ L} =$  0.3 ml

So:

Add **5ml** of NaOH to the **0.3 ml** of concentrate  $\text{H}_3\text{PO}_4$ , mix ; then add sufficient water to bring the final volume to 0.1 liters (100 ml), and check the pH.

Remember that the two ionic species involved in the buffer are:



**b) From solid  $\text{NaH}_2\text{PO}_4$  and solid  $\text{NaOH}$ .**

Remember that the two ionic species involved in the buffer are:




**Calculations:**

Start with **0.0045 mole** of  $\text{NaH}_2\text{PO}_4$  (HA) and add **0.003 moles** of  $\text{NaOH}$  to convert  $\text{NaH}_2\text{PO}_4$  to give  $\text{Na}_2\text{HPO}_4$  ( $\text{A}^-$ ):

→ Weight in grams of  $\text{NaH}_2\text{PO}_4$  needed = no. of moles x MW =  $0.0045 \times 119.98 = \underline{0.54 \text{ g}}$

→ Weight in grams of  $\text{NaOH}$  needed = no. of moles x MW =  $0.003 \times 40 = \underline{0.12 \text{ g}}$

So: Dissolve the **0.548g** of  $\text{NaH}_2\text{PO}_4$  and **0.12g** of  $\text{NaOH}$  in some water, mix ; then add sufficient water to bring the final volume to 0.1 liters (100 ml), and check the pH.





# Practical Part

# Objective:

---

- To learn how to prepare a buffer by different laboratory ways.

# Method:

□ Prepare 0.1 liters of 0.045 M sodium phosphate buffer, pH=7.5, [pKa1= 2.12, pKa2 = 7.21 and pKa3 = 12.30]:

**a) From concentrated (15M)  $\text{H}_3\text{PO}_4$  and solution of 1.5 M NaOH :**

Add **5ml** of NaOH to the **0.3 ml** of concentrate  $\text{H}_3\text{PO}_4$ , mix ; then add sufficient water to bring the final volume to 0.1 liters (100 ml), and check the pH.

**b) From solid  $\text{NaH}_2\text{PO}_4$  and solid NaOH :**

Dissolve the **0.584g** of  $\text{NaH}_2\text{PO}_4$  and **0.12g** of NaOH in some water, mix ; then add sufficient water to bring the final volume to 0.1 liters (100 ml), and check the pH.

# Homework:

□ Prepare 0.1 liters of 0.045 M sodium phosphate buffer, pH=7.5, [pka1= 2.12, pka2 = 7.21 and pka3 = 12.30]:

c) You are provided with solid  $\text{Na}_2\text{HPO}_4$  and 2M HCl.

d) You are provided with solid  $\text{Na}_3\text{PO}_4$  and 2 M HCL .

Now, to prepare the required buffer:

c) You are provided with solid  $\text{Na}_2\text{HPO}_4$  and 2M HCl?

Calculations:

Start with **0.0045 mole** of  $\text{Na}_2\text{HPO}_4$  add **0.0015 moles** of HCl to titrate  $\text{Na}_2\text{HPO}_4$  to give  $\text{NaH}_2\text{PO}_4$  (HA):

No. of moles needed of HCl= **0.0015 moles**

→ Volume of HCl needed= no.of moles / M = 0.0015/ 2 = 0.00075 L = **0.75 ml**

→ Weight of  $\text{Na}_3\text{PO}_4$  needed = no.of moles x MW = 0.0045 x 380.12 = **1.71 g**

So: Dissolve 1.71 g of  $\text{Na}_3\text{PO}_4$  in some water, mix ; then add 0.75 ml of HCl. Finally, add sufficient water to bring the final volume to 0.1 liters (100 ml), and check the pH.

d) You are provided with solid  $\text{Na}_3\text{PO}_4$  and 2 M solution of HCL .

Calculations:

Start with **0.0045 mole** of  $\text{Na}_3\text{PO}_4$  and add **0.0045 moles** of HCl to convert  $\text{Na}_3\text{PO}_4$  completely to give  $\text{Na}_2\text{HPO}_4$  (A-), then add **0.0015 moles** of HCl to convert  $\text{Na}_2\text{HPO}_4$  to give  $\text{NaH}_2\text{PO}_4$  (HA):

No. of moles needed of HCl= 0.0045+0.0015= **0.006 moles**

→ Volume of HCl needed= no.of moles / M = 0.006/ 2= 0.003 L = **3 ml**

→ Weight of  $\text{Na}_3\text{PO}_4$  needed = no.of moles x mwt = 0.0045 x 380.12 = **1.71 g**

So: Dissolve 1.71 g of  $\text{Na}_3\text{PO}_4$  in some water, mix ; then add 3 ml of HCl. Finally, add sufficient water to bring the final volume to 0.1 liters (100 ml), and check the pH .

Remember that the two ionic species involved in the buffer are:

