

# Buffer Capacity

# Buffers:

- Buffer solutions are solutions that can resist changes in pH upon addition of small amounts of acid/base.
- Common buffer mixtures contain two substances: conjugate acid and a conjugate base .
- Together the two species (conjugate acid and conjugate base) resist large changes in pH by absorbing the  $H^+$  ions or  $OH^-$  ions added to the system.

# How buffers resist the change in pH:

1. When **H<sup>+</sup> ions** are added to the buffer system they will react with the **conjugate base** in the buffer as following:



2. When **OH<sup>-</sup> ions** are added they will react with the **conjugate acid** in the buffer as following:



➔ **Thus the buffer is effective as long as it does not run out of one of its components.**

( There are enough conjugated base and conjugated acid to absorb the H<sup>+</sup> ions or OH<sup>-</sup> ions added to the system respectively).

# Buffer Capacity:

- Quantitative measure of this resistance to pH changes is called **buffer capacity**.
- Buffer capacity can be defined in many ways, **it can be defined as:**
  - ➔ The number of moles of  $\text{H}^+/\text{OH}^-$  ions that must be added to one liter of the buffer in order to decrease /increase the pH by one unit respectively.
- The buffer capacity is **expressed as  $\beta$**  and can be derived from Henderson Hasselbalch equation:

$$\beta = \frac{2.3 K_a [\text{H}^+][\text{C}]}{(K_a + [\text{H}^+])^2}$$

From the equation ➔ the buffer capacity is **directly proportional** to the buffer concentration.

- **Where :**

$\beta$  = the buffer capacity ,  $[\text{H}^+]$  = the hydrogen ion concentration of the buffer ,  $[\text{C}]$  = concentration of the buffer and  $K_a$  = acid dissociation constant

# Practical Part

# Objective:

---


- To understand the concept of buffer capacity.
- To determine the maximum buffer capacity of a number of buffer solutions.
- To establish the relationship between buffer capacity and buffer concentration.

# Method:

- You are provided with two acetate buffer (pH=5) ; 0.1 M acetate buffer and 0.2 M acetate buffer .
- In one beaker add 8ml of the 0.1 M acetate buffer buffer, and in another beaker add 8ml of 0.2 M acetate buffer.
- Start the titration by adding 0.5 ml of 0.1 M HCl from the burette and determine the pH of the solution after each addition.
- Continue adding acid in until pH falls to about 2 pH units from your starting pH .
- Plot a Curve of pH against ml of HCl added.
- Calculate the buffer capacity (which one has higher buffer capacity. why?)

What do you conclude finally about the relationship between,  
Buffer concentration and buffer capacity ?

# Results:





What do you conclude finally about the relationship between,  
Buffer concentration and buffer capacity ?

# Results:

