BCH312 [Practical]

### **Buffer Capacity**



• Quantitative measure of buffer resistance to pH changes is called **<u>buffer capacity.</u>** 

Buffer capacity can be defined in many ways, it can be defined as:

The number of moles of  $H^+/OH^-$  ions that must be added to <u>one liter</u> of the buffer in order to decrease /increase the pH by <u>one unit</u> respectively.

The instantaneous buffer capacity is **expressed as**  $\beta$  and can be derived from Henderson Hasselbalch equation:

#### □ Where :

 $\beta$  = the buffer capacity , [H+] = the hydrogen ion concentration of the buffer , [C] = concentration of the buffer and Ka= acid dissociation constant.

# **Practical buffer capacity:**

Buffer capacity of acid and alkaline direction:

→ Buffer capacity  $_{a}(BC_{a}) =$  the number of moles of H<sup>+</sup> that must be added to one liter of the buffer in order to decrease the pH by one unite.

This called buffer capacity in the ACID direction.

$$BC_{a} = \frac{9[HA] [A^{-}]}{10 [HA] + [A^{-}]}$$


→ Buffer capacity  $_{b}(BC_{b}) =$  the number of moles of OH<sup>-</sup> that must be added to one liter of the buffer in order to <u>increase the pH by one unite</u>.

This called buffer capacity in the ALKAILNE direction.


$$BC_{b} = \frac{9[HA] [A^{-}]}{10 [A^{-}] + [HA]}$$

### Buffer capacity in acid and base direction:

 $CH_3COOH(aq) + H_2O(l) = H_3O^+(aq) + CH_3COO^-(aq)$ 



# Buffer capacity in acid and base direction:



# Titration of buffer with acid and base:




Figure 6 (a) Ethanoic acid buffer with a strong base added (b) Ethanoic acid buffer with a strong acid added. The pH changes quickly once all of the available buffer is depleted.

# **Example:** Calculate the instantaneous ( $\beta$ ) and the practical buffer capacity in both directions of a 0.05 M Tricine buffer, pH 7.5, pKa = 8.15.

First calculate the concentration of the weak base and its conjugated acid that make up the buffer with 0.05M:  
7.5 = 8.15 +log 
$$\frac{[A-]}{[HA]}$$
  $\rightarrow$  -0.65= log  $\frac{[A-]}{[HA]}$   $\rightarrow$  Anti log for both sides  $\rightarrow$  0.224 =  $\frac{[A-]}{[HA]}$   
SO:  $[A-] = \frac{0.224}{1.224} \times 0.05 \equiv 0.009 \text{ M}$ ,  $[HA-] = \frac{1}{1.224} \times 0.05 = 0.041 \text{ M}$   
\* since the pH < pKa, the [HA] will be higher than [A<sup>-</sup>].

#### Calculate the instantaneous buffer capacity:


$$\beta = \frac{2.3 \text{ Ka [H^+][C]}}{(\text{Ka + [H^+])^2}} \Rightarrow \frac{2.3 \text{ x } 7.08 \text{ x } 10^{-9} \text{ x } 3.16 \text{ x } 10^{-8} \text{ x } 0.05}}{(7.08 \text{ x } 10^{-9} + 3.16 \text{ x } 10^{-8})^2} \Rightarrow \beta = 0.017 \text{ M}$$

$$\beta = \frac{2.3 \text{ [HA] [A^-]}}{[\text{HA}] + [\text{A}^-]} \Rightarrow \frac{2.3 \text{ x } 0.041 \text{ x } 0.009}{0.041 + 0.009} \Rightarrow \beta = 0.017 \text{ M}$$

#### **Example cont':** Calculate the instantaneous ( $\beta$ ) and the practical buffer

capacity in both directions of a 0.05 M Tricine buffer, pH 7.5, pKa = 8.15.

<u>Calculate the practical buffer capacity in both directions:</u>



#### Prestieel Pert



□ To understand the concept of buffer capacity.


□ To determine the buffer capacity in alkaline and acid directions.



- □ You are provided 0.1 M acetate buffer (pH=5).
- □ In two beakers add 8ml of the 0.1 M acetate buffer.
- Titrate the first beaker by adding 0.5 ml of 0.1 M HCl and the second one by 0.1 M NaOH from the burette and determine the pH of the solution after each addition.
- □ Continue adding the acid/base until you record a notable change in the pH.
- □ Record the titration table.



- Plot a curve of pH against the volume (ml) of HCl and NaOH added. calculate pH after addition of 0.5 ml, 2ml of HCl.
- □ Calculate the buffer capacity in both direction from the graph and the formula.
- Determine the buffering region.





- Compare between the value of the buffer capacity you got from the curve and formula.
- Did your buffer have a larger capacity for acid or base? Why?
- □ How can you relate your results with the buffer pH.