




BCH 312 Experiment (5)

Preparation of Buffer Solutions by Different Ways



+ Objectives




- 1) To learn how to prepare a buffer by different ways.



Introduction

For example: prepare **0.1 liters** of **0.045 M** sodium phosphate buffer, **pH 7.5** solution

Dissociation of phosphoric acid



- The pH of this buffer is a little above the pK_2 of H_3PO_4 , consequently, the two major ionic species present are H_2PO_4^- (conjugate acid) and HPO_4^{2-} (conjugate base) with the HPO_4^{2-} predominating { since the pH of the buffer is slightly basic }

+ The buffer can be prepared in any one of several ways :

1. By mixing NaH_2PO_4 (**conjugate acid**) and Na_2HPO_4 (**conjugate base**) in the proper proportions ,
2. By starting with H_3PO_4 and converting it to NaH_2PO_4 plus Na_2HPO_4 by adding the proper amount of NaOH ,
3. By starting with NaH_2PO_4 and converting a portion of it to Na_2HPO_4 by adding NaOH ,
4. By starting with Na_2HPO_4 and converting a portion of it to NaH_2PO_4 by adding a strong acid such as HCl ,
5. By starting with Na_3PO_4 and converting it to NaH_2PO_4 plus Na_2HPO_4 by adding HCl ,and
6. By mixing Na_3PO_4 and NaH_2PO_4 in the proper proportions .



- **Regardless of which method is used , the first step involves calculating the proportion and amounts of the two ionic species in the buffer .**

■ Total no. of moles of phosphate buffer= $M \times V = 0.1 \times 0.045 = 0.0045$ moles

- $\text{pH} = \text{pK}_{\text{a}_2} + \log \left[\frac{\text{HPO}_4^{-2}}{\text{H}_2\text{PO}_4^-} \right]$
- Assume $[\text{A}^-] = y$, $[\text{HA}] = 0.045 - y$
- $7.5 = 7.2 + \log (y / 0.045 - y)$
- $7.5 - 7.2 = \log (y / 0.045 - y)$
- $0.3 = \log (y / 0.045 - y) \rightarrow \text{antilog of } 0.3 = 2 = y / 0.045 - y$
- $Y = 0.09 - 2 y \rightarrow 3 y = 0.09$
- $Y = 0.9 / 3 = 0.03 \text{ M} = [\text{HPO}_4^{-2}] = [\text{A}^-]$
- $[\text{H}_2\text{PO}_4^-] = [\text{HA}] = 0.045 - 0.03 = 0.015 \text{ M}$

■ No. of moles of $\text{A}^- = M \times V = 0.03 \times 0.01 = 0.003$ moles


■ No. of moles of $\text{HA} = M \times V = 0.015 \times 0.01 = 0.0015$ moles



You are provided with concentrated (15M) H_3PO_4 and solution of 1.5 M NaOH .

Calculations:

- Start with **0.0045 mole** of H_3PO_4 and add 0.0045 moles of NaOH to titrate H_3PO_4 completely to give H_2PO_4^- (HA) , then add 0.003 moles of NaOH to titrate H_2PO_4^- to give HPO_4^{2-} (A^-)



- No. of moles needed of NaOH= $0.0045+0.003= 0.0075$ moles
 - Volume of NaOH needed= no.of moles / M = $0.0075/ 1.5 = 0.005 \text{ L} = \mathbf{5 \text{ ml}}$
 - Volume of H_3PO_4 needed =no.of moles / M = $0.0045/ 15 =0.0003 \text{ L} = \mathbf{0.3 \text{ ml}}$
- Add 5ml of NaOH to the 0.3 ml of concentrate H_3PO_4 , mix ; then add sufficient water to bring the final volume to 0.1 liters (100 ml), and check the pH



You are provided with solid $\text{NaH}_2\text{PO}_4(\text{HA})$ and solid NaOH

Calculations

- Start with **0.0045 mole** of NaH_2PO_4 and add **0.003** moles of NaOH to titrate NaH_2PO_4 to give $\text{Na}_2\text{HPO}_4 (\text{A}^-)$
- Wt of NaH_2PO_4 needed = no. of moles x mwt = $0.0045 \times 141.98 =$
0.638 g
- Wt of NaOH needed = no. of moles x mwt = $0.003 \times 40 =$ **0.12 g**
- Dissolve the NaH_2PO_4 and NaOH in some water, mix ; then add sufficient water to bring the final volume to 0.1 liters (100 ml), and check the pH



You are provided with solid Na_3PO_4 and 2 M solution of HCl

- Start with **0.0045 mole** of Na_3PO_4 and add 0.0045 moles of HCl to titrate Na_3PO_4 completely to give Na_2HPO_4 (A^-), then add **0.0015** moles of HCl to titrate Na_2HPO_4 to give NaH_2PO_4 (HA)

- No. of moles needed of HCl = $0.0045 + 0.0015 = 0.006$ moles

- Volume of HCl needed = $\text{no. of moles} / \text{M} = 0.006 / 2 = 0.003 \text{ L} = \mathbf{3 \text{ ml}}$

- Wt of Na_3PO_4 needed = $\text{no. of moles} \times \text{mwt} = 0.0045 \times 380.12 = \mathbf{1.71 \text{ g}}$


→ Dissolve 1.71 g of Na_3PO_4 in some water, mix ; then add 3 ml of HCl.

Finally, add sufficient water to bring the final volume to 0.1 liters (100 ml), and check the pH



You are provided with solid NaH_2PO_4 (HA) and Na_3PO_4

- The NaH_2PO_4 (HA) and Na_3PO_4 react to form Na_2HPO_4 (A^-). The **NaH_2PO_4 acts as an acid** and the **Na_3PO_4 acts as a base**.



- **Note that each mole of NaH_2PO_4 and Na_3PO_4 yields 2 moles of Na_2HPO_4 →** Thus to produce 0.003 mole of Na_2HPO_4 , 0.0015 mole of NaH_2PO_4 and 0.0015 mole of Na_3PO_4 are required.

** But, in addition to the 0.003 mole of Na_2HPO_4 , the final solution contains 0.0015 mole of NaH_2PO_4 .

Therefore, dissolve 0.0030 mole of NaH_2PO_4 and 0.0015 mole of Na_3PO_4 in water.

- Of the original 0.003 mole of NaH_2PO_4 , 0.0015 mole reacts with the Na_3PO_4 to produce 0.0030 mole of Na_2HPO_4 , leaving 0.0015 mole as NaH_2PO_4 .

-You need 0.0015 mole of K_3PO_4 .

- Wt of Na_3PO_4 = no. of moles \times mwt = $0.0015 \times 380.12 = 0.57 \text{ g}$

- Wt of NaH_2PO_4 = no. of moles \times mwt = $0.003 \times 141.96 = 0.425 \text{ g}$

- Dissolve the NaH_2PO_4 and Na_3PO_4 in some water, mix, then add sufficient water to make 0.1 liters (100 ml) and check the pH