

Lab sheet #2

-Preparation of solutions-

Objectives:

- To learn how to prepare solutions with different concentration expressions

Method:

(1) You are provided with solid NaOH, Prepare 50ml with 0.08M NaOH solution.

Calculation:

.....
.....
.....
.....
.....
.....

- To prepare the 0.08M NaOH solutiong of solid NaOH should be dissolved in a little volume of water then the volume made up toml, by the addition of water.

(2) You are provided with solid NaCl, Prepare 50ml with 1.5 w/v% solution of NaCl.

Calculation:

.....
.....
.....
.....
.....
.....

- To prepare the 1.5 w/v% solutiong of NaCl should be dissolved in little water and the volume made up toml by the addition of water.

(3) Prepare 100ml with 0.4 M HCl solutions starting with the concentrated HCl solution you are provided with: (w/w%= 36 , S.Gr =1.15).

Calculation:

.....

.....

.....

.....

.....

.....

- To prepare the 100ml of 0.4M HCl solutionml of stock (i.e. concentrated HCl) solution is needed and the volume made up toml by the addition of water.
- Measure and record the pH value of the acid you prepared.....
- Calculate the pH of the acid (pH= - log [H+])
- Determine your accuracy.....

Note: Atomic weights: Na = 23, Cl= 35.5, O = 16, H = 1

Formulas

$$\text{Molarity} = \frac{\text{moles of solute (mole)}}{\text{volume of solution in (L)}}$$

$$\text{Moles} = \frac{\text{weight (g)}}{\text{Molecular weight } \left(\frac{\text{g}}{\text{mole}}\right)}$$

$$\text{W/V}\% = \frac{\text{weight of solute in (g)}}{\text{volume of solution in (ml)}} \times 100$$

$$\text{W/W}\% = \frac{\text{weight of solute in (g)}}{\text{weight of solution in (g)}} \times 100$$

$$\text{Weight (wt)} = \text{volume (ml)} \times \text{SG} \times \text{w/w}\% \text{ (as decimal)}$$