Lab sheet #7

-Titration of a weak acid with a strong base-

Objectives:

- To study titration curves.
- Determine the pKa value of a weak acid.
- Calculate the pH value at a given point.
- Reinforce the understanding of buffers.

Method:

You are provided with **0.1M CH₃COOH** (weak acid solution) and **0.1M NaOH** (strong base solution):

- 1. Fill up the Burette with **0.1M NaOH** solution using the funnel.
- 2. To a beaker, add 10 ml of **0.1M CH₃COOH** solution and measure its pH value.
- 3. Start the titration: slowly add 0.5ml of **0.1M NaOH** (drop-wise) to the weak acid solution and <u>mix</u>, then record the pH value.
- 4. Keep on titration (step 3) until the <u>pH reaches 9.</u>

Results:

ml of 0.1M NaOH	pН	ml of 0.1M NaOH	pН
0		8.5	
0.5		9	
1		9.5	
1.5		10	
2		10.5	
2.5		11	
3		11.5	
3.5		12	
4		12.5	
4.5		13	
5		13.5	
5.5		14	
6		14.5	
6.5		15	
7		15.5	
7.5		16	
8		16.5	

- **1.** Record the pH values in the titration Table
- **2.** Plot the titration curve (pH versus ml of 0.1M NaOH added.
- 3. The calculations
 - a. <u>Calculate</u> the pH of the weak acid HA solution at 0 ml of NaOH and after the addition of 3ml, 5ml, and 10ml of 0.1M NaOH. (using the theoretical Pka = 4.76)
 - b. <u>Determine</u> pH value from the curve (measured)

(In the discussion: compare between pH values from a&b)

0.1M NaOH (ml)	Calculated pH	pH from titration curve
0		
3		
5		
10		

c. <u>Determine</u> the pKa value of weak acid from the curve (In the discussion: compare the calculated Pka value with the theoretical one)

In the Discussion

At what pH-range did the acid show the **best** buffering behavior? What are the chemical species at that region, what are their proportions?