

Blood Biochemistry BCH 471[Practical]

**Lab (2) Determination of Non-functional Plasma Enzymes  
in Serum**



# Blood Enzymes

---

- Plasma, serum or **blood proteins**, are proteins present in blood plasma which have several functions.
- Some blood proteins also act as **enzymes**.
- **Enzymes** are biocatalysts that increase the rate of the chemical reaction.
- **Clinical enzymology** refers to measurement of enzyme activity in body fluids for the diagnosis and treatment of diseases.
- Most clinical enzyme measurements using **serum** or **plasma**, occasionally other fluids, such as urine and gut secretions are also investigated.

# Differences Between Plasma Enzymes

## Plasma Enzymes

### 1. Plasma-specific Enzymes (Functional)

Enzymes that are normally present in the plasma and perform their primary function in the blood.

### 2. Non-plasma specific Enzymes (Non functional)


Intracellular enzymes that are normally present in very small amount in blood and perform no known function in blood.

	Functional plasma enzymes	Non functional plasma enzymes
<b>Their substrate</b>	Always present in the blood	Absent from the blood
<b>Site of synthesis</b>	Liver	Different organs e.g. liver, heart, muscles, and brain
<b>Effect of diseases in its plasma levels</b>	Decrease in liver diseases	Different enzymes increase in different organ diseases
<b>Examples</b>	Thrombin Plasmin Ceruplasmin	ALT LDH Acid Phosphatase Amylase

 **Pause and Think** Which of these enzymes is a better diagnostic indicator ? Why?

# Sources of Non functional Plasma Enzyme

1. **Cell damage** with the release of its content of enzymes into blood e.g. **Myocardial infarction** and **viral hepatitis**.
2. **Block in the secretory pathway** e.g. elevation of blood pancreatic amylase and lipase in **pancreatitis**.
3. **Increase enzyme synthesis** e.g. elevation of serum alkaline phosphatase in **bone cancer**.
4. **Increased permeability of cell membrane** as in hypoxia.



So estimation of the plasma concentration of these enzymes in blood is useful for the diagnosis of disease depending on their tissue origin.

# Clinical Significance of Non-Functional Plasma Enzymes


---

**Measurement of non-functional enzymes is important for:**

- 1. Diagnosis of diseases.**
- 2. Prognosis of the disease:** following up of the treatment by measuring plasma enzymes before and after treatment.

# Lactate Dehydrogenase (LDH)

- LDH is a hydrogen transfer enzyme which catalyzes the **interconversion of pyruvate and lactate** with the mediation of **NAD<sup>+</sup>** as hydrogen acceptor, eventually **converting pyruvate to glucose**.
- The optimum pH for lactate pyruvate (L→P) reaction is **8.8 – 9.8**, While for pyruvate to lactate (P→L) is **7.7 – 7.8**.
- The enzyme is inhibited by **sulphydryl reagents** and **mercuric ions**.



# Lactate Dehydrogenase (LDH)


- **Lactic acid dehydrogenase (LDH)** is a cytoplasmic enzyme that helps produce energy.
- LDH is most often measured to evaluate the presence of tissue damage (diagnostic).
- The enzyme LDH is found in many body tissues, especially the **heart, liver, kidney, skeletal muscle, brain, blood cells and lungs.**

High LDH Level in the Plasma in Different Diseases	Diseases	Examples
	Myocardial infarction	
Liver Disease		Toxic jaundice
		Viral hepatitis
		Obstructive jaundice
Anemia		Pernicious anemia
		Megaloblastic anemia
Renal Diseases		Tubular necrosis
		Pyelonephritis
Malignant Disease		Lung Cancer
		Hodgkin's disease

# LDH Isoenzymes

- LDH exists in **5 forms** (isoenzymes), which differ slightly in structure.
- All of these isoenzymes can be measured in the blood, and can be separated by electrophoresis.


Isoenzyme	Cellular level	Diseases associated
LDH-1	<b>Heart tissue</b>	Myocardial infarction
LDH-2	<b>WBC (monocytes), RBC</b>	<ul style="list-style-type: none"> <li>➤ Megaloblastic anemia</li> <li>➤ leukemia</li> </ul>
LDH-3	<b>Lung tissue</b>	Pulmonary embolism
LDH-4	<b>kidneys, placenta and pancreas</b>	Pancreatitis
LDH-5	<b>liver and skeletal muscle</b>	<ul style="list-style-type: none"> <li>➤ Toxic hepatitis with jaundice</li> <li>➤ Muscular dystrophy</li> </ul>





# Alanine Transaminase (ALT)

- ALT is an enzyme that catalyzes a type of reaction (**transamination**) between an amino acid and  $\alpha$ -keto acid.
- It is important in the production of various amino acids.
- Also called alanine transferase (**ALT**), serum glutamate-pyruvate transaminase (**SGPT**).



- **Transamination reaction** is the process by which **amino** groups are removed from **amino acids** and transferred to acceptor **keto-acids** to generate the amino acid version of the keto-acid and the keto-acid version of the original amino acid.

# ALT Diagnostic Importance

---

- Normally, high concentrations of ALT occur in the liver, and relatively low concentrations are found in the heart, muscle, and kidney.
- In the **serum**, only **low level** (10–35 U/L) of the ALT is found, thus an elevated level is a sensitive index of acute hepatocellular injury.
- Elevated serum ALT (SGPT) level are found in **hepatitis, cirrhosis** and **obstructive jaundice**.
- Levels of ALT are only slightly elevated in patient following a **myocardial infarction**.

# Practical Part

---



# Objectives

---

- To determine the level of **Lactate Dehydrogenase (LDH)** in serum
- To determine the level of **Alanine Transaminase (ALT)** in serum.
- To evaluate the presence of tissue damage.

# Lactate Dehydrogenase Assay

---

## Principle

LDH catalysis the following reaction:



- The rate of NADH formation is indicated by **increase the absorbance at 340 nm** and it is directly proportional to serum LDH activity.

If:

- NADH is **product** → **increase** the absorbance/min
- NADH is **reactant** → **decrease** the absorbance/min

# Method

---

Pipette into cuvettes at 30 °C:

**Serum Sample**

20  $\mu$ l

**Buffer**

1000  $\mu$ l

\*Containing (TRIS buffer, Pyruvate)

Mix and incubate at 30 °C for 3 minutes

**Substrate**

250  $\mu$ l

\*Containing (NADH)

Mix and incubate at 30 °C for 1 minute, then read the absorbance at 340 nm against distilled water (blank) every minute for 2 minutes. Then determine  $\Delta A/\text{min}$ .

Measure enzyme kinetics using UV-visible spectroscopy:

2) Applications  $\rightarrow$  2) Simple Kinetics  $\rightarrow$  wavelength (340 nm)  $\rightarrow$  1) Seconds  $\rightarrow$  Duration (120 sec)

$\rightarrow$  Intervals (60 sec)  $\rightarrow$  Print Data Table (off)  $\rightarrow$  Press start (2 times)

# Results and Calculations

---

## Results

	Time (min)	Absorbance at 340 nm
A <sub>1</sub>	0	
A <sub>2</sub>	1	
A <sub>3</sub>	2	


## Calculations

- $\Delta A_1 = A_1 - A_2$        $\Delta A_2 = A_2 - A_3$
- $\Delta A/\text{min} = (\Delta A_1 + \Delta A_2) / 2$
- **LDH Activity (U/L) =  $\Delta A/\text{min} \times 10080$**

**Normal Values:** Adults: 160 – 320 (U/L)

# Alanine Transaminase Assay

## Principle



- The rate of NAD<sup>+</sup> formation is indicated by **decreased the absorbance at 340 nm** and it is indirectly proportional to serum LDH activity.

### If:

- NADH is **product** → **increase** the absorbance/min
- NADH is **reactant** → **decrease** the absorbance/min



# Method

---

Pipette into cuvettes at 37 °C:

**ALT reagent**

1000  $\mu$ l

\*Containing (L-Alanine, Oxoglutarate LDH, NADH, buffer)

Pre-warm at 37 °C for 3 minutes and add

**Sample (serum)**

100  $\mu$ l

Mix and incubate at 37 °C for 1 minute, then read the absorbance at 340 nm against distilled water (blank) every minute for 2 minutes. Then determine  $\Delta A/\text{min}$ .

---

Measure enzyme kinetics using UV-visible spectroscopy:

2) Applications  $\rightarrow$  2) Simple Kinetics  $\rightarrow$  wavelength (340 nm)  $\rightarrow$  1) Seconds  $\rightarrow$  Duration (120 sec)

$\rightarrow$  Intervals (60 sec)  $\rightarrow$  Print Data Table (off)  $\rightarrow$  Press start (2 times)

---

# Results and Calculations

---

## Results

	Time (min)	Absorbance at 340 nm
A <sub>1</sub>	0	
A <sub>2</sub>	1	
A <sub>3</sub>	2	

## Calculations

- $\Delta A_1 = A_1 - A_2$        $\Delta A_2 = A_2 - A_3$
- $\Delta A/\text{min} = (\Delta A_1 + \Delta A_2) / 2$
- **ALT Activity (U/L) =  $\Delta A/\text{min} \times 1768$**

**Normal Values:** Males: 10-40 (U/L). Female: 7- 35 (U/L).