

BCH 462- Biotechnology & Genetic engineering [Practical]
Lab (1) Plasmid Isolation and Purification



Let's assume we want to make **insulin** for the **treatment of diabetes**.

- Which is better to use **cow** OR **bacteria** as a biological factory? Why?

→ How would the cell produce **human insulin** ?




Figure 1. Schematic representation of recombinant insulin production

DNA cloning techniques

Are techniques used to create copies of certain DNA fragments.


1- **PCR** (in vitro)


[polymerase chain reaction].

2- **Cell-based** (in vivo)

[using a vector e.g. **plasmid** carrying the DNA of interest, which eventually inserted to a host cell “usually bacteria” and self replicate].

Remember

Vector  e.g. **Plasmid**

Host  e.g. **Bacteria**

Plasmid

- The DNA of most bacteria is contained in a single circular molecule, called the **bacterial chromosome**.
- Many bacteria contain an extrachromosomal element of DNA, termed a **plasmid**.
- **Plasmid** is a relatively small, covalently closed circular molecule that replicate independently from a bacterial chromosome. (Why ?)
- Every plasmid has its own origin of replication (**replicon**) and use the enzymes and proteins that encoded by its host for its replication and transcription.
- Plasmid found in a wild variety of bacterial species and they are not essential for the bacterium but benefit the survival of the organism (**Symbiotic** relationship with the host ?).




Figure 2. Illustration of a bacterium showing chromosomal DNA and plasmids

Plasmid cont.

❑ Plasmids classes:

- I. **Virulence** plasmids encoding toxin genes.
- II. **Drug-resistance** plasmids that confer resistance to antibiotics.
- III. Plasmids encode **genes** required for bacterial conjugation.
(which can be advantageous for host cell)

❑ Plasmids applications:

- i. **Molecular cloning**
- ii. Gene therapy
- iii. Drug production
- iv. Making a large amount of proteins.




Figure 3. Illustration of *E.coli* showing chromosomal DNA and plasmids

Plasmid cont.




Figure 5. TEM picture of plasmid DNA molecules in solution




Figure 6. Electron micrograph of an *E. coli* cell ruptured to release its DNA.

Plasmid as a vector

- Plasmids are widely used as **vectors** in **molecular cloning**, serving to drive the replication of **recombinant DNA** sequences within host organisms (It is used to provide a “vehicle” in which to insert a desired DNA fragment).
- Recombinant DNA**, molecules of DNA from **two different species** (human/bacteria) that are inserted into a **host organism** (bacteria) to produce new genetic combinations (human insulin).
- In the laboratory, the modified plasmids (**recombinant DNA**) are usually reintroduced into a host cell for replication via process called *transformation*.




Figure 7. Recombinant DNA





Figure 8. Transformation

Plasmid vectors should contain three important parts:


1. Origin of replication (Ori)
2. antibiotic resistance
3. gene cloning site

The **Ori** is a DNA sequence which allows initiation of replication of the plasmid by cellular enzymes.



Cloning site is a short segment of DNA which contains several restriction sites allowing for the easy insertion of DNA (A place to insert foreign DNAs).

Antibiotic resistance gene allows for selection of plasmid-containing bacteria. (?)



Plasmid isolation and purification

- Is an essential step for many molecular biology procedures.
- In general, plasmid purification involved three steps:
 1. **Growth** of the bacterial culture.
 2. **Harvesting** and **lysis** of bacteria.
 3. **Purification** of plasmid DNA.

1. Growth of the bacterial culture

Depending upon nutritional status, bacteria exhibit **different growth patterns** which include:

- I. **Lag phase:** in this phase bacteria adapt themselves to growth conditions and synthesis its own DNA, RNA and proteins.
- II. **Log phase:** it is exponential phase, the bacterial cells divide and the production of new cells is proportion to increased time.
- III. **Stationary phase:** the growth rate slows as nutrients become limited, waste products accumulate and the rate of cell division equals the rate of death.
- IV. **Death phase:** due to continuous accumulation of toxic metabolites and the lack of nutrients, death occurs of the bacteria.




Figure.9. Bacterial culture growth curve.

💡 **Pause and Think** on which phase should we purify plasmid?

2. Harvesting and lysis of bacteria

1. Bacteria are recovered by **centrifugation**.
2. **Cell lysis** by any one of many methods, including:
 - Treatment with **detergents, alkali, organic solvents, and heat**.
 - **The choice among these methods depends on three factors:**
 - The size of plasmid.
 - The bacterial strain.
 - The technique used to subsequently purify the plasmid DNA.

3. Purification of plasmid DNA

- The plasmid purification procedures, unlike the procedures for purification of genomic DNA, should involve removal of not only protein but also another major impurity **bacterial chromosomal DNA**.
- **There are basic methods of plasmid preparation:**
 1. Chemical base lysis methods.
 2. Application of affinity matrices for plasmid or proteins.

Practical Part


Practical part

▪ Aim:

- To isolate pure plasmid DNA from **E. coli** using **alkaline lysis method**.

▪ Principle:

- In the alkaline lysis method, cells are lysed and DNA denatured by **SDS** and **alkaline pH**.
- The **SDS** will lyse the bacterial cell membrane and denature the proteins.
- **Alkaline pH** will denature the genomic DNA and the proteins too.
- **Neutralization** of the solution.
- **Precipitation** of protein-SDS complexes.
- Subsequently both complexes, DNA and protein, are removed by **centrifugation** leaving native plasmid molecules in the supernatant.






Figure.10. Alkaline lysis purification method performing steps

Practical part

- **Results:**

- Concentration of plasmid DNA (ng/ μ l) = _____
- Plasmid purity: A260/A280 = _____

- **Methodology:**

1- Centrifuge the bacterial samples at 4 °C, maximum speed for 5 minutes, using microcentrifuge device.

Bacterial sample was centrifuged at 4°C, maximum speed for 5 minutes using microcentrifuge.

- **References:**

Endnote, Mendeley or Cite This For Me: Web Citer (*extension in Google Chrome*).