


KINETICS ANALYSIS OF B-FRUCTOFURANOSIDASE ENZYME

2-The effects of PH value on the rate of a reaction catalyzed by β - fructofuranosidase.

- Enzymes are affected by changes in pH. The most favorable pH value - the point where the enzyme is most active is known as the optimum pH.
- Extremely high or low pH values generally result in complete loss of activity for most enzymes. pH is also a factor in the stability of enzymes. As with activity, for each enzyme there is also a region of pH optimal stability.





Optimum pH: is the pH at which the rate of reaction is maximum.

At **higher or lower** pH, the rate of an enzymatic reaction **decrease**.


For most enzymes, the optimum pH lies in the range from pH 5 to pH 9


The bell shape of pH activity curve is resulted from the following factors:

- Enzyme denaturation at extremely high or low pH
- Effects on the charged state of the substrate or enzyme.



**Different Enzymes
with Different
optimum PH**





For the majority of enzymes, the relationship between the rate of an enzymatic reaction and pH takes form of a **bell-shape**.


Objective

- a) To establish the relationship between pH and the rate of an enzyme catalyzed reaction.

- b) To determine the optimum pH for β -Fructofuranosidase enzyme.

Principle

Within acidic environment using acetate buffer (PH= 4.7) **β -fructofuranosidase** enzyme cleavage its substrate (Sucrose) non reducing sugar to mixture of reducing sugar glucose and fructose, using 3,5,dinitrocylic acid .



Material

Solutions :-

- 0.05M Sodium Acetate buffer , pH 4.7 .
- 0.18 M Sucrose ,
- Reducing sugar (0.005M glucose + 0.005M fructose)
- Beta Fructofuranosidase (Invertase) enzyme extract from yeast.
- DNS (dinitrosalicylic acid) Reagent .
- Sodium Bicarbonate .

Method:

1- Prepare 6 tubes of different by following the table provided. (1) :

| Tube | Buffer 1.0 (ml) | Sucrose (ml) |
|---------|-----------------|--------------|
| A | pH 3.0 | 2.0 |
| Blank A | pH 3.0 | 2.0 |
| B | pH 4.7 | 2.0 |
| Blank B | pH 4.7 | 2.0 |
| C | pH 7.0 | 2.0 |
| Blank C | pH 7.0 | 2.0 |
| D | pH 8.0 | 2.0 |
| Blank D | pH 8.0 | 2.0 |
| E | pH 12.0 | 2.0 |
| Blank E | pH12.0 | 2.0 |

Table (1).

- 2- Mix each tube properly then incubate all tubes at 40°C for 5min .
- 3- Start the reaction by adding 0.05ml of diluted enzyme to all test tubes except for the blanks add 0.05ml of distilled water instead , mix and start the stop clock immediately , incubate each tube for 10min , then stop the reaction by adding 2.0ml of the DNS reagent to each tube and mix well .(follow table 2 for adding enzyme and DNS to tubes) .
- Note : Mix each tube frequently during the incubation time .

Method:

| Tube | Start Time (min) | Stop by adding 2.0ml DNS . (min) |
|---------|------------------|----------------------------------|
| Blank A | 0 | 10 |
| Blank B | 1.0 | 11 |
| Blank C | 2.0 | 12 |
| Blank D | 3.0 | 13 |
| Blank E | 4.0 | 14 |
| A | 5.0 | 15 |
| B | 6.0 | 16 |
| C | 7.0 | 17 |
| D | 8.0 | 18 |
| E | 9.0 | 19 |

Table(2)

Method:

- 4- Mix properly , cover each tube by aluminum foil and place in a boiling water bath for 5min to allow the color to develop .
- 5- Remove from water bath cool under tap water , add 20ml of distilled water to each tube , mix properly then measure the absorbance at 540nm .
- 6- Record the absorbance of each test tube in the following table (3),
- 7- Convert the Absorbance reading obtained to micromoles of sucrose hydrolyzed making use of the standard reducing sugars calibration curve , determine the initial velocity v_i for each tube and record all in table 3 .
- 8- Obtain the relationship between the initial velocity v_i and pH , by drawing a graph between the initial velocity v_i and pH . Determine the optimum pH for your enzymatic reaction reaction.

Result

Plot velocity against PH value. Describe the shape of this curve and discuss the reasons for its shape.

| Tube | Absorbance 540nm | μ moles of sucrose hydrolyzed | μ moles of sucrose hydrolyzed/min(vi) |
|------|------------------|-----------------------------------|---|
| A | | | |
| B | | | |
| C | | | |
| D | | | |
| E | | | |

Table (3)

Discussion

Comment on the curve shape and conclude the relationship between PH value and the rate of an enzyme catalyzed reaction.



Thank You