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
READING THE API 20

TESTS	SUBSTRATE	REACTION TESTED	- RESULTS	+ RESULTS
ONPG	ONPG	beta-galactosidase	colorless	yellow
ADH	arginine	arginine dihydrolase	yellow	red/orange
LDC	lysine	lysine decarboxylase	yellow	red/orange
ODC	ornithine	ornithine decarboxylase	yellow	red/orange
CIT	citrate	citrate utilization	pale green/yellow	blue-green/blue
H2S	Na thiosulfate	H2S production	colorless/gray	black deposit
URE	urea	urea hydrolysis	yellow	red/orange
TDA	tryptophan	deaminase	yellow	brown-red
IND	tryptophan	indole production	yellow	red (2 min.)
VP	Na pyruvate	acetoin production	colorless	pink/red (10 min.)
GEL	charcoal gelatin	gelatinase	no diffusion of black	black diffuse
GLU	glucose	fermentation/oxidation	blue/blue-green	yellow
MAN	mannitol	fermentation/oxidation	blue/blue-green	yellow
INO	inositol	fermentation/oxidation	blue/blue-green	yellow
SOR	sorbitol	fermentation/oxidation	blue/blue-green	yellow
RHA	rhamnose	fermentation/oxidation	blue/blue-green	yellow
SAC	sucrose	fermentation/oxidation	blue/blue-green	yellow
MEL	melibiose	fermentation/oxidation	blue/blue-green	yellow
AMY	amygdalin	fermentation/oxidation	blue/blue-green	yellow
ARA	arabinose	fermentation/oxidation	blue/blue-green	yellow
OX	oxidase	oxidase	colorless/yellow	violet

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