

# Copper Reduction Method (Somogyi–Nelson Method)

Done by: Sahar AL-subaie

# Principle:

- ❖  $\text{Cu}^{+2} + \text{glucos} \longrightarrow \text{Cu}^{+1} + \text{oxidation products of glucose}$
- ❖  $\text{Cu}^{+1} + \text{Phosphomolybdic acid}(\text{MO}^{+6}) \longrightarrow \text{MO}^{+4} \text{ (blue color)}$
- ❖ Then the concentration of  $\text{MO}^{+4}$  is measured using spectrophotometer, at 520 nm,
- ❖ N.B. the 1<sup>st</sup> test is done under alkaline pH and heat conditions.

# Procedure:

unknown	Blank	St. 1	St. 2	St. 3	St. 4	St. 5
Distilled water						
x	1 ml	0.8 ml	0.7 ml	0.5 ml	0.3 ml	x
Diluted Glucose standard (0.2 mg/ml)						
x	x	0.2 ml	0.3 ml	0.5 ml	0.7 ml	1.0 ml
Protein-free filtrate						
1 ml	x	x	x	x	x	x
Copper reagent						
1 ml	1 ml	1 ml	1 ml	1 ml	1 ml	1 ml
Incubate in boiling water for 20 min., then cool.						
Arsenomolybdate Reagent						
1 ml	1 ml	1 ml	1 ml	1 ml	1 ml	1 ml
Let it stand for 1 min.						
Distilled water						
7 ml	7 ml	7 ml	7 ml	7 ml	7 ml	7 ml

N.B.

Reoxidation of cuprous ion by oxygen from the air is prevented by adding **Sodium Sulphate** in the reagent to decrease the solubility of oxygen.



## N.B.

- ❖ Due to the interference of proteins in this method, the proteins are precipitated by the addition of **Barium hydroxide and Zinc sulphate.**
- ❖ Protein is removed as Zinc proteinate

# Just for your information

Sulphydryl compounds as Zinc salts and the remaining zinc and barium ions as **Zinc hydroxide** and **Barium sulphate**.



N.B.

Recently, rapid colorimetric procedures using O-toluidine or enzymes (such as glucose oxidase and peroxidase) have replaced the Somogyi–Nelson method.



# Calculations:

$$C_1 \times V_1 = C_2 \times V_2$$