



322 BCH

EXP (7)

THE EFFECT OF SUBSTRATE
CONCENTRATION ON THE
RATE OF AN ENZYME
CATALYZED REACTION



In this experiment, we will continue to study acid phosphatase kinetics.




OBJECTIVES

- To establish the relationship between substrate concentration and the rate of an enzyme catalyzed reaction.
- To determine the K_m and V_{max} of the enzyme for a particular substrate.

THE EFFECT OF SUBSTRATE CONCENTRATION


- One of the important parameters affecting the rate of a reaction catalyzed by an enzyme is the substrate concentration, $[S]$.
- During enzyme substrate reaction, the initial velocity V_0 gradually increases with increasing concentration of the substrate. Finally a point is reached, beyond which the increase in V_0 will not depend on the $[S]$.



THE EFFECT OF SUBSTRATE CONCENTRATION


- At relatively low concentration of substrate, the rate of reaction increase linearly with an increase in substrate concentration.

→ The catalytic site of the enzyme is empty, waiting for substrate to bind, for much of the time, and the rate at which product can be formed is limited by the concentration of substrate which is available.



THE EFFECT OF SUBSTRATE CONCENTRATION


- At higher substrate concentration, the rate of reaction increases by smaller and smaller amounts in response to increase in substrate concentration.



THE EFFECT OF SUBSTRATE CONCENTRATION

- *However beyond a particular substrate concentration, the velocity remains constant without any further increase. This plateau is called the maximum velocity, V_{max}*
- This is because as the concentration of substrate increases, the enzyme becomes saturated with substrate.

So there is usually a **hyperbolic** relationship between the rate of reaction and the concentration of substrate



MICHAELIS-MENTEN EQUATION

- Michaelis-Menten equation give the relationship between $[S]$ and velocity of enzymatic reaction.
- The hyperbolic shape of this curve can be expressed algebraically by the Michaelis – Menten equation:

$$V = \frac{V_{\max} [S]}{K_m + [S]}$$


V_i = intial velocity

V_{\max} = maxiumm velocity

$[S]$ = substrate concentration


K_m = Michaelis constant

Michaelis Menten Plot



MICHAELIS CONSTANT (K_m)


- K_m is the substrate concentration at half V_{max} .
- The relationship between rate of reaction and concentration of substrate depends on the affinity of the enzyme for its substrate. This is usually expressed as the **K_m** of the enzyme, an inverse measure of affinity
- **The larger the k_m , the weaker the binding and the larger the $[S]$ needed to reach the half the maximum rate**
- The K_m can vary greatly from enzyme to enzyme, and even for different substrates of the same enzyme



LINEWEAVER – BURK EQUATION



- The Michaelis -Menten equation can be algebraically transformed into forms that are useful in the practical determination of K_m and V_{max} .
- One common transformation is derived simply by taking the reciprocal of both sides of the Michaelis -Menten equation to give Lineweaver – Burk equation:

$$\frac{1}{V_i} = \frac{1}{V_{max}} + \frac{K_m}{V_{max}} \cdot \frac{1}{[S]}$$



- By plotting $1/v$ against $1/[S]$ a straight line plot, Lineweaver – Burk plot is obtained.
- Both V_{max} and K_m can be obtained accurately from intercepts of the straight line with the y – axis and x-axis

PRINCIPAL



Its concentration can be measured at **405 nm**.

Method:


Place in a water bath maintained at 37 °C for 5 minutes.

[S] mM 0 0.5 1 2.5 5 10 25 50



Add to each tube

- 0.5 ml Corresponding pNPP



Add to each tube

- 0.5 ml of buffer
- 0.5 ml MgCl_2
- 5 ml water


All the factors that affect enzyme kinetics are constant [S] where it varies in
each tube

Time = 5 min Temp= 37 °C pH= 5.7

Method:

To start the reaction add add 0.5ml of enzyme

To stop the reaction → add 0.5ml of KOH



Start at	0	0	2	4	6	8	10	12
Stop at	0	5	7	9	11	13	15	17

For Blank: Add KOH ((FIRST)) then E, **to prevent the reaction from happening.**

After all the reactions have been terminated, determine the absorbance at **405 nm** for each sample against blank.



Tube	A	B	C	D	E	F	G
Start at	0	2	4	6	8	10	12
Stop at	5	7	9	11	13	15	17

Time (min)	Tube	Addition
0	A	Enzyme
2	B	Enzyme
4	C	Enzyme
5	A	KOH
6	D	Enzyme
7	B	KOH
8	E	Enzyme
9	C	KOH
10	F	Enzyme
11	D	KOH
12	G	Enzyme
13	E	KOH
15	F	KOH
17	G	KOH

To convert the time table to an easier way try the following



RESULTS

Record the absorbance and calculate the velocity:

Tube	Absorbance at 405 nm	Velocity $\mu\text{M}/\text{min}$	[S] mM
A			0
B			0.5
C			1
D			2.5
E			5
F			10
G			25
H			50

In order to draw using Lineweaver – Burk:
calculate the following:

[illegible]

CALCULATIONS

$$\text{Velocity (V)} = (A \times 10^6) / (E \times \text{time}) = \text{ }\mu\text{mole of PNP/min}$$

A= absorbance

E= extension coefficient= 18.8×10^3

Time = 5 min

RESULTS

- Draw the curve using Michaelis -Menten and determine V_{max} and K_m for acid phosphatase.
- Prepare the double –reciprocal plot of Lineweaver and Burk and determine the K_m and V_{max} from the x and y intercepts.

DISCUSSION

- Describe the curve of the effect of substrate concentration on the enzymatic activity
- Comment on the value of V_{max} and K_m and define each of them