

Isolation of Glycogen from Camel's liver

BCH 447


Objective:

- ▶ To illustrate the method for isolating glycogen.




Glycogen

- ▶ Glycogen is the main polysaccharide energy reserve in animals which is analogous to the starch in plants.
- ▶ It is the storage form of glucose.
- ▶ It is synthesized and stored mainly in **liver and muscles**.
- ▶ Its concentration **in the liver is higher** than in muscles.



Structure of glycogen

- It is a very large, branched polymer of glucose residues.
- Various samples of glycogen have been measured at 1,700-600,000 units of glucose.



Why does the body store glycogen?

- ▶ The breakdown of glycogen and release of glucose increase the amount of glucose that is available between meals.
- ▶ Hence, glycogen **serves as a buffer** to maintain blood-glucose levels.
- ▶ Maintaining blood-glucose levels is especially important because glucose is virtually the only fuel used by the *brain*, except during prolonged starvation.
- ▶ Moreover, the glucose from glycogen is readily metabolized and is therefore a good source of energy for sudden activity.



Principle of the experiment:

Grinding with TCA and sand

- Trichloroacetic acid (TCA) is used to precipitate large molecules by **denaturation**
- Sand helps in the grinding process

Centrifugation


- Proteins and nucleic acids are precipitated
- Glycogen remain in the supernatant with other sugars and water soluble substances

Alcohol Precipitation

- Glycogen is separated by precipitation with alcohol, **polysaccharides are less soluble than sugars in aqueous alcohol.**

Procedure:

1



5.0 g cold liver +
0.5 g sand +
5ml 10%TCA

Centrifuge

at 3,000 rpm for 5min



Transfer
into 50ml



Why cold? because low temp. inhibit the enzymatic hydrolysis of glycogen

2




Rinse out the
rest with 5ml
of 5% TCA

Centrifuge

at 3,000 rpm for 5min



Transfer
into 50ml



Discard




CAUTION

*TCA causes severe burns;
wash accidental spills on skin
with plenty of running tap
water.*

3


- Record total volume.
- Add equal volume of 95% ethanol, slowly with stirring.
- Allow to stand while precipitate settles.

Procedure:



Calculations:

Glycogen content (g) =



How to make glycogen solution?

Dissolve glycogen pellet in **phosphate buffer** (32mg glycogen in 4 ml phosphate buffer)

- Example:

I have 3.52 g glycogen ??

A- $3.52\text{g} \times 1000 = 3520\text{ mg}$

B- $32\text{mg} \rightarrow 4\text{ ml}$

$3520\text{mg} \rightarrow ???\text{ ml}$



Results:

Record total **yield** and **glycogen content** per 100 g liver.

Example:

Liver weight 5 g - The glycogen content 1.5 g

** the glycogen **yield** was 1.5 g

** the **glycogen content** per 100 g is

1.5 g → 5 g liver

??g → 100 g liver

= $1.5 \times 100 / 5 = 30\text{g} / 100 \text{ g liver}$

factors affecting glycogen yield:

-**Food state:** Fed animals get more glycogen yield

-**Stress:** Any stress for animal will lead to glycogen consumption .

