



The effect of incubation time on the rate of an enzyme catalyzed reaction




Objectives

- To monitor the progress of an enzyme catalyzed reaction (Acid phosphatase).
- To determine the initial rate of the reaction (V_i).
- **Important terms and points:**



Enzyme kinetics

- The central approach for studying the mechanism of an enzyme-catalyzed reaction is to study enzyme kinetics.
- It **determines** the rate of the enzymatic reaction (**velocity**) and its changes in response with the changes in **parameters** such as **substrate concentration**, **enzyme concentration**, **pH**, **temperature**.




Measuring the rate of reaction (velocity)

The progress of an enzyme catalyzed reaction may be followed by measuring either the

- amount of substrate consumed,
 - or Amount of product formed
- per unit time


Velocity (V) = rate of reaction = change in [P] or [S] per unit time

Unit : $\mu\text{moles/minute}$



Initial velocity


- **The initial rate of reaction (V_i)** is used in the study of enzyme kinetics
- It is measured as the **slope** of the curve at the origin (time= 0) measured as the slope at the origin (time= 0).



Progress curve is a graphic representation of a enzyme-catalyzed reaction in which the product or the substrate conc. is plotted against time

Why we measure initial velocity ?

- **The rate of the reaction is highest at time zero** and decreases with increasing time, eventually falling to zero itself, reaching a **plateau**.
- This usually occurs either when all the substrate is used up or when equilibrium is reached.
- So V_i reflects the enzyme activity in the beginning of the reaction at optimum parameters where the substrate conc. is maximum.



Acid phosphatase




- In this experiment, you will measure the velocity of the reaction catalyzed by purified acid phosphatase (**ACP**) from [wheat germ](#).
- Acid phosphatase is a **phosphatase** that acts on **monoesters** of orthophosphoric acid. It does **not act** on phosphoric diesters or triesters.
- Acid phosphatase is an important enzyme that plants use to **obtain their energy** from hydrolyze phosphate esters.

Enzyme Assay of ACP





1. The enzyme catalyzes the hydrolysis of **p-nitrophenyl phosphate** to **inorganic phosphate** and **p-nitrophenol**, under **acidic conditions (pH=5.7)**, with **optimum temperature 37°C**

- Mg^{++} ion act as an activator for the enzyme



Principal of the enzyme assay in vitro



2. Both p-nitrophenyl phosphate and p-nitrophenol are colorless at acidic pH values.
3. Addition of **alkaline solution**, p-nitrophenol is converted to a **p-nitrophenolate** (yellow color) and concentration can be measured at **405 nm**.

Note: Since ACP works under acidic condition, adding alkaline will cause the **enzymatic reaction to stop**.

Method: The effect of incubation Time

- The reaction is stopped after 5 minutes intervals (by Addition of KOH), hence called **Fixed-time assay**




Method:

- Prepare a series of **seven reaction tubes** labeled 0 through 30 minutes at 5 minute intervals (Blank, 5, 10 ... minutes).
- Follow the following addition protocol for **all the tubes**:

Add in each tube:


- 0.5 ml of buffer
- 0.5 ml of MgCl_2
- 0.5 ml of pNPP (S)
- 5 ml of water



Place in a water bath maintained at 37 °C for 5 minutes

To start the reaction → add 0.5ml of E

To stop the reaction → add 0.5ml of KOH



Start at (min)	0	0	2	4	6	8	10
Stop at (min)	0	5	12	19	26	33	40



Add KOH to blank ((FIRST)), to prevent the reaction from happening.

After all the reactions have been terminated, determine the absorbance at **405 nm** for each sample against blank.

Tube	5	10	15	20	25	30
Start at (min)	0	2	4	6	8	10
Stop at (min)	5	12	19	26	33	40



Time (min)	Tube	Addition 0.5 ml of
0	5	Enzyme
2	10	Enzyme
4	15	Enzyme
5	5	KOH
6	20	Enzyme
8	25	Enzyme
10	30	Enzyme
12	10	KOH
19	15	KOH
26	20	KOH
33	25	KOH
40	30	KOH



To convert the time table to an easier way try the following




Results

Time (min)	Absorbance at 405nm	[P] (μmolar)
0		
5		
10		
15		
20		
25		
30		

$[P] = (A \times 10^6) / E \times L = \mu\text{molar}$
E= extension coefficient= 18.8×10^3
A= absorbance
L= path length (1 cm)

Results

- Draw the graph [do not forget the title and the units]
 - X axis = Time
 - Y axis = [P]
 - Calculate from the graph the value initial velocity



Discussion

- introductory statement
- Comment on the graph that you get [in details](#)
- Mention the initial velocity value you got with **unit**, and why we use V_i



Questions

- How is the initial velocity of an enzymatic reaction determined?
- What are the factors that influence the rate of enzyme catalyzed reactions?