



The effect of temperature on the rate of an enzyme catalyzed reaction



- In this experiment, we will continue to study acid phosphatase kinetics




Objectives

- To establish the relationship between temperature and the rate of an enzyme catalyzed reaction.
- To determine the optimum temperature for such a reaction.

Optimum Temperature


- The rate of an enzyme catalyzed reaction is affected by changes in temperature.
- Each enzyme has a temperature that it works optimally in called (**Optimum Temperature**).
- For most enzymes, the optimum temperature is at the temperature of the cells in which the enzyme is found in **vivo**.
- For example, humans' enzymes works optimally **at around 37** degrees Celsius.
- However, some enzymes from hot spring bacteria work really well in higher temperatures like 95 degrees Celsius.



What is The effect of Temperature on the activity of Enzyme?


1-At suboptimal temperature, increasing temperature increase the kinetics energy of the reactant. As they move faster, they move likely to collide and interact with each other and the enzyme.

2-The rate continues to rise until reaches a peak at the optimum temperature




What is The effect of Temperature on the activity of Enzyme?

3-Above this temperature, the rate usually falls dramatically. This is because the increase energy cause bonds that maintain the enzymes shape to break , and the enzyme becomes **denatured**.(three dimensional structure will change)



- The changed shape means that the substrate can no longer fit into the active site, and enzyme activity is lost.



Method:

- This exercise illustrates the effect of increasing enzyme concentrations on reaction rate.
- You will perform a series of **5-minutes assays**, in which different temperature will be used

| In order to detect the effect of Temperature you must fix all the component except the temperature | |
|--|-----------------|
| Time (5 minutes) | constant |
| Enzyme concentration | constant |
| Temperature (degree) | Variable |
| pH (5.7) | constant |

Method:

- The effect of the following temperature will be studied on Acid phosphatase:

| Desired temperature (°C) | Method of preparation |
|--------------------------|-------------------------------------|
| 4°C | Ice plus tap water in an ice bucket |
| 25°C | Room temperature |
| 37°C | Thermos-tatted water bath |
| 65°C | Thermos-tatted water bath |
| 80°C | Thermos-tatted water bath |
| 100°C | Thermos-tatted water bath |

Method:

- Label 12 assay tubes as the following:

| Each should have its own Blank | |
|--------------------------------|-------------|
| Blank-4C° | Test-4 C° |
| Blank-25 C° | Test-25 C° |
| Blank -37 C° | Test-37 C° |
| Blank-65 C° | Test-65 C° |
| Blank-80 C° | Test-80 C° |
| Blank-100 C° | Test-100 C° |

Method:

- You must prepare the following for each tube (Blank + Test):

| Chemical | Volume |
|--|--------|
| 1.0M sodium acetate buffer (pH 5.7) | 0.5 ml |
| 0.1M MgCl₂ | 0.5 ml |
| p-nitrophenyl phosphate | 0.5 ml |
| Water | 5 ml |

- Place the tubes in the labeled temperature and let the temperature equilibrate for 5 min

Method

- For the blank you must first add 0.5 ml KOH then add the enzyme (we do not want any product to appear in the blank)
- For **TEST** , Add **0.5 ml of enzyme** extract to TEST and allow the reaction to proceed for **5 min**, using all the water bath temperatures described in the previous table
- Stop the reaction by the addition of **0.5ml of KOH** .
- **When all of the reaction mixtures have returned to room temperature**, determine the **absorbance at 405 nm** of each experimental tube against its own blank.

Results :

| Temperature | Absorbance 405 nm | Velocity (μ mole of PNP/min) |
|-------------|----------------------|--------------------------------------|
| 0 | | |
| 25 | | |
| 37 | | |
| 65 | | |
| 80 | | |
| 100 | | |

Plot a graph illustrating the effect of different temperatures on the rate of the reaction.

Calculations:

$$\text{Velocity (V)} = (A \times 10^6) / (E \times \text{time}) = \text{ }\mu\text{mole of PNP/min}$$

A= absorbance

E= extension coefficient= 18.8×10^3

Time = 5 min

Discussion

- Discuss the shape of the graph. And From the curve, explain and discuss the relationship between the activity of acid phosphatase and temperature.
- Define the optimum temperature and determine it from the curve.