



قسم الكيمياء الحيوية  
Biochemistry Department  
College of Science - King Saud University

## BCH 471


# **Determination of plasma enzymes**

Determination of LDH in serum

# Objectives

- To determine the level of LDH in serum.
- To evaluate the presence of tissue damage.

Most clinical enzyme measurements using serum or plasma, occasionally other fluids, such as urine and gut secretions, are investigated.




# Differences of Functional and Nonfunctional plasma enzymes

|                           | Functional plasma enzymes              | Nonfunctional plasma enzymes                              |
|---------------------------|--|---|
| <b>Their substrate</b>    | Always present in the blood            | Absent from the blood                                     |
| <b>Site of synthesis</b>  | Liver                                  | Different organs<br>e.g. liver, heart, muscles, and brain |
| <b>Effect of diseases</b> | Decrease in liver diseases             | Different enzymes increase in<br>different organ diseases |
| <b>Examples</b>           | Clotting factors<br>Lipoprotein Lipase | ALT<br>LDH<br>Acid Phosphatase<br>Amylase                 |

# Sources of Nonfunctional Plasma Enzyme

- Cell damage with the release of its content of enzymes into blood e.g. Myocardial infarction and viral hepatitis
- Obstruction of normal pathways e.g. Obstruction of bile duct increases alkaline phosphatase
- Increase of the enzyme synthesis e.g. bilirubin increases the rate of synthesis of alkaline phosphatase in obstructive liver disease
- Increased permeability of cell membrane as in hypoxia



# Medical Importance of Non Functional Plasma Enzymes

Measurement of non functional enzymes is important for:


- Diagnosis of diseases
- Prognosis of the disease: following up of the treatment by measuring plasma enzymes before and after treatment.

# Lactate Dehydrogenase (LDH)

- Lactic acid dehydrogenase (LDH) is an enzyme that helps produce energy.
- LDH is most often measured to evaluate the presence of tissue damage (diagnostic)
- The enzyme LDH is in many body tissues, especially the heart, liver, kidney, skeletal muscle, brain, blood cells, and lungs.

# LDH Reaction

- LDH is a hydrogen transfer enzyme which catalyzes the interconversion of pyruvate and lactate with the mediation of NAD as hydrogen acceptor, eventually converting pyruvate to glucose.
- The optimum pH for lactate pyruvate (L→P) reaction is 8.8 – 9.8
- While for pyruvate to lactate (P→L) is 7.7 – 7.8.
- The enzyme is inhibited by sulfhydryl reagents and mercuric ions.





# LDH Isoenzymes

- LDH exists in 5 forms (isoenzymes), which differ slightly in structure.
- All of these isoenzymes can be measured in the blood, and can be separated by electrophoresis.

| LDH isoenzyme | Tissues   |
|---------------|---|
| LDH-1         | is found primarily in heart muscle and red blood cells. |
| LDH-2         | is concentrated in white blood cells.                   |
| LDH-3         | is highest in the lung                                  |
| LDH-4         | is highest in the kidney, placenta, and pancreas        |
| LDH-5         | is highest in the liver and in skeletal muscle          |


**↑LDH**  
**in**  
**plasma**

| Diseases              | Examples             |
|-----------------------|----------------------|
| Myocardial infarction |                      |
| Liver Disease         | Toxic jaundice       |
|                       | Viral hepatitis      |
|                       | Obstructive jaundice |
| Anemia                | Pernicious anemia    |
|                       | Megaloblastic anemia |
| Renal Diseases        | Tubular necrosis     |
|                       | Pyelonephritis       |
| Malignant Disease     | Lung Cancer          |
|                       | Hodgkin's disease    |

# LDH Assay

## Principle:

LDH catalysis the following reaction:



The rate of NADH formation is indicated by **increase the absorbance** at 340nm and it is **directly proportional to serum LDH activity**.

If:

NADH is **product** : **increase** the absorbance /min

NADH is **reactant**: **decrease** the absorbance /min

# Method

|   | Tube  |
|---|-------|
| <b>LDH reagent</b>  | 1 ml  |
| Pre-warm at 37 °C for 3 minutes and add   |       |
| <b>Sample (serum)</b>   | 25 µl |
| Mix and incubate at 37 °C for 1 minutes, then read the absorbance at 340 nm against <u>distilled water (blank)</u> every minute for 3 minutes and determine $\Delta A/\text{min}$ . |       |

**2) Applications → 2) Simple Kinetics → wave length (340 nm) → 1) Seconds → Duration (180 sec = 3 min) → Intervals (60 sec = 1 min) → Print Data Table (off) → Press start (2 times)**

# Results

|    | Time (min) | Absorbance at 340 nm |
|----|------------|----------------------|
| A1 | 1          |                      |
| A2 | 2          |                      |
| A3 | 3          |                      |

# Calculations

$$\Delta A_1 = A_2 - A_1$$

$$\rightarrow \Delta A/\text{min} = (\Delta A_1 + \Delta A_2) / 2$$

$$\Delta A_2 = A_3 - A_2$$

$$\text{LDH (U/L)} = \Delta A \times 6592$$

- **Normal Values** 109 to 245 U/L