

Lab# 1


# Scanning spectrophotometry and spectrophotometric determination of concentration

---

BCH 333 [PRACTICAL]

## Objectives:

1. What is absorption spectrum and determination of  $\lambda_{\text{max}}$ .
2. Standard curve of known concentration and determine the concentration for a solution “with unknown concentration”.



How a spectrophotometer works.  
(spectrophotometer components)

**Spectrophotometer:** it can be used to measure the amount of light absorbed by a solution. How?

-By using the spectrophotometer, we can quantitatively measure the absorbance, and this information can be used to determine the concentration of the absorbing molecule.

-More concentrated solution will absorb more light and transmits less.[why?]



So, the more concentrated solution → high absorbance value.

And Less concentrated solution → less absorbance value.

# Absorption spectrum


## Absorption spectrum:

The curve that display the action or behaviour of absorption of molecule [solute] at different weave lengths, Known as Absorption spectrum curve.


-Every molecule has its own Absorption spectrum, So it's considered as fingerprint for each molecule.

**$\lambda_{max}$** : it is the weave length, at which the molecule has the maximum absorbance at this weave length.[best absorbance]. It can be determined from absorption spectrum curve.

[The wavelength of maximum absorption ( $\lambda_{max}$ )]



Absorption spectrum curve and determination of  $\lambda_{max}$ .



Absorption spectrum curve for different molecules.



## What is the benefit of studying the Absorption spectrum :

- Determine  $\lambda_{\text{max}}$ .
- Used to identify substances.
- It also can be used to know if there is any contamination with another molecule.

How does a spectrophotometer work?

<http://www.youtube.com/watch?v=pxC6F7bK8CU>



# **Standard curve of Concentrations**

# Principle Standard curve of Concentrations:

## Beer-Lambert law:

The absorption of light by a solution is described by the Beer-Lambert law as:

There is linear relationship between absorbance and concentration of an absorbing species.

$$A = \epsilon lc$$

A= is the absorbance.

$\epsilon$  = extinction(absorption) coefficient.

l = length of the light path through the solution.

c = concentration of the absorbing substance.

## So, what does standard curve of concentrations mean?

## A standard curve for concentrations:

- It is a graph that shows the relationship between different known concentrations of a substance and the absorbance at a specific wave length.

-Standard curve are most commonly used to determine the concentration of a substance, using serial dilution of solutions of known concentrations.

-So, what is the principle of a standard curve concentrations???


Wait..... what is standard solutions???

## Standard Solution:



is a solution containing a precisely known concentration of an element or a substance.

**How to determine unknown concentration of a solution with known absorbance value?**

If the unknown concentration of a solution has absorbance value =0.45, the conc. From the curve will be .....??



## To determine the concentration of a solution with “an unknown concentration”:



### From standard curve of concentration:

-Measure the absorbance of the “solution with unknown concentration” in order to determine the concentration, from the curve.

### Beer-Lambert law:

Using available information of any standard solution to determine the “ $\epsilon$ ”, then using these information to get the unknown concentration. Using:  $A = \epsilon lc$

Note: “ $\epsilon$ ” will changed when the weave length changed.








Principles and Techniques of  
**Biochemistry** and  
**Molecular Biology**

SEVENTH EDITION

Edited by  
**Keith Wilson**  
and  
**John Walker**



## **Spectrophotometry Introduction:**

[http://www.youtube.com/watch?v=qbCZbP6\\_j48](http://www.youtube.com/watch?v=qbCZbP6_j48)



## **Spectrophotometry Example**

[http://www.youtube.com/watch?v=VqAa\\_cmZ7OY&feature=relmfu](http://www.youtube.com/watch?v=VqAa_cmZ7OY&feature=relmfu)

