


Molecular Biology



جامعة
الملك سعود
King Saud University



Serial Dilution

Vortex mixer

Spectrophotometer



MIC.



MIC.

How to do serial of dilution?

A dilution definition is:

a process that reduces the concentration of a substance in a solution.

A serial dilution definition is:

the repeated dilution of a solution to amplify the dilution factor quickly

In which experiment do we often need dilutions?

- It's commonly performed in experiments requiring highly dilute solutions with great accuracy.
such as involving experiments to determine density of bacteria.

How to do serial dilution?


1- Determine the proper dilution liquid
distilled water or culture media

2- Prepare several test tubes with 9 mL of dilution liquid

- You will be adding your undiluted sample to the first tube and then serially diluting into the following tubes.
- label all of your tubes
- Each tube will be a 10-fold dilution starting from the undiluted tube
- The first tube will be a 1:10 dilution.....etc.

3- Prepare a test tube with at least 2 mL of your undiluted solution


- to perform this serial dilution is 1 mL of undiluted solution
- Thoroughly mix your solution before starting any dilutions



How to do serial dilution?

4- Perform the first dilution. Draw 1 mL of undiluted solution from test tube with a **micropipette** and transfer it to the test tube labeled $1:10$ containing 9 mL of the dilution liquid and mix thoroughly
The solution, therefore, has been diluted by a factor of 10

5- Perform the second dilution. For the second serial dilution, you will take 1 mL of solution from tube $1:10$ and add it to the 9 mL of dilution liquid in the tube $1:100$.
Thoroughly mix tube $1:10$ before adding to the next tube.



How to do serial dilution?

6- Calculating Final Dilution Factor and Concentration

•Concentration definition:

the amount of solute dissolved in a volume of solution.


- (1) volume of solution increases from V_{initial} to V_{final}
- (2) concentration of solution decreases from c_{initial} to c_{final}

•Dilution factor:

refers to the ratio of the volume of the initial (concentrated) solution to the volume of the final (dilute) solution, that is the ratio of V_1 to V_2 or, $V_1 : V_2$

•A dilution factor, DF, can be calculated:

$$DF = V_2 \div V_1$$



How to do serial dilution?

•Example:

- prepare a 1:50 dilution of the solution.
 - take a known volume of the stock solution (V_{initial}) and add enough solvent to it so that the solution has a new volume, V_{final} of $50 \times V_{\text{initial}}$.
 - $V_{\text{final}} = 50 \times V_{\text{initial}}$
 - The "1:50" tells you the dilution factor, the ratio of volumes, to use to prepare the new solution.
 - $V_1:V_2$
 - 1:50
 - In this case it tells us that $V_1 = 1$ and $V_2 = 50$
- so the dilution factor, $DF, = V_2 \div V_1 = 50 \div 1 = 50$
- That is, the new, diluted solution will have a volume 50 times greater than the volume of the original, undiluted, solution:

$$V_{\text{final}} = DF \times V_{\text{initial}}$$

How to do serial dilution?

The table below gives you a number of different options for preparing a 1:50 dilution:

V_{initial} (initial volume)	1 mL	1 L	0.1 mL	2 mL	25 μL
V_{final} (final volume)	50 mL	50 L	5 mL	100 mL	1250 μL
ratio of volumes used $V_{\text{initial}} : V_{\text{final}}$ $V_1 : V_2$	1 : 50 1 : 50	1 : 50 1 : 50	0.1 : 5 1 : 50	2 : 100 1 : 50	25 : 1250 1 : 50
dilution factor $\text{DF} = V_2 \div V_1$	$50 \div 1$ = 50	$50 \div 1$ = 50	$50 \div 1$ = 50	$50 \div 1$ = 50	$50 \div 1$ = 50

How to do serial dilution?

7- Determine the concentration of the solution following dilution

The equation is $C_{final} = C_{initial}/D$ where C_{final} is the ending concentration of the diluted solution, $C_{initial}$ is the starting concentration of the original solution and D is the dilution ratio previously determined


Example: If you started with a solution of cells with a concentration of 1,000,000 cells per mL and your dilution ratio is 1,000, what is the final concentration of your diluted sample?

Using the equation:

$$C_{final} = C_{initial}/D$$

$$C_{final} = 1,000,000/1,000$$

$$C_{final} = 1,000 \text{ cells per mL.}$$



Vortex Mixer

Vortex: - is a simple device used commonly in laboratories to mix small vials of liquid. As the motor runs the rubber piece oscillates rapidly in a circular motion.

- When a test tube or other appropriate container is pressed into the rubber cup (or touched to its edge) the motion is transmitted to the liquid inside and a vortex is created.
- **Vortex mixers** are quite commonplace in bioscience laboratories.
- In cell culture and microbiology laboratories they may be used to suspend cells. In a biochemical or analytical laboratory they may be used to mix the reagents of an assay or to mix an experimental sample and a dilutant.



Vortex mixer



Heated Microplate Vortexer

Spectrophotometer


A spectrophotometer is

- commonly used for the measurement of transmittance or reflectance of solutions, transparent or opaque solids, such as polished glass, or gases.
- to measure the diffusivity on any of the listed light ranges that usually cover around 200 nm - 2500 nm using different controls and calibrations.

Calibration in measurement technology and metrology is the comparison of measurement values delivered by a device under test with those of a calibration standard of known accuracy.



Roua AL Kufeidy



Spectrophotometer


Estimating Bacterial Numbers by Indirect Methods

A spectrophotometer is used to determine turbidity ("cloudiness") by measuring the amount of light that passed through a suspension of cells.

More cells = more turbidity; more turbidity = less light passing through the suspension

%T is percent transmission - fewer cells present (less turbidity) will allow more light to pass through, the %T is higher when the cell number is lower.

Absorbance is the opposite of %T. More light is absorbed when more cells are present - some people like this measure better because absorbance goes up as turbidity (or cell number) goes up.



sources

- 1- <http://abacus.bates.edu/~ganderso/biology/resources/dilutions.html>
- 2- <http://study.com/academy/lesson/serial-dilution-in-microbiology-calculation-method-technique.html>
- 3- <http://study.com/academy/lesson/serial-dilution-in-microbiology-calculation-method-technique.html>
- 4- <http://biology.kenyon.edu/courses/biol09/tetrahymena/serialdilution2.htm>
- 5- <http://scienceprimer.com/serial-dilution>
- 6- <http://abacus.bates.edu/~ganderso/biology/resources/dilutions.html>
- 7- <http://www.usetute.com.au/dfactor.html>
- 8- <http://www.usetute.com.au/dfactor.html>
- 9- J Buie (September 7, 2011 (acc. January 18, 2014)). *"Evolution of Biological Shakers and Stirrers."*
- 10- ISO 12647-2: Graphic technology — Process control for the production of halftone colour separations, proof and production prints — Part 2: Offset lithographic processes. Geneva: International Organization for Standardization. 2013. p. 13.
- 11- Allen, D., Cooksey, C., & Tsai, B. (2010, October 5). Spectrophotometry. Retrieved from
- 12- <http://www.nist.gov/pml/div685/grp03/spectrophotometry.cfm>
- 13- <http://www.usetute.com.au/dfactor.html>
- 14- <http://www.wikihow.com/Do-Serial-Dilutions>
- 15- <http://classes.midlandstech.edu/carterp/courses/bio225/chap06/Microbial%20Growth%20ss5.htm>
- 16- <https://en.wikipedia.org/wiki/Calibration>