
SPECTROPHOTOMETER I

LAB 3


DEFINITION

- A spectrophotometer is an instrument that measures the amount of light absorbed by a sample.
- Spectrophotometer techniques are used to measure the concentration of solutes in solution by measuring the amount of the light that is absorbed by the solution in a cuvette placed in the spectrophotometer

COMPONENTS OF SPECTROPHOTOMETER

COMPONENTS OF SPECTROPHOTOMETER

- **Light Source** : to provide a sufficient of light which is suitable for marking a measurement.
- 1- Tungsten Lamp
- 2- Hydrogen Lamp
- 3- Xenon Lam
- Tungsten Lamp :It is the most common light source used in spectrophotometer wavelength , range of about 330 to 900 nm

DISPERSION DEVICES

Dispersion devices

*Monochromator Accepts polychromatic input light from a lamp and outputs monochromatic light.

- Dispersion devices causes a different wavelength of light to be dispersion at different angles.
- Types of Dispersion devices used are
 - 1- Prism
 - 2- Filters
- Prism is used to isolate different wavelength.
- Filters separate different parts of the electromagnetic spectrum by absorbing or reflecting certain wavelengths and transmitting other wavelengths.

CUVETTES

- Absorption cells(Cuvettes) A cuvette is a kind of cell (usually a small square tube) sealed at one end, made of Plastic, glass or optical grade quartz and designed to hold samples for spectroscopic experiments


DETECTORS

- **Detectors** Any photosensitive device can be used as a detector of radiant energy .The photocell and phototube are the simplest photo detectors, producing current proportional to the intensity of the light striking them.

DISPLAY DEVICE

- **Display devices** The data from a detector are displayed by a readout device, such as an analog meter, a light beam reflected on a scale, or a digital display , Or liquid crystal display(LCD) .The output can also be transmitted to a computer.

COMPONENTS OF SPECTROPHOTOMETER

THE WORK OF SPECTROPHOTOMETER

- 1.A sample solution is placed inside the spectrophotometer.
- 2.A light source shines light toward the sample.
- 3.A monochromator splits the light into each color, or rather, individual wavelengths. An adjustable slit allows only one specific wavelength of light through to the sample solution.
- 4.The wavelength of light hits the sample, which is held in a little container called a cuvette.
- 5.Whatever light passes through the sample is read and displayed on the output screen.

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