




Triglyceride determination



Introduction:





- Triglycerides are esters of fatty acids and are hydrolyzed to glycerol and free fatty acids (by lipase)
- Triglyceride determinations when performed in conjunction with other lipid assays are useful in the diagnosis of **primary and secondary hyperlipoproteinemia**.
- They are also of interest in following the course of diabetes mellitus, nephrosis, biliary obstruction, and various metabolic abnormalities due to endocrine disturbances.



- **Hyperlipoproteinemia:** abnormally elevated of fat in blood (disorder in lipid metabolism).
- Standard methods for the measurement of triglyceride concentrations involved either **enzymatic** or **alkaline hydrolysis** to liberate glycerol.


- Principle:

The enzymatic reaction sequence employed in the assay of Triglycerides is as follows:



- The present procedure involves hydrolysis of triglycerides by lipase.
- The glycerol concentration is then determined by enzymatic assay coupled with Trinder reaction that terminates in the formation of a quinoneimine dye.
- **The amount of the dye formed**, determined by its absorption at 505 nm, is directly proportional to the concentration of triglycerides in the samples.

- Specimen collection and storage:

1. Fresh, non-hemolyzed serum from fasting patients is recommended.
 2. Triglycerides in serum appears stable for three days when stored at 2-8 °C.
 3. Prolonged storage of the samples at room temperature is not recommended since other glycerol containing compounds may hydrolyze, releasing free glycerol with an apparent increase in total triglycerides content.
- 

- Method :

- By Triglyceride reagent kit.

-Follow the table:

	Blank	Standard	Test
Reconstituted Reagent	1 ml	1 ml	1 ml
Pre-worm at 37°C for 2 min and add:			
Standard	---	0.01 ml (10 µl)	---
Sample	---	---	0.01 ml (10 µl)
Mix and incubate at 37°C for 10 min			
↓			
Read the absorbance of standard and sample at 505 nm against blank			

-Calculation:

$$\text{Conc. of TG} = \frac{\text{Ab Test}}{\text{Ab Std.}} \times \text{conc. of Std. (200 mg/dl)}$$

- Normal range: 10 -190 mg/dl





HDL-Cholesterol determination



- Introduction:

- Cholesterol is a fatty substance found in blood, bile and brain tissue.
- It serves as a precursor to **bile acids, steroids and vitamin D**.
- **In the plasma, cholesterol is transported by three lipoproteins:** high density lipoprotein (HDL-Cholesterol), low density lipoprotein (LDL-Cholesterol), and very low density lipoprotein (VLDL- Cholesterol).
- The concentration of **total cholesterol** in serum has been associated with metabolic, infectious and coronary heart diseases.



- The concentration of HDL-cholesterol in serum has important in diagnosis of the how the level of **risk to get coronary heart diseases**.
 - (More HDL-Chol. That indicate low risk to get coronary heart disease.)
 - Castelli and co-workers have indicated that an **inverse relationship** exists between serum HDL-Cholesterol and the risk of coronary heart disease.
 - The measurement of **HDL Cholesterol and triglyceride** provides valuable information for the prediction of coronary heart disease and for lipoprotein phenotyping.
- 



- Specimen collection:

1. Specimen should be serum and free from hemolysis.
2. Patient should be fasting for 12-14 hours.

- Principle:

- HDL cholesterol determination

- Enzymatic methods, involving cholesterol esterase and oxidase and Trinders color system.
- **The enzymatic reaction sequence employed in the assay of cholesterol is as follows:**



- Cholesterol Esters are hydrolyzed to produce cholesterol, Hydrogen peroxide is then produced from the oxidation of cholesterol by cholesterol oxidase. In a coupled reaction catalyzed by peroxidase, quinoneimine red colored dye is formed from 4-aminoantipyrine, phenol and hydrogen peroxide. The absorption of light at 505 ± 5 nm of the solution of this dye is proportional to the concentration of cholesterol in the sample.

Method :

- HDL Cholesterol:

- Follow the Table:

	Blank	Calibrator	Assay
Reagent R1	300 µl	300 µl	300 µl
Calibrator		3 µl	3 µl
Mix vigorously, let stand for 5 min at 37°C . Read absorbance A1 at 600 nm against blank.			
Add Reagent R2	100 µl	100 µl	100 µl
Mix vigorously, let stand for 5 min at 37°C . Read absorbance A2 at 600 nm against blank.			

- Calculation :

* Determine the HDL Cholesterol conc.

$$\Delta \text{ Abs.} = (A2 - 0.75 A1)$$

$$\text{Conc. of HDL} = \frac{\Delta \text{ Ab Assay}}{\Delta \text{ Ab Calibrator}} \times \text{conc. of calibrator (50 mg/dl)} = \text{mg/dl}$$

- Normal value of :

- HDL-Cholesterol :

- Low level (risk factor) $< 40 \text{ mg/dl}$
- High HDL (protector factor) $\geq 60 \text{ mg/dl}$