




**BCH 447**

# **Triglyceride Determination in Serum**

# Introduction:

- Triglycerides are esters of fatty acids and are hydrolyzed by lipase to glycerol and free fatty acids.







- Triglyceride determinations when performed in conjunction with other lipid assays are useful in the diagnosis of **primary and secondary hyperlipidemia** (abnormally elevated of fat in blood )



# Principle:

- Standard methods for the measurement of TG concentrations involved either enzymatic or alkaline hydrolysis to liberate glycerol.
- The glycerol concentration is then determined by enzymatic assay coupled with Trinder reaction that terminates in the formation of a quinoneimine dye.
- **The amount of the dye formed is** determined by its absorption at 505 nm, it is directly proportional to the concentration of triglycerides in the samples.



Trinder reaction: It is the reaction between hydrogen peroxide and the phenol and aminoantipyrine (AAP) to form quinoneimine (red-violet dye), catalyzed by the presence of a peroxidase

## Specimen collection and storage:

1. Fresh, non-hemolyzed serum from fasting patients is recommended.
2. Triglycerides in serum appears stable for three days when stored at 2-8 °C.
3. Prolonged storage of the samples at room temperature is not recommended since **other glycerol containing compounds may hydrolyze**, releasing free glycerol with an apparent increase in total triglycerides content.

## Method :

- By Triglyceride reagent kit.


**-Follow the table:**

	Blank	Standard	Test
Reconstituted Reagent	1 ml	1 ml	1 ml
Pre-warm at 37°C for 2 min and add:			
Standard	---	0.01 ml (10 µl)	---
Sample	---	---	0.01 ml (10 µl)
Mix and incubate at 37°C for 10 min			
↓			
Read the absorbance of standard and sample at 505 nm against blank			

# Calculation:

$$\text{Conc. of TG} = \frac{\text{Ab Test}}{\text{Ab Std.}} \times \text{conc. of Std. ( 200 mg/dl)}$$

**Normal range:** 10 -190 mg/dl




TG → **Hyperthyroidism**

**Malnutrition**

**Low-fat diet**

**Malabsorption**



TG → **Hyperlipidemia**

**Hypothyroidism**

**Obesity**

**Kidney Disease**

**Diabetes**




# **HDL-Cholesterol determination**



# Introduction:


- Cholesterol is a fatty substance found in blood, bile and brain tissue.
- It serves as a precursor to **bile acids, steroids and vitamin D**.



## •In the plasma, cholesterol is transported by three lipoproteins:




- high density lipoprotein (HDL-Cholesterol),
  - low density lipoprotein (LDL-Cholesterol),
  - and very low density lipoprotein (VLDL- Cholesterol).
- 
- The concentration of **total cholesterol** in serum has been associated with metabolic, infectious and coronary heart diseases.



- The concentration of HDL-cholesterol in serum has important in diagnosis of the how the level of **risk to get coronary heart diseases**.
  - There is an **inverse relationship** between serum HDL-Cholesterol and the risk of coronary heart disease
    - **More HDL-Chol** → indicate **low risk** of coronary heart disease.
  - The measurement of **HDL Cholesterol and triglyceride** provides valuable information for the prediction of coronary heart disease and for lipoprotein phenotyping.
- 

# Principle:

- Enzymatic methods, involving cholesterol esterase and oxidase and Trinders color system.
- The enzymatic reaction sequence employed in the assay of cholesterol is as follows:



- The amount of the dye formed is** determined by its absorption at  $505 \pm 5 \text{ nm}$ , it is directly proportional to the concentration of cholesterol in the samples.

## Method :

- By HDL-Cholesterol reagent kit.

**-Follow the table:**

	Blank	Calibrator	Assay
Reagent R1	300 $\mu$ l	300 $\mu$ l	300 $\mu$ l
Calibrator /Sample		3 $\mu$ l	3 $\mu$ l
Mix vigorously, let stand for 5 min at 37°C . Read absorbance A1 at 600 nm against blank.			
Add Reagent R2	100 $\mu$ l	100 $\mu$ l	100 $\mu$ l
Mix vigorously, let stand for 5 min at 37°C . Read absorbance A2 at 600 nm against blank.			

## Calculation :

$$\Delta \text{ Abs.} = (A2 - 0.75 A1)$$

$$\text{Conc. of HDL} = \frac{\Delta \text{ Ab Assay}}{\Delta \text{ Ab Calibrator}} \times \text{conc. of calibrator (50 mg/dl)} = \text{mg/dl}$$

### - Normal value of HDL-Cholesterol :

- Low level (risk factor) < 40 mg/dl
- High HDL (protector factor)  $\geq$  60 mg/dl