

CLS 281

Basic Biochemistry and Biomolecules


جامعة  
الملك سعود  
King Saud University




Experiment 6

# Color Tests for Specific Carbohydrates: Ketoses, Pentoses, and Polysaccharides.


# Carbohydrate Classification




# Disaccharide Examples



Sucrose  
glucose + fructose




Lactose  
galactose + glucose



Maltose  
glucose + glucose

# Aldose vs Ketose




**Aldose** has a **carbonyl group** at the end of the carbon chain.  
**Ketose** has a **carbonyl group** in the middle of the carbon chain.

# Classification of Important Sugars


	<b>Aldoses</b>	<b>Ketoses</b>
Trioses ( $C_3H_6O_3$ )	Glycerose (glyceraldehyde)	Dihydroxyacetone
Tetroses ( $C_4H_8O_4$ )	Erythrose	Erythrulose
Pentoses ( $C_5H_{10}O_5$ )	Ribose	Ribulose
Hexoses ( $C_6H_{12}O_6$ )	Glucose, galactose, mannose	Fructose
Heptoses ( $C_7H_{14}O_7$ )	—	Sedoheptulose

# Most Common Polysaccharides


Starch





Glycogen



Cellulose



# Dehydration of Monosaccharides into Furfural



# Today's experiments

1. Seliwanoff's Resorcinol Test.
2. Bial's Orcinol Test.
3. Iodine Test.




# 01 Seliwanoff's Resorcinol Test

- **Aim**
  - This test is used to distinguish between ketoses and aldoses monosaccharides.
- This test gives **fast**, clear **positive** results with ketohexoses. Thus, it's a specific test for ketohexoses.
- **Why ketohexoses?**
  - **The Ketohexoses** are about 20 to 25 times faster than aldohexoses in producing furfural derivatives. And form considerably more furfural derivatives.
- If you end the reaction within 1 minute, you get a positive result with the fastest sugar (ketohexoses).
- **Note:**
  - Free or bounded ketohexoses can respond to the test. e.g. sucrose that contains fructose.

## 01

# Seliwanoff's Resorcinol Test Principle

- Seliwanoff's Resorcinol Reagent
  - HCL
  - Distilled Water (D.W)
  - Phenolic compound (Resorcinol)
- Principle



▲  
Ketohehexoses Hexoses + HCL → Hydroxymethyl furfural + 3 H<sub>2</sub>O

▲  
Hydroxymethyl furfural + Resorcinol → **Red product** + 3 H<sub>2</sub>O

This aldehyde undergoes **condensation** along with Resorcinol.

▲ = Heat

01

# Seliwanoff's Resorcinol Test Procedure


Note: swirl the samples and reagent bottles before use to aspirate homogenous solution.

Steps	Tube No.	Tube 1	Tube 2	Tube 3	Tube 4	Tube 5	Tube 6
1	Sample	1% Fructose	1% Glucose	1% Sucrose	1% Sorbose	1% Xylose	D.W
	Volume	2 drops	2 drops	2 drops	2 drops	2 drops	2 drops
2	Resorcinol Reagent	5 ml	5 ml	5 ml	5 ml	5 ml	5 ml
3	<b>Mix</b> and Incubate in a boiling water bath for 1 min. Observe the color >>> record your result as (result after 1 min)						
4	Continue incubation for 4 min. Observe the change in color. >>> record your result as (result after 4 min)						
	Result 1 min						
	Result 4 min						

## 01

# Seliwanoff's Resorcinol Test Result


- **Positive Result**
  - Ketohexoses → **red complex**
- **Interference**
  1. Aldohexoses (glucose) → **light yellow to faintly pink color**
  2. Pentose → **blue to green color**
- **How to solve this interference?**
  - Make the concentration of **HCL** less than 12%.
  - Make the **incubation** period shorter than 1 minute.
  - If **aldohexose**, e.g., glucose, is present, it must not be in a concentration greater than 2%.




# 02 Bial's Orcinol Test

- **Aim**
  - It is a simple, rapid qualitative test for pentoses.
- **Usage:**
  - It can be used For quantitative assay of pentoses (e.g., **Ribonucleic acid**) in the absence of interfering substances.
- **Bial's Reagent contains:**


- HCL
- Phenolic compound (Orcinol)
- Ferric chloride  
(used to increase the sensitivity of the test).



## D-Aldoses



## D-Ketoses



## 02

# Bial's Orcinol Test Principle

- Principle



Carbohydrate undergoes **dehydration** upon the introduction of concentrated Hydrochloric acid, resulting in the formation of **furfural**.



**Furfural** undergoes **condensation** along with orcinol.

## 02

# Bial's Orcinol Test Procedure


Note: swirl the samples and reagent bottles before use to aspirate homogenous solution.

Steps	Tube No.	Tube 1	Tube 2	Tube 3	Tube 4	Tube 5	Tube 6
1	Sample	1% Xylose	1% Glucose	1% Fructose	1% Lactose	1% Starch	D.W
	Volume	1 drop	1 drop	1 drop	1 drop	1 drop	1 drop
2	Bial's Reagent	3 ml	3 ml	3 ml	3 ml	3 ml	3 ml
3	Incubate in a boiling water bath for 3-5 mins. Observe the change in color formed.						

# 02

## Bial's Orcinol Test Result


- **Result**
  - **Blue-green color product**
- **Interference in this test:**
  - Hexoses generally react to form **green**, **red**, or **brown products**.
  - However, all of these hexoses colors are considered negative results.






# 03 Iodine Test for Polysaccharide

- **Aim**
  - The iodine test is used for the detection of polysaccharides.
- **Reagent**
  - Iodine
- **Principle**
  - Iodine forms colored adsorption complexes with polysaccharides.
  - Color is due to the coordination complex between the helically coiled polysaccharide chains and the iodine centrally located within the helix.




(b)



(a)


# 03 The Starch-Iodine Complex

- In the presence of **starch**, the iodine will fit into the center of the coiled polysaccharides chain and form a chain of 6 iodine molecules ( $I_6$ ) inside the helix, and the color turns into an intense **blue complex**.
- A similar complex is formed with other polysaccharides but will show a different color.



- (a) Schematic structure of the starch-iodine complex. The amylose chain forms a helix around  $I_6$  units.
- (b) View down the starch helix, showing iodine inside the helix.

(b)



(a)

## 03

# Iodine Test Procedure

Note: swirl the samples and reagent bottles before use to aspirate homogenous solution.


Steps	Tube No.	Tube 1	Tube 2	Tube 3	Tube 4	Tube 5
1	Sample	1% Starch	1% Glycogen	1% Dextrin	1% Cellulose	D.W
2	Volume	2 ml	2 ml	2 ml	2 ml	2 ml
3	Iodine Reagent	2 drops	2 drops	2 drops	2 drops	2 drops
4	Observe the color and record the result.					

# 03 Iodine Test Result

● Starch + Iodine → **Blue to black color**

● Dextrin + Iodine → **Red to violet color**

● Glycogen + Iodine → **Red to brown color**



# Summary of Color Tests for Specific Carbohydrates

Test	Detect	Reagent	Principle	Positive Result	Negative Result	Interference	Note
<b>Seliwanoff's Resorcinol Test</b>	Specifically for the detection of <u>Ketohexoses</u>	<b>Seliwanoff's Resorcinol Reagent:</b> HCL, Distilled Water (D.W), Phenolic compound (Resorcinol).	<p>Hexoses + <u>HCL</u> → <u>Hydroxymethyl furfural</u> + 3 H<sub>2</sub>O</p> <p>Hydroxymethyl furfural + <u>Resorcinol</u> → <b>Red product</b> + 3 H<sub>2</sub>O</p>	<b>Red complex</b>	Others	<p>1- Aldohexoses (glucose) → <b>light yellow to faintly pink color</b></p> <p>2-Pentose → <b>blue to green color.</b></p>	<p>This test is used to distinguish between <u>ketoses</u> and <u>aldoses</u> monosaccharides.</p> <p>This test gives <b>fast</b>, clear <b>positive</b> results with <u>ketohexoses</u>.</p>
<b>Bial's Orcinol Test</b>	<b>Pentoses</b>	<b>Bial's Reagent contains:</b> HCL, Phenolic compound (Orcinol), Ferric chloride (used to increase the sensitivity of the test).	<p>Pentose + <u>HCL</u> → furfural + 3 H<sub>2</sub>O</p> <p>furfural + <u>orcinol</u> → <b>Blue-green color product</b> + 3 H<sub>2</sub>O</p>	<b>Blue-green color product</b>	Others	<p><u>Hexoses</u> generally react to form <b>green, red, or brown products.</b> However, all are considered negative.</p>	<p>It can be used For <u>quantitative assay of pentoses</u> (e.g., <b>Ribonucleic acid</b>) in the absence of interfering substances</p>

# Summary of Color Tests for Specific Carbohydrates

Test	Detect	Reagent	Principle	Positive Result	Negative Result
Iodine Test	Polysaccharide	Iodine	<ul style="list-style-type: none"><li>• Iodine <u>forms colored adsorption complexes</u> with polysaccharides.</li><li>• Color is due to the <u>coordination complex</u> between the helically coiled polysaccharide chains and the iodine centrally located within the helix.</li></ul>	<ul style="list-style-type: none"><li>• Starch → <b>Blue to black color</b></li><li>• Dextrin → <b>Red to violet color</b></li><li>• Glycogen → <b>Red to brown color</b></li></ul>	Others

# Guideline for writing the lab report

## Total: 5 marks

All the following information should be included in your report:

- a) Course # (CLS 281)
- b) Experiment title
- c) Date of the experiment
- d) Student's names and university ID#
- e) Section #

The lab report is broken down into 6 sections:

1. Experiment **title**
2. The **aim** of the experiment (objective, or what the test detects specifically) (1 mark)
3. **Principle** (chemical reaction) (1 mark)
4. **Methodology** (written in **steps**, **not in tables**)
5. **Result** (1 mark)
6. **Interpretation or Comment** (2 mark)

Deadline: Next lab    Submission: via email