

CLS 281

Basic Biochemistry and Biomolecules

جامعة  
الملك سعود  
King Saud University




Experiment 7

# Determination of Reducing Sugars by Somogyi-Nelson Method

# Reducing Sugars

- A reducing sugar is a sugar with a free or potentially free aldehydic or ketonic group.




- Oxidation is the loss of electrons
- Reduction is the gain of electrons

# Review of Reducing Sugars Test


- *Benedict's Test Result*

Color of the Precipitate	% of Reducing Sugar
Green	0.5%
Yellow	1%
Orange	1.5%
Red	2% or more



***Semi-quantitative***

- *Barfoed's Test Result*



Blue Solution

Carbohydrates absent

Red Precipitation

Within few minutes - monosaccharides  
After 3 minutes- disaccharides


***Qualitative***

# Type of Testing

- **Qualitative** examinations measure the presence or absence of a substance.
- **Semi-quantitative** examinations provide an estimate (e.g. %) of how much of the measured substance is present.
- **Quantitative** examinations are used for determining the amount of an analyte in a sample. The amount is always expressed as a number with appropriate units.


# Somogyi-Nelson Method

- The Nelson-Somogyi method is one of the classical and widely used methods for the **quantitative** determination of reducing sugars.




- The reducing sugars, when heated with alkaline copper tartrate, reduce the copper from the cupric to the cuprous state, and thus cuprous oxide is formed.
- When cuprous oxide is treated with Phosphomolybdic acid, the reduction of molybdic acid to molybdenum **blue** takes place.
- How do we analyze this test quantitatively and measure the concentration of glucose?** Using a spectrophotometer to measure the concentration of MO+4 at 520 nm.

# The Electromagnetic Spectrum



# Spectrophotometer



- If the light is sent through a solution ( $I_o$ ) and only some of the light passes through ( $I_t$ ), this suggests that the rest of the light has been absorbed by the molecules in the solution.

# Converting light intensity to concentration


- To convert light intensity to concentration the Beer-Lambert Law is usually applied.
- It is commonly applied to determine the concentration of various molecules in solution.

## Beer's law

- States that the **absorbance is directly proportional to the concentration** of a solution.
- If you plot absorbance versus concentration, the resulting graph yields **a straight line**.

$$A = \epsilon b C$$

**A** → Absorbance  
 **$\epsilon$**  → Molar absorptivity → L/(mol cm)  
**b** → Path length → cm  
**C** → Concentration Mol/L



$$C_1 \times V_1 = C_2 \times V_2$$

(Start reading) (final reading)



# Procedure steps

Step 1: Construction of glucose standard curve:


Use Beer's law to measure standards of glucose.

Step 2: Measure the unknown glucose sample:

Calculate the concentration using data obtained from the standard curve.

# Materials

- Spectrophotometer
- Cuvettes
- micropipette
- Tubes
- Rack
- Water bath




# Preparation of a Standard Curve

- A standard curve is generated to be used to determine the concentration of an unknown sample.

## How to make a standard curve?

1. Create a series of solutions (4-7 tubes) in known increasing concentrations of your analyte of interest (in this case, glucose of known concentration)
2. Use the formula  $C_1V_1 = C_2V_2$  to measure the concentration of each standard tube.
3. Measure the absorbance of all your standard at 520 nm ( $\lambda_{\text{max}}$ ) for glucose).
4. Used the values of concentration and absorbance to create a standard curve.

Now you can move to **step 2** and use the standard curve to determine the concentration of glucose in an unknown sample.



# Step 1 Construction of Glucose std. Curve

Note: **Accurate pipetting** is very important in performing your standard curve.

Total Volume (V<sub>2</sub>) = 1 ml

Steps	Tube No.	Blank	Std.1	Std.2	Std.3	Std.4	Std.5	Unknown sample
1	Volume of Glucose ([C1] 0.2mg/ml) (ml) [V1]	-	0.2	0.4	0.6	0.8	1	-
2	D.W (ml)	1	0.8	0.6	0.4	0.2	-	-
	Unknown sample (ml)	-	-	-	-	-	-	1
3	Copper reagent (ml)	1	1	1	1	1	1	1
4	Mix and incubate in a boiling water bath for 20 mins and cool.							
5	Arsenomolybdate reagent (ml)	1	1	1	1	1	1	1
6	Incubation Room Temperature for 1 minute.							
7	D.W (ml)	10	10	10	10	10	10	10
8	<ul style="list-style-type: none"> <li>• Add each into a cuvette.</li> <li>• Set the spectrophotometer and read the absorbance of standards against the blank at wavelength 520 nm.</li> </ul>							

# Step 1

## Standard Curve Calculation

1- Calculate the final concentration C2 for each standard using  $C1 \times V1 = C2 \times V2$  :

- C1 = Concentration of standard glucose = 0.2mg/ml
- V2=total volume of glucose + D.W dilution = 1ml in all tubes
- V1 = volume of std. Glucose was added to each tube.
- C2 = ? → The final concentration of standards. Glucose after dilution

2- Create The Result Table.

Read Abs. at 520 nm.

Tube No.	Blank	Std.1	Std.2	Std.3	Std.4	Std.5
Glucose (0.2m g/ml) V1	-	0.2	0.4	0.6	0.8	1
D.W	1	0.8	0.6	0.4	0.2	-


Record your result here	Abs	C2
Blank		
Std.1		
Std.2		
Std.3		
Std.4		
Std.5		
Unknown		

# Step 1


## Standard Curve plot

3- Draw The Curve of absorbance against concentration

- Standard curves are graphs of light absorbance versus solution concentration which can be used to figure out the solute concentration in unknown samples.



## Step 2 *Finding the concentration of unknown*



# Definitions

- **Standards:** a series of tubes containing increasing concentrations of the substance to be assayed (glucose) are treated with the reagent.
- **Blank:** a blank is prepared to contain all reagents except the compound to be measured (glucose), and then Absorbance is measured for each of the standard tubes against the blank.
- **Unknown sample:** the sample to be measured of unknown concentration of the analyte under investigation.



**FYI**

# Using Excel to plot standard curve

- [https://www.youtube.com/watch?v=KeTFzJUG\\_rA](https://www.youtube.com/watch?v=KeTFzJUG_rA)

# Note

- Reoxidation of cuprous ions by oxygen from the air is prevented by adding Sodium Sulphate in the reagent to decrease the solubility of oxygen.
- The Somogyi-Nelson method has been replaced recently using rapid colorimetric procedures like O-toluidine or enzymes (such as glucose oxidase and peroxidase).

# Report Criteria **Total: 5 marks**

- Course # (CLS 281), Experiment title, Date of the experiment, Student's names and university ID#, Section #, Experiment title.
- The **aim** of the experiment (objective, or what the test detects specifically) (1 mark).
- **Principle** (chemical reaction) (1 mark).
- **Result (2 mark):**
  - Table of standard concentration and absorbance value, write detailed calculation with units.
  - Standard curve – do not forget the units.
  - The absorbance of the unknown sample.
  - The concentration of the unknown sample + explains how you obtained the concentration value.
- **Interpretation or Comment (1 mark)** - Describe the observed relationship between the absorbance and the concentration.

Deadline: Next lab

Submission: Handout