

Lab sheet #2**-Estimation of reducing sugars by Dinitrosalicylic acid method-****Method:**

1- In clean test tubes add:

	Glucose solution (100mg/dl) (ml)	Sample (ml)	Water (ml)	DNS reagent (ml)		Sodium potassium tartrate (ml)
Blank	--	--	1	3	Cover the tubes (with aluminum foil) And heat for 5 min in a boiling water bath	1
1	0.1	--	0.9	3		1
2	0.2	--	0.8	3		1
3	0.3	--	0.7	3		1
4	0.4	--	0.6	3		1
5	0.5	--	0.5	3		1
6	0.6	--	0.4	3		1
7	0.7	--	0.3	3		1
8	0.8	--	0.2	3		1
Honey	--	1	---	3		1
	--	0.6	0.4	3		1
Milk	--	1	---	3		1
	--	0.6	0.4	3		1

- Mix the contents.
- Cool by immersing in cold water and read at 540 nm.
- Plot the standard curve, calculate the amount in the sample from standard curve and calculate the contents.

Results:

Tube	Absorbance at 540 nm	CHO content (mg/dl) ($C_1 \times V_1 = C_2 \times V_2$)
B	--	--
1		Ex: $C_1 \times V_1 = C_2 \times V_2$ $100 \times 0.1 = C_2 \times 1 \rightarrow C_2 = 10$
2		
3		
4		
5		
6		
7		
8		
Honey		
Milk		

Calculations:

- Dilution factor (DF) = final volume / aliquot volume
.....
.....
- Concentration of sample A = Concentration obtained from standard curve x 100
= A x DF (if used)
- The amount of carbohydrate in the sample=mg/dl x DF
.....
.....