




## BCH 447

# Oral Glucose Tolerance Test (OGTT)




## *Objectives:*

- Use OGTT in diagnosis of diabetes mellitus.



# Introduction

- There are two types of glucose tolerance test (Oral and IV)
- The most common glucose tolerance test is **the oral glucose tolerance test (OGTT)**.
- The test reveals how quickly glucose is metabolized from the bloodstream for use by cells as energy source.
- Serial measurement of plasma glucose before and after glucose is given orally should provide a standard method to evaluate individuals and establish values for normal and disease states.



# How does stress, caffeine and steroids affect OGTT ?

## Stress:

When the body is under stress, the adrenal glands trigger the release of glucose stored in various organs, which often leads to elevated levels of glucose in the bloodstream.

## Caffeine:

caffeine intake can acutely lower insulin sensitivity and increase glucose concentrations

## Corticosteroid:

These drugs have also been called "glucocorticoids" because of their effects on glucose metabolism: Increases in blood glucose are common among people taking prednisone and other steroids.

## How is the test performed?

- When an OGTT is ordered, the following conditions should be met:

(1) Omit medications known to affect glucose tolerance.

(2) Perform the test in the morning after 3 days of unrestricted diet and activity.

(3) Perform the test after a 10-16 hours fast (**12 hour is best**).

- Oral dose : For adults, the recommended load is 75g and for children, 1.75 g/kg,

- Plasma glucose should be measured **fasting** then every 30 min for 2h **after** an oral glucose load

- **Note:** the time of collection is different, it is depend on the situation.

## This test is recommended for :

- Generally most healthcare providers recommend that all pregnant women should be screened for ***gestational diabetes***.
- Experts recommend this test to pregnant women who are between 24 and 28 weeks of pregnancy .
- This test is also recommended for anyone suspected of developing adult **diabetes**.

## **OGTT side effect:**

- Some people feel sweaty, light-headed, or may even feel short of breath or faint after drinking the glucose.
- However, serious side effects of this test are very uncommon.



## Normal and abnormal results :

### Normal :

Fasting: 60 -128 mg/dL      1 hour: less than 200 mg/dL      2 hours: less than 140mg/dL

-Higher-than-normal levels of glucose may mean you have **prediabetes** , **diabetes (type 2)**, or **gestational diabetes**.


*Or impaired glucose tolerance.*

It means you are at increased risk for developing diabetes.



## Principle:

- Several reagents can be used to assay reducing sugars such as 3, 5 dinitrosalicylic acid in one of the compounds.
- In alkaline solution it is reduced to 3-amino-5-nitro salicylic acid, which is orange-red.
- Absorbance is determined at 540 nm.



## Method:

	Plasma	Standard	dH2O	DNS reagent
Test (a1) (Fasting plasma)	0.1	-	-	2 ml
Test (a2) (Fasting plasma)	0.1	-	-	2 ml
Test (b1) Two- hour	0.1	-	-	2 ml
Test (b2) Two- hour	0.1	-	-	2 ml
Standard (1)	-	0.1	-	2 ml
Standard (2)	-	0.1	-	2 ml
Blank	-	-	0.1	2 ml

↓  
Mix the contents of each tube and cover each tube by Aluminum foil

↓  
Boiling water bath for 5 minutes

↓  
cool the tubes for 1-3 min

↓  
Read absorbance at 540 nm

## - RESULT:

<b>Tubes</b>	<b>Absorbance at 540 nm</b>
Test (a1)	
Test (a2)	
Test (b1)	
Test (b2)	
Standard (1)	
Standard (2)	

# Calculations:

- Conc. Of Std. = 0.1 g/dl.

- Sample A = Fasting plasma glucose

- Sample B = Two hour plasma glucose

- Amount of glucose in plasma =  $\frac{\text{Mean Abs Test}}{\text{Mean Abs Std.}} \times \text{conc. of Std} = Z \text{ g/dl}$

-  $Z \text{ g/dl} \times 1000 = \underline{Y \text{ mg/dl}}$

Calculate the glucose in fasting and two hrs plasma glucose ..

Then discuss your results ..