Photometric Law of Distance

1 Objective

- The luminous intensity *I* emitted by a punctual source is determined as a function of distance from the source.
- The photometric law of distance is verified by plotting illuminance as a function of the reciprocal value of the square of the distance.

2 Prelab Questions

- 1. What is the difference between radiometry and photometry?
- 2. What is a *solid angle*? State the solid angle of a sphere.
- 3. State the Photometric Law of Distance and explain it.
- 4. Define the following terms: luminous flux, luminous intensity, illuminance, luminance and luminous existance, stating their respective laws and SI units.

3 Principles

The luminous intensity emitted by a punctual source is determined as a function of distance.

4 Apparatus

- Luxmeter.
- Luxmeter probe.
- Optical profile bench.
- Several slide mounts.
- Filament lamp 6 V/5 A.
- Universal clamp.
- Power supply.

5 Precautions

- 1. The input voltage for the lamp is 6 V.
- 2. The luxmeter must be calibrated before actually carrying out measurements. Possible background luminance must be determined with the lamp switched off and must be taken into account during evaluation.
- 3. Notice that the luxmeter probe is offset from the centre of the clamp.

6 Experimental Steps

- 1. Switch on the luxmeter and measure the background illuminance E_0 .
- 2. Set the separation r between the probe and lamp to 55 cm.
- 3. Switch on the lamp and measure the illuminance E.
- 4. Decreasing the separation in steps of 5cm, record E until you reach a separation of r = 25 cm.

7 Evaluation

- 1. Calculate r^2 and $\frac{1}{r^2}$.
- 2. Plot E vs. r^2 , confirming the Photometric Law of Distance.
- 3. Plot E vs. $\frac{1}{r^2}$ and find the slope.
- 4. Using the slope, calculate the luminous intensity I of the lamp, knowing that:

$$E = \frac{I}{r^2} \tag{1}$$

8 Postlab Questions

- 1. What is the difference between photopic and scotopic vision?
- 2. A surface oriented perpendicularly is positioned 2 m away from a lightbulb emitting 100 W of radiant flux.
 - (a) Calculate the irradiance E_e at the surface.
 - (b) If all 100 W is emitted from a red bulb at $\lambda = 650$ nm, calculate the illuminance E_v at the surface.

Hint: use the following equation: $photometricunit = K(\lambda) \times radiometricunit$. Where $K(\lambda)$ is the luminous efficacy given by: $K(\lambda) = 685V(\lambda)$ and $V(\lambda)$ is the luminous efficiency obtained from a CIE luminous efficiency curve.