

# ***Biochemistry 2 Booklet***

## ***CLS 331***

***Prepared by :***

***M.s Shaikhah Al- Subie***

***2019-1440***

## Chemical list biochemistry2

### CLS 331

Name	Unt	EXPIRATION DATE
Phosphate buffer (ph:6.7)	gm	
Starch	gm	
Iodine	Gm	
Sodium chloride	gm	
Phosphate buffer(ph:6.7)	Gm	
Starch	Gm	
Sodium chloride	Gm	
Sodium hydroxid	Gm	
maltose	Gm	
3,5Dinitrosalicylic acid	gm	
Sodium potassium tartarate	Gm	
Disodium hydrogen phosphate	Gm	
Potassium dihydrogen phosphate	Gm	
Yeast suspension	Gm	
Glucose	Gm	
Sodium nitroprusside	Gm	
Ammonium hydroxide	Gm	
Ammonium sulphate	Gm	
Sodium sulphite	Gm	
Sodium hydroxide	gm	
2,4dinitrophenylhydrazine hydrochloride	Gm	
Hydrochloric acid	Gm	
piperidine	MI	
Trichloroacetic acid	ml	
Silica gel (20×20cm)	Box	
Butanol	MI	
Ethanol	MI	
Alanine	Gm	
argnine	Gm	
proline	Gm	
Tyrosine	Gm	

Ninhydrin	gm	
Acetic acid	gm	
Hydrochloric acid	ml	

Name	Unt	EXPIRATION DATE
Egg albumin	gm	
Copper sulphate	Gm	
Sodium potassium tartarate	Gm	
Sodium hydroxide	Gm	
Potassium iodide	Gm	
Sodium carbonate	Gm	
Folin-ciocalteau reagent	ml	
Frish milk	MI	
Saturated litmus solution	ml	
Pancreatin	Gm	
Calcium hydroxide	Gm	
Olive oil	MI	
Sodium cholate	Gm	
Pancreatin	Gm	
Phenolphthalein	Gm	
Sodium hydroxide	Gm	
Urease(type 111 powder,sigma)	Gm	
Sodium edetate	Gm	
Phenol	Gm	
Sodium nitroprusside	Gm	
Sodium hydroxide	Gm	
Sodium hypochlorite	gm	
Ammonium molybdate	Gm	
Copper sulphate	Gm	
Sodium acetate	Gm	
Acetic acid	ml	
Metol	Gm	
Sodium sulphite	gm	
Trichloracetic acid	Gm	
Ppotassium dihydrogen phosphate	Gm	
Phosphorus acid	gm	

## Experiment:1

### The hydrolytic activity of salivary amylase on starch

#### Reagent :

1) phosphate Buffer( 0.1M , ph ( 6)

Ready solution Or a powder ( Dissolve in large Erlenmeyer flask(1L)

Dilute to mark by D.W .

2) Buffered starch substrate (0.5% in phosphate buffer )

a)5ml phosphate buffer (ph 6)+0.5 gm starch (**Max**)

b)100ml boiling phosphate buffer (ph 6)+ (a) ( **Very max +boiling** )

c)27.5 ml (b) 50ml phosphate buffer (ph 6) ( **Very max**)

3) Iodine solution (0.1N in 3% potassium iodide (KI )

a) 30gm potassium iodide + 1L D.W

b) (0.1N)Iodine :  $N = \frac{V}{m.wt} * 1L$

$$0.1 = \frac{V}{127} * 1L$$

$$V = 0.1 * 127 = 12.7 \text{ gm}$$

**(12.7gm Iodine + 1L potassium iodide)**

4) Iodine solution (0.005N in 3% KI) :

**( 1ml(0.1N)Iodine + 20ml (3%) KI )**

5)Nacl (1%):

**(1gm Nacl +100ml D.W)**

**Glass Ware :**

**2) pipette :**

plastic ( 2)

1ml : (2)

2ml : ( 2)

5ml : ( 1)

**3) Tubes ( 4 ).**

**4) Wooden sticlts .**

**5) Spottory plates.**

**6) Mechnecl Pipette .**

**Machine:**

**1) Water bath : a) 37c      b) 95c**

## Expermint:2

### Quantitative determination of amylase activity

#### **Reagent:**

#### **1) Alkaline reagent :**

1gm (3,5 dinitrosalicylic acid )+ 1gm sodium potassium tartarate +  
30gm sodium hydroxide + 1L D .W

#### **2) 1% Maltose :**

**( 1gm Maltose + 100ml D.W )**

#### **4) Phosphate buffer (0.1M) ph : 7 .**

( See Expermint:1)

#### **5) Buffered starch substrate (0.5% in phosphate buffer )**

( See Expermint:1)

#### **6) Nacl (1%):**

**(1gm Nacl +100ml D.W)**

#### **7) (2N) NaOH :**

$$N = \frac{V}{m.wt} \times 1L$$

$$2 = \left( \frac{V}{23+16+1} \right) \times 1L$$

$$2 = \left( \frac{V}{40} \right) \times 1L$$

$$V = 2 \times 40 = 80 \text{ gm}$$

**(80 gm NaOH + 1L D.W )**

***Glass Ware :***

- 1) pipette : 1ml (3)\ 10ml( 2) \5ml ( 2)
- 2) Tubes ( 8 ).
- 3) Mechnecl Pipette .
- 4) Rax.

***Machine:***

- 1) Water bath :a) 95c
- 2) Vortex .
- 3) Spectrophotometers.

### Expermint:3

## The production of pyruvate and acetaldehyde during the fermentation of glucose by yeast.

### Reagent:

1) (0.5 M) Disodium hydrogen phosphate (  $\text{Na}_2\text{HPO}_4$ ).

**(14.19gm $\text{Na}_2\text{HPO}_4$  + 100ml D.W )**

2)( 0.5M) potassium dihydrogen phosphate (  $\text{KH}_2\text{PO}_4$ ).

**(13.19 gm ( $\text{KH}_2\text{PO}_4$ )+100ml D.W)**

3) Yeast suspension (10% in  $\text{Na}_2\text{HPO}_4$ ) . (Fresh)

**( 10 gm yeast suspension +100ml (0.5M)  $\text{Na}_2\text{HPO}_4$ )**

4) Yeast suspension (10% in  $\text{KH}_2\text{PO}_4$ ) . (Fresh)

**( 10 gm yeast suspension +100ml (0.5M)  $\text{KH}_2\text{PO}_4$ )**

5) Yeast suspension ( 10% in Water) .(Fresh)

**( 10gm Yeast suspension + 100m mlD.W )**

6) Glucose (10%) .

**( 10gm Glucose + 100 ml D.W )**

7)Sodium nitroprusside ( 5%) .( Fresh )



**(5 gm Sodium nitroprusside + 100ml D.W)**

**(In Dark bottal)**

**8)Sodium hydroxide ( 10%) .**

**( 10gm NaOH + 100ml D.W)**

**9)Piperidine ( 3%) .(Fresh)**

**(3gm Piperidine +100mlD.W)**

**10) Trichloracetic acid (10%).**

**( 10 gm TCA + 100 ml D.W)**

**11) Ammonium sulphate ( 1\2 inch ).**

**( 1\2 inch = 0.5 gm Ammonium Sulphate )**

**In 2 tubes ( 15ml size ).**

**12) Ammonium hudroxide .**

**( Ready in bottal )**

**13)2,4-Dinitrophenylhydrazine hydrochloride solution in 2N hydrochloric acid .**

**a) ( 2N) Hcl : (wright on bottal(Hcl)M=63.01 \1L= 1.19 Kg)**

$$\text{Wt} = \text{Vol} \times \text{Sp.qvarity} \times 36.46$$

$$= 1000 \times 1.19 \times 0.364$$

$$= 433.16$$

$$M = Wt \backslash Mwt = 433.16 \backslash 36.46 = 11.88$$

$$M \times V = M' \times V'$$

$$11.88 \times V = 2 \times 1000$$

$$V = 2 \times 1000 \backslash 11.88$$

$$V = 168.3$$

**( 168.3 ml + 1L D.W )**

**b) (4Gm)2,4-Dinitrophenylhydrazine + ( 2N) Hc**

**14) Sodium Sulphite : take abawder .**

***Glass Ware :***

**1) pipette : plastic (1)**

**2) Tubes 15ml ( 2 ) \with centrifuge ( 2).**

**3) Mechnecl Pipette .**

**4) Rax.(2).**

**Machine:**

- 1) Water bath :a) 37c
- 3) centrifuge.
- 2) Spectrophotometers.
- 4) Vortex .

## Expermint:4

### Identification of amino acids by TLC using silica gel plates

#### Reagent:

1)Solvent :

**(7ml Butanol – 2ml acetic acid – 2ml Water )**

2) Silica gel .( *Fresh* )

**(Placed in oven When (105C) For 45 Minutes).**

3) Amino acid in (0.1M) Hcl . **(wright on bottal(Hcl M=63.01 \ 1L= 1.19 Kg)**

$$Wt = Vol \times Sp.qvarity \times 36.46$$

$$= 1000 \times 1.19 \times 0.364$$

$$= 433.16$$

$$M = Wt \backslash Mwt = 433.16 \backslash 36.46 = 11.88$$

$$M \times V = M ' \times V'$$

$$11.88 \times V = 0.1 \times 1000$$

$$V = 0.1 \times 1000 \backslash 11.88$$

$$V = 8.41$$

**( 8.42 ml Hcl + 1L D.W )**

**A) Alanine (5M)in (0.1M)Hcl**

$$M=(V \setminus m.wt) \times 1L$$

$$5 = (V \setminus 89.09) \times 1L$$

$$V= 89.09 \times 5 = 445.45gm$$

**( 445.45gm Alanine + 1L (0.1M)Hcl)**

**B) Arginine (5M) in (0.1M)Hcl**

$$M=(V \setminus m.wt) \times 1L$$

$$5 = (V \setminus 174.2) \times 1L$$

$$V= 5 \times 174.2 = 871 gm$$

**( 871 gm Arginine + 1L (0.1M)Hcl)**

**C) Proline (5M) in (0.1M)Hcl**

$$M=(V \setminus m.wt) \times 1L$$

$$5 = (V \setminus 115.3) \times 1L$$

$$V = 5 \times 115.3 = 576.5gm$$

**( 576.5 gm Proline + 1L (0.1M)Hcl)**

**D) Tyrosine ( 5M) in (0.1M)Hcl.**

$$M=(V \setminus m.wt) \times 1L$$

$$5 = (V \setminus 181.19) \times 1L$$

$$V = 5 \times 181.19 = 905.95 \text{ gm}$$

**( 905.95 gm Tyrosine +1L (0.1M)HCl)**

**4) Unknown: (1)**

**( 1ml alanine + 1ml arginine +1ml tyrosine ).**

**Unknown: (2)**

**( 1ml proline + 1ml arginine )**

**Unknown: (3)**

**( 1ml proline + 1ml alanine + 1ml tyrosine )**

**Unknown: (4)**

**( 1ml alanine + 1ml tyrosine )**

**5) Ninhydrine (0.3%)( N4) .**

**( 0.3 gm (N4) + 100ml Butanol + 3ml acetic acid )**

**( In Location spray )**

***Glass Ware :***

- 1) Micro pipette (1-5).
- 2) Tips .

***Machine:***

- 1) Hair dresser.

## Experiment:5

### Quantitative Determination of protein by the biuret & the lowry Reactions .

#### Reagent:

1) Protine Standard ( 5mg \ml ) .

**( 5 gm albumin + 100ml D.W )**

2) Biuret reagent .

a)(3gm) Copper sulphate + (9gm)Sodium potassium tartarate.

b) Sodium hydroxide ( 0.2N)

$$M=( V \setminus m.wt) \times 1L$$

$$0.2 = (V \setminus 23+16+1) \times 1L$$

$$0.2 ( V \setminus 40 ) \times 1L$$

$$V= 0.2 \times 40 = 8 \text{ gm}$$

**( 8 gm NaOH + 1L D.W)**

c) Potassium Iodide ( 5 gm )

d) 250 ml (0.2N)NaOH + (a)+ mix well in Machine + (c) + mix well in Machine + 250 ml NaOH+ placed in Flask standard 1L + Dilute to mark by D.W .

3) Unknown: (1)

**In Class :( 9.6 ml Protine standard +30.4ml D.W )**



In *finel exam* :(160ml protine standard+ 40ml D.W )

4) Alkaline reagent (1) : ( 2% Na<sub>2</sub>CO<sub>3</sub> in 0.1N NaOH )

a) NaOH (0.1N)  $N = (V \setminus m.wt) \times 1L$

$$0.1 = (V \setminus 1+16+23) \times 1L$$

$$V = 0.1 \times 40 = 4 \text{ gm} \quad (4 \text{ gm} + 1L \text{ D.W})$$

b) 2 gm (Na<sub>2</sub>CO<sub>3</sub>) + 100 ml (0.1M) NaOH

5) Alkaline reagent (2)(0.5% CuSO<sub>4</sub> in 1% Na<sub>1</sub>K tartrate) .(Fresh )

a) Copper sulphate 0.5% ( 0.5 gm CuSO<sub>4</sub> + 100ml D.W )

b) Sodium potassium tartrate 1% ( 1 gm Na<sub>1</sub>K tartrate + 100ml D.W )

c) 50 ml (a) + 50 ml ( b) (Fresh)

6) Alkaline reagent (3) (Fresh)

50 ml Alkaline reagent (1) + 1ml Alkaline reagent ( 2- c)

7) Folin Ciocalteau reagent .

(100ml Folin Ciocalteau reagent+ 100ml D.W)

8) Protine standard .

( 0.2 gm Albumin + 100ml D.W)

9) Unknown: (2)

(50 ml Protine standard + 50 ml D.W )

**Glass Ware :**

- 1) pipette : 5ml (3) \2ml (4) \1ml (4)
- 2) Tubes 15ml (8 ) \small( 8).
- 3) Mechnecl Pipette .
- 4) Rax .

**Machine:**

- 1) Water bath :a) 37c
- 2) Spectrophotometers.
- 3) Vortex .

## Experiment:6

### Enzymatic Digestion of fat by pancreatic lipase .

#### **Reagent:**

1) Homogenised whole milk preferably (*fresh.*)

2) saturated Litmus solution :

(0.5gm Litmus + 100ml D.W ) Or (ready in bottle)

3) Pancreatic solution :

a) unboiled : 50ml Ethanol absolute + 50ml ml water + 0.8 gm  
pancreatin ( pH = 8)

b) boiled : 50ml (a) in water bath about 10 minutes.

4) Calcium hydroxide (1%) .

1gm Calcium hydroxide + 100ml D.W .

**Glass Ware :**

- 1) Pipette : 1) plastic pipette (2) .2ml ( 2) . 1ml (1) .
- 2) Test tubes ( 2) .
- 3) Mechanical Pipette .
- 4) Rack .

**Machine:**

- 1) Water bath.45c

## Expermint:7

### Effect of Bile salts on pancreatic Lipase activity .

#### **Reagent:**

- 1) Vegetable oil (olive ): or any oil
- 2) Sodium cholate solution .  
10 gm sodium cholate + 100ml D.W
- 3) D.W
- 4) Pancreatin solution (Lipase ) .  
See Expermint (6) .(unboiled)
- 5) Phenolphthalein indicator (1%) .  
2gm phph + 200ml ethanol + 120ml D.W
- 6) Sodium hydroxide ( 0.05N) .  
 $M=(V \setminus m.wt) \times 1L$   
 $0.05 = (V \setminus 23+16+1) \times 1L$   
 $0.05 = (v \setminus 40) \times 1L$   
 $V = 0.05 \times 40 = 2gm$   
( 2gm NaOH + 1L D.W )

**Glass Ware :**

- 1) Pipette : 1) plastic pipette (1) .2ml ( 1) . 1ml (3) .
- 2) Test tubes ( 2) .
- 3) Mechanical Pipette .
- 4) Flask (2) .
- 5) burette .
- 6) Funnel .

**Machine:**

- 1) Water bath 45 c.

## Experiment:8

### Determination of serum urea

#### **Reagent:**

- 1) Buffered urease solution :
  - a) In large Erlenmeyer flask 100ml 1gm sodium edetate+ 80 ml D.W(pH: 6.5) .
  - b) (a) + 0.1 gm urease . Dilute to mark by D.W.
- 2) Phenol-sodium nitroprusside solution.  
In large Erlenmeyer flask 1L 50 gm phenol + 0.25 gm sodium nitroprusside + Dilute to mark by D.W .
- 3) Sodium hydroxide – sodium hypochlorite (NaOCl)  
25gmNaOH + 2.1 gm sodium hypochlorite + 1L D.W
- 4) Urea solution .  
0.12 gm urea + 100ml .....

#### **Glass Ware :**

- 1) Pipette : 1) plastic pipette (1) .2ml ( 1) . 1ml (3) .
- 2) Test tubes ( 2) .
- 3) Mechanical Pipette

**Machine:**

- 1) Water bath 37 c.
- 2) Spectrophotometers.



## Expermint:9

### The colorimetric estimation of inorganic phosphate .

#### Reagent:

1)Ammonium molybdate (5%) .

5gm Ammonium molybdate + 100gm D.W .

2) Copper acetate buffer ( PH: 4).

a) (2N) acetic acid : (*wright on bottal( HC<sub>2</sub>H<sub>3</sub>O<sub>2</sub>)M=60.05 \ 1L= 1.05 Kg*).

$$\begin{aligned} \text{Wt} &= \text{Vol} \times \text{Sp.qvarity} \times 60.05 \\ &= 1000 \times 1.05 \times 60.05 \ / \ 100 \\ &= 630 \end{aligned}$$

$$M = \text{Wt} \ / \ \text{Mwt}$$

$$M = 630 \ / \ 60 .05 = 10.49$$

$$M \times V = M' \times V'$$

$$10.49 \times V = 2 \times 1000$$

$$V = 2 \times 1000 \ / \ 10.49 = 190.6$$

**( 190.6 ml HC<sub>2</sub>H<sub>3</sub>O<sub>2</sub> + 1L D.W )**

b)2.5gm copper sulphate + 46 gm sodium acetate .

c) ( a) + ( b) **3)Metol – sodium sulphate solution .(Fresh)**

a) sodium sulphite (10%) .

10gm sodium sulphite + 100ml D.W .

b) 2gm Metol + (a) . (in dark bottle )

4) TCA (5%) :

(5 gm +100ml D.w ) .

5) stock phosphate solution .

a) 0.438 gm potassium dihydrogen phosphate (pdp) + 100ml D.W

b) 0.1 gm phosphorus acid + (a) .

6) Working phosphate .

1ml stock phosphate solution + 100ml 5% TCA .

7) Unknown :

In class : (50 ml working phosphate + 50 ml 5% TCA)

In Finel exam : (80 ml working phosphate+20ml 5%TCA)

**Glass Ware :**

1) Pipette : 1ml ( 6) 2ml ( 1)

2) Test tubes (7) .

**Machine:**

1) Spectrophotometers



كلية العلوم الطبية التطبيقية

قسم علوم المختبرات الاكلينيكية