Experiment 6: Oral Glucose tolerance test (OGTT)

Procedure:

- Set up 9 test tube the Please follow the table:

| Chemical | Fasting | | After 1 hour | | After 2 hour | | standard | | Blank |
|---------------------|---------|--------|--------------|--------|--------------|--------|----------|--------|--------|
| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | |
| Fasting | 0.1 ml | 0.1 ml | - | - | - | = | - | - | - |
| After 1 hour | - | - | 0.1 ml | 0.1 ml | - | = | - | - | - |
| After 2 hour | - | - | - | - | 0.1 ml | 0.1 ml | | | - |
| Working standard | - | - | - | - | - | - | 0.1 ml | 0.1 ml | |
| Distilled water | - | - | - | - | - | - | - | - | 0.1 ml |
| O-toluidine reagent | 7 ml | 7 ml | 7 ml | 7 ml | 7 ml | 7 ml | 7 ml | 7 ml | 7 ml |

 Cover all tube by aluminum foil then use Vortex to mix the contents, then place all tube in water bath at boiling degree for 5 min.

Read the absorption of 6 tubes against blank at 630 nm. And record your result **Results:**

| Test t | tube | Absorption at 630nm | Mean |
|--------------|--------|---------------------|------|
| Fasting | Test 1 | | |
| | Test 2 | | |
| After 1 hour | Test 3 | | |
| | Test 4 | | |
| After 2 hour | Test 5 | | |
| | Test 6 | | |
| Working | Test 7 | | |
| standard | Test 8 | | |

Calculation:

| Amount of glucose in plasma = | Mean Ab Test X C Mean Ab Standard Concentration of Stand | _ | = mg/dl |
|-------------------------------|--|---|---------|
| Mean Absorbance of Standa | rd = | | |
| Mean Absorbance of the Fas | ting sample = | | |
| Mean Absorbance of the Aft | er 1 hour sample = | | •••••• |
| Mean Absorbance of the Aft | er 2 hour sample = | | •••••• |
| Glucose concentration in | Fasting plasma= | | |
| | | | |
| Glucose concentration af | ter 1 hour= | | |
| Glucose concentration af | ter 2 hour= | | |