




# The effect of incubation time on the rate of an enzyme catalyzed reaction



# Objectives

- To monitor the progress of an enzyme catalyzed reaction (Acid phosphatase).
- To determine the initial rate of the reaction ( $V_i$ ).
- **Important terms and points:**





# Enzyme kinetics

- The central approach for studying the mechanism of an enzyme-catalyzed reaction is to study enzyme kinetics.
- It **determines** the rate of the enzymatic reaction (**velocity**) and its changes in response with the changes in **parameters** such as **substrate concentration**, **enzyme concentration**, **pH**, **temperature**.


# Measuring the rate of reaction (velocity)

The progress of an enzyme catalyzed reaction may be followed by measuring either the

- amount of substrate consumed,
  - or Amount of product formed
- } per unit time


**Velocity (V)** = rate of reaction = change in [P] or [S] per unit time

Unit :  $\mu\text{moles/minute}$




# Initial velocity

- The initial rate of reaction is used in the study of enzyme kinetics
- The initial rate of reaction ( $V_i$ ) is measured as the slope at the origin (time= 0).



# Why we measure initial velocity ?

- **The rate of the reaction is highest at time zero** and decreases with increasing time, eventually falling to zero itself, reaching a **plateau**.
- This usually occurs either when all the substrate is used up or when equilibrium is reached.
- **The Solution is** to Measure  $V$  at very early times in reaction, before  $[S]$  decreases significantly
- The initial rate of reaction,  $v_i$ , measured as the **slope** to the above curve at the origin (time= 0)



# Acid phosphatase




- In this experiment, you will measure the velocity of the reaction catalyzed by purified acid phosphatase (**ACP**) from [wheat germ](#).
- Acid phosphatase is a **phosphatase** that acts on **monoesters** of orthophosphoric acid. It does **not act** on phosphoric diesters or triesters.
- Acid phosphatase is an important enzyme that plants use to obtain their energy from phosphate. Acid phosphatase's function in plants is to hydrolyze phosphate esters during energy metabolism.

# Acid phosphatase from wheat germ





1. The enzyme catalyzes the hydrolysis of **p-nitrophenyl phosphate** to **inorganic phosphate** and **p-nitrophenol**, under **acidic conditions (pH=5.7)**, with **optimum temperature 37°C**

- $Mg^{++}$  ion act as an activator for the enzyme



# Principal of the enzyme assay in vitro



2. Both p-nitrophenyl phosphate and p-nitrophenol are colorless at acidic pH values.
3. Addition of **alkaline solution**, p-nitrophenol is converted to a **p-nitrophenolate** (yellow color) and concentration can be measured at **405 nm**.

**Note:** Since ACP works under acidic condition, adding alkaline will cause the **enzymatic reaction to stop**.

## Method: The effect of incubation Time

- The reaction is stopped after 5 minutes intervals (by Addition of KOH ), hence called **Fixed-time assay**




# Method:

- Prepare a series of **seven reaction tubes** labeled 0 through 30 minutes at 5 minute intervals (Blank, 5, 10 ... minutes).
- Follow the following addition protocol for **all the tubes**:

Add in each tube:


- 0.5 ml of buffer
- 0.5 ml of  $\text{MgCl}_2$
- 0.5 ml of pNPP (S)
- 5 ml of water



Place in a water bath maintained at 37 °C for 5 minutes

To start the reaction → add 0.5ml of E

To stop the reaction → add 0.5ml of KOH



Start at (min)	0	0	2	4	6	8	10
Stop at (min)	0	5	12	19	26	33	40



Add KOH to blank ((FIRST)), to prevent the reaction from happening.

After all the reactions have been terminated, determine the absorbance at **405 nm** for each sample against blank.

<b>Tube</b>	<b>5</b>	<b>10</b>	<b>15</b>	<b>20</b>	<b>25</b>	<b>30</b>
<b>Start at (min)</b>	<b>0</b>	<b>2</b>	<b>4</b>	<b>6</b>	<b>8</b>	<b>10</b>
<b>Stop at (min)</b>	<b>5</b>	<b>12</b>	<b>19</b>	<b>26</b>	<b>33</b>	<b>40</b>



<b>Time (min)</b>	<b>Tube</b>	<b>Addition 0.5 ml of</b>
0	5	Enzyme
2	10	Enzyme
4	15	Enzyme
5	5	<b>KOH</b>
6	20	Enzyme
8	25	Enzyme
10	30	Enzyme
12	10	<b>KOH</b>
19	15	<b>KOH</b>
26	20	<b>KOH</b>
33	25	<b>KOH</b>
40	30	<b>KOH</b>



To convert the time table to an easier way try the following



# Results

Time (min)	Absorbance at 405nm	[P] (μmoles)
0		
5		
10		
15		
20		
25		
30		


$$[P] = (A \times 10^6) / E \times L = \mu\text{mole}$$

E= extension

coefficient= $18.8 \times 10^3$

A= absorbance

L= path length (1 cm)



Progress Curve

$$V_1 = \frac{y_2 - y_1}{x_2 - x_1}$$

# Discussion

- Draw the graph [do not forget the title and the units]
  - X axis = Time
  - Y axis = [P]
  - Calculate from the graph the value initial velocity
- Principle
- Comment on the graph that you get in details



# Questions

- How is the initial velocity of an enzymatic reaction determined?
- What are the factors that influence the rate of enzyme catalyzed reactions?