[BCH 322]

Method of Enzyme Assay

1-	Discontinuous	assav ((col	lorimetric	assav)
_	2 000001000000		,	0.0.000.00	

7.4	-41	۱ م	1	
/VI	ett	าก	П	

Prepare each tube as follows:

	BLANK	SAMPLE		
ALT Reagent	0.5 ml	0.5 ml		
Pre-warm at 37 °C for	r <u>5 minutes</u> and ad	d:		
Distilled Water	0.1 ml	-		
Serum Sample	-	0.1 ml		
Mix, and incubate at 3	37 °C for exactly 3	0 minutes, and add:		
Color Reagent	0.5 ml	0.5 ml		
(DNPH)				
Mix, and return at 37 °C for exactly 10 minutes, then add:				
Color Developer	5.0 ml	5.0 ml		
(NaOH)				
Mix, and return to 37 °C for exactly <u>5 minutes</u> . Read absorbance of all tubes at 546nm against blank.				

Result

Absorbance at 546 nm = \dots

ALT (SGPT) activity of serum sample (<u>from graph</u>)=U/L

2-	Continuous	assav	(UV/	kinetics)
----	-------------------	-------	------	----------	---

Method

Type of cuvette used -----

Prepare the following:

ALT Reagent	3 ml			
Pre-warm at 37°C for 3 minutes and add				
Serum Sample	0.2 ml			
Mix and incubated at 37 °C for 1 minute, then read absorbance (at 340 nm against				
distilled water) every minute for 3 minutes) and determine $\Delta A/min$				

Choose the following on the spectrophotometer:

2) Applications → 2) Simple Kinetics → wave length (340 nm) → 1) Seconds → Duration (180 sec = 3 min) → Intervals (60 sec= 1 min) → Print Data Table (off) → Press start (2 times)

Result

Time	Absorbance 340nm		$\Delta A/min = ((A1-A2)+(A2-A3))/2$
1 min	A1		
2 min	A2		
3 min	A3		

Calculation

ALT Activity (U/L) = $\Delta A/min x 1768$

 $ALT\ Activity =U/L$