


Antibiotic sensitivity Testing

Antimicrobial Susceptibility Testing for Bacteriology Laboratories



Manal AL khulaifi

- **Antibiotic:** Chemical produced by a microorganism that kills or inhibits the growth of another microorganism

- **Types of antibiotic:-**
 - **Narrow spectrum antibiotic** active against either Gram +ve bacteria only or gram -ve bacteria only

 - **Broad spectrum antibiotic** active against both Gram's +ve and -ve bacteria

Antibiotic Sensitivity Testing

- Antibiotic sensitivity testing is used to determine the susceptibility of the microorganism to various antimicrobial agents.

Sensitivity testing techniques:

Disc Diffusion Method.


Serial Dilution Method

MIC (Minimum inhibitory concentration).

MECHANISMS OF ACTION OF ANTIBACTERIAL DRUGS

■ Mechanism of action include:

- Inhibition of cell wall synthesis
- Inhibition of protein synthesis
- Inhibition of nucleic acid synthesis
- Inhibition of metabolic pathways
- Interference with cell membrane integrity



Classification of antibiotics based on mode of action

- **Inhibit Cell synthesis**
- Penicillin (P) Bacitracin (Ba) –Vancomycin (Va)
- **Inhibit Protein synthesis**
- Erythromycin (E) - Gentamycin (GN)
- **Inhibit nucleic acids**
- DNA/RNA Rifampicin (RF) –Ciprofloxacin (CIP)
- **Inhibit of selective permeability of cell membrane**
- Polymyxins
- **Inhibit bacterial metabolism**
- SXT) sulphamethoxazole/trimethoprin(

Routine Susceptibility Tests

- **Disk diffusion**
 - (Kirby Bauer)
 - (Modified Stock Technique)
- **Broth micro-dilution MIC / MBC**
- **Etest**

Disc Diffusion Method

- The effectiveness of an antibiotic in this technique is based on the size of the zone of inhibition that surrounds a disc that has been impregnated with a specific concentration of the agent.


- **Advantages:**

- ✓ Rapid
- ✓ Accurate
- ✓ Inexpensive



ANTIMICROBIAL SUSCEPTIBILITY TESTING

- Probably the most widely used testing method is the disk-diffusion method, also known as the Kirby-Bauer test.



- **The recommended medium →**
Mueller Hinton agar

The inoculum:

The turbidity of a broth culture / saline suspension of the test organism has to match a defined standard

“0.5 McFarland” (a barium sulphate standard)

Disc Diffusion method

■ Materials:

- Cultures of bacteria
- Filter paper disc
- Antibiotic solutions (Vancomycin, Augmentin, Ceftazidime, Azactam)
- 20 ml melted nutrient agar
- Petri-dish
- Sterile 1ml pipette


Prepare inoculum suspension



Select colonies



Mix well



**Standardize inoculum
suspension**

Manal AL-Khulafi

Swab plate






Add disks

Incubate overnight




Procedure:




Results:

- Measure the diameter of each inhibition zone
- * The diameter of the inhibition zones are directly proportional to the susceptibility of the microorganism to the antibiotics.




Microbial Antibiotic Susceptibility Test

(The agar diffusion test)




Measure Zones


Transmitted
Light



Manal AL khulaifi




Reflected
Light



Susceptibility Testing Methods






**Innoculate
MH plate**



**Place disks
on agar plate**




**Incubate plate
18-24 hr, 35 C
Measure and
record zone of
inhibition around
each disk**



Disk Diffusion Susceptibility Testing


Improper agar & disk placement



Use Mueller Hinton agar

Manal AL khulaifi

Mueller Hinton agar &
good disk placement



MIC:

It is the **lowest concentration of the antimicrobial agent that inhibits the growth of the test organism** but not necessarily kills it.

MBC (minimum bactericidal conc.):

It is the **lowest concentration of the antimicrobial agent that kills the test organism.**


Clinical Conditions when MICs are Useful


- Endocarditis
- Meningitis
- Septicemia
- Osteomyelitis
- Immuno suppressed patients (HIV, cancer, etc.)


Methods for determination of MIC

1. **Serial dilution** method
2. **Agar diffusion** method
3. **E- test**

■ Procedure:









Incubate the plate **Uninverted** at 37°C for 24h

Results:

- Measure the diameter of each inhibition zone
- * The diameter of the inhibition zones are directly proportional to the conc. of the antibiotic.


Zone diameter

Well diameter




Determination of MIC by the serial dilution method


■ Procedure:



Prepare inoculum suspension



Microdilution MIC tray




E-test

Antimicrobial Gradient Testing E-test®

A plastic coated strip contains a gradient of antibiotic concentrations and the minimal inhibitory concentration is read from a scale printed on the strip


**Read plates
after
recommended
Incubation**



**Read MIC
where
intersects
scale**

the E test


configuration






*A Penicillin Resistant Pneumococcus (PRP)
with an MIC of 3 µg/ml.*



E-Test



Manal AL Khulafni



E-test for *P.aeruginosa*

