Cloning and rDNA (I)

Dr. Aws Alshamsan
Department of Pharmaceutics

Office: AA87


Tel: 4677363


aalshamsan@ksu.edu.sa


Objectives of this lecture

By the end of this lecture you will be able to:


- Understand the concept of recombinant DNA technology
- 2. Define terms related to cloning and DNA recombination
- 3. Realize the value of recombinant DNA technology









Human Insulin Production



Isolation of Gene of Interest

- The gene of interest is a small segment in the large DNA molecule
- Restriction enzymes (endonucleases) cut the DNA double helix at specific sites
- Different restriction enzymes have different sequence specificities (palindrome)





Examples of Palindromes

Noon

Murder for a jar of red rum

مودته تدوم لكل هول *** وهل كل مودته تدوم

Enzymes with wobbling cuts give sticky ends

• HindIII - leaves 5' overhangs


• Kpnl leaves 3' overhangs


 Enzymes that cut at same position on both strands give blunt ends

• Smal

5' --cccggg-- 3'
3' --gggccc-- 5'

3' --ggg ccc-- 5'







DNA Ligase

 DNA ligation is the act of joining together DNA strands with covalent bonds with the aim of making new viable DNA or plasmid

T4 DNA ligase has the unique ability to join sticky and blunt

ended fragments





DNA Ligase

- Sticky ends that are complementary can be ligated together even if they are produced by different restriction enzymes.
- Sticky ends that are not complementary cannot be ligated together.

GATCC-BamHI-G

> G--CCTAG

BgIIIGATCT-**-**A


> -TCTAG **A**-

Result -GGATCT-


No longer palindromic, so -CCTAGA- not cut by BamHI or Bg/III


DNA Ligase can also join blunt ends


DNA fragments with blunt ends generated by different enzymes can be ligated together (with lower efficiency), but usually cannot be re-cut by either original restriction enzyme.



- Ligations that re-constitute a Smal or Dral site (CCCGGG or AAATTT) can be re-cut by Smal or Dral.
- Mixed ligation products (CCCTTT + AAAGGG) <u>cannot</u> <u>be re-cut by Smal or Dral.</u>







You are now able to:

- ✓ Understand the concept of recombinant DNA technology
- ✓ Define terms related to cloning and DNA recombination
- ✓ Realize the value of recombinant DNA technology