Titration curve

of polyprotic acids

Dissociation of polyprotic acids

- A polyprotic acid has more than one proton per molecule, thus it ionizes in successive steps.
- Example: H₂A a "polyprotic acid", diprotic acid


$$H_2A \xrightarrow{K_{a1}} H^+ + HA^- \xrightarrow{K_{a2}} H^+ + A^{-2}$$

$$K_{a1} = \frac{[H^+][HA^-]}{[H_2A]}$$
 $K_{a2} = \frac{[H^+][A^{-2}]}{[HA^-]}$

Dissociation of polyprotic acids cont'ed

- $ightharpoonup K_{a1}$ is larger than the K_{a2}
- The pH of H₂A solution is mainly dependant on the first ionization step.

Titration curve of H₃PO₄ a poly protic acid



Titration curve of H₃PO₄ a poly protic acid cont'ed

$$H_{3}PO_{4} \xrightarrow{K_{1}} H^{+} + H_{2}PO_{4}^{-}$$

$$H_{2}PO_{4}^{-} \xrightarrow{K_{2}} H^{+} + HPO_{4}^{-2}$$

$$HPO_{4}^{-2} \xrightarrow{K_{3}} H^{+} + PO_{4}^{-3}$$

$$K_{a1} = \frac{[H^{+}][H_{2}PO_{4}^{-}]}{[H_{3}PO_{4}]}$$

$$K_{a2} = \frac{[H^{+}][HPO_{4}^{-2}]}{[H_{2}PO_{4}^{-}]}$$

$$K_{a3} = \frac{[H^{+}]PO_{4}^{-3}}{[HPO_{4}^{-2}]}$$

Titration curve of H₃PO₄ a poly protic acid cont'ed

$$K_{a1} > K_{a2} > K_{a3}$$

Example

How many ml of 0.1 M NaOH are required to titrate 100 ml of 0.1 M H_3PO_4 ?

No. of moles of $H_3PO_4 = V \times M = 0.1 \times 0.1 = 0.01$ mole

 \triangleright No. of moles of H⁺ = V \times M = 0.01 \times 3 = 0.03 mole

Titration curve of H₃PO₄ a polyprotic acid cont'ed

Thus: we need 0.03 moles of NaOH to titrate the acid

```
M = no. of moles / V
V = no. of moles / M
V = 0.03 / 0.1
V = 0.3 L = 300 ml
```

Example 1

- ▶ 1.025 g of anhydrous sodium acetate is dissolved in 100 ml of 0.25 M acetic acid CH₃COOH. Calculate:
 - A) The pH of the final solution.
 - B) The molarity of the final solution (resulting buffer)

Mwt of anhydrous sodium acetate = 82; $pK_a = 4.7$

```
A) pH = pK<sub>a</sub> + log \frac{[A^-]}{[HA]}
no. of moles of A<sup>-</sup> in buffer = wt / Mwt
= 1.025 / 82
= 0.0125 moles
```

no. of moles of HA in buffer =
$$M \times V$$

= 0.25×0.1
= 0.025 moles

$$pH = 4.7 + log (0.125 / 0.25)$$

 $pH = 4.39$

B) The molarity of the buffer = the molarity of HA + the molarity of A⁻
 Buffer molarity = 0.25 + 0.125 = 0.375 M

<u>OR</u>

```
No. of moles of buffer = no. of moles of HA + the no. of moles of A<sup>-</sup>
No. of moles of buffer = 0.025 + 0.0125 = 0.0375 moles
Molarity of buffer = no. of moles / V in L
Molarity of buffer = 0.0375 / 0.1 = 0.375 M
```

Example 2

5 L of 0.1 M phosphate buffer with a pH = 12.32 was prepared from Na₃PO₄ and Na₂HPO₄.

Calculate the weight in grams of each component which was used to prepare the buffer, $pK_a = 12$.

SELF SOLVE!

Titration Curve

of amino acids

Titration curves

- It is a curve that monitors the pH of a solution as amounts of alkali or acid is added.
- Amino acids are simple weak polyprotic acids.
- Neutral amino acids (as alanine, glycine) are treated as diprotic acids.
- Acidic amino acids (as aspartic or glutamic acid) and Basic amino acids (as lysine or histidine) are treated as triprotic acids.


Titration curve of alanine

At point a:

- Before titration
- NH⁺₃CHRCOOH
- The net charge = +1

At point b:

- \circ Pk₁ \cdot = pH
- Here it has buffering capacity
- NH⁺₃CHRCOOH = NH⁺₃CHRCOO⁻
- The net charge = +0.5



Titration curve of alanine cont'ed

At point c:

- Isoelectric point
- pI = pH; to calculate pI: from the plot or pI = (pK₁`+pK₂`) \ 2
- NH⁺₃CHRCOO⁻ a zwitter ion
- The net charge = 0

At point d:

- Pk₂` = pH
- Here it has buffering capacity
- NH⁺₃CHRCOO⁻ = NH₂CHRCOO⁻
- The net charge = -0.5

Titration curve of alanine cont'ed

At point e:

- End of titration
- NH₂CHRCOO⁻
- The net charge = −1
- This curve has two flat zones, at point b and d; meaning it has two ionized groups.