

# Blood Biochemistry BCH 220 [Practical] Lab (8) Clotting Time, Bleeding Time and Prothrombin Time

## Coagulation

- **Coagulation** is a complex process by which <u>blood forms clots</u>.
- It is an important part of **hemostasis** (the cessation of blood loss from a damaged vessel).
- Disorders of coagulation can lead to an increased risk of bleeding (hemorrhage) or clotting
  - (thrombosis). Bleeding Clotting Hemostasis

### Hemostasis

- Hemostasis is a complex of homeostatic reactions, which result in <u>arrest of bleeding</u> from damaged blood vessels.
- Maintained in the body via three mechanisms:
  - 1. Vascular spasm
  - Damaged blood vessels constrict.
    - 2. Platelet plug formation

Platelets adhere to damaged endothelium to form platelet plug (primary hemostasis).

#### **3. Blood Coagulation**

Clots form upon the conversion of fibrinogen to fibrin (secondary hemostasis).





(a) The general steps of clotting

## **Clotting Cascade**

 A cascade is a mechanism in which enzymes activate other enzymes sequentially usually leading to an amplification of an initial signal.

### **Pathways**

- **1.** Extrinsic, which normally is triggered by trauma.
- 2. Intrinsic, which begins in the bloodstream and is triggered by internal damage to the wall of the vessel.
- Initially independent, then they converge on common pathway leading to the formation of a **fibrin clot**.
- Each of these pathways leads to the conversion of factor X (inactive) to factor Xa (active).

## What triggers extrinsic and intrinsic pathways?

#### Extrinsic

- Damage to tissue **outside** the blood vessel.
- This pathway acts to clot blood that has **escaped** from the vessel.

#### Intrinsic

- Damage to blood vessel wall.
- Coagulates the blood **within** the damaged vessel.
- Triggered by elements that **lie within** the blood itself.

### **Clotting Cascade**



**Pause and Think why extrinsic pathway is completed** 

in seconds whereas intrinsic pathway takes few minutes?



## **Clotting Time**

- Test for <u>intrinsic system.</u>
- Simple test but takes time and rarely done now.



#### **Method**:

- Venous blood is taken and placed on glass test tube at 37°C and it observed at time intervals until clotting occurs.
- Normal blood takes **5 10 min** to clot.
- Longer periods  $\rightarrow$  Coagulation defects (e.g. Hemophilia).

### **Clotting Time - Capillary Method**





# **Bleeding Time (BT)**

Provides assessment of platelet count and function.

#### Method:

• It is determined by noting time at which blood coming out a small cut, no longer forms a spot on a

piece of filter paper placed in contact with cut surface.

• The normal range from **2-4 min**.



# **Prothrombin Time (PT)**

Measures effectiveness of the <u>extrinsic pathway</u>

#### Method:

- An excess of tissue factor and Ca<sup>2+</sup> ions are added to diluted plasma containing citrate (anticoagulant) and then the time taken for the mixture to clot is measured.
- Normal value  $\rightarrow$  10-15 secs
- High  $PT \rightarrow low levels of thrombin$
- Results from liver disease due to deficiency of prothrombin, fibrinogen,
  V, VII and X factors.

#### **Prothrombin Time [PT]**

